

# **Agilent N77xx Series**

N7751A 1-Channel Optical Attenuator and Power Meter N7752A 2-Channel Optical Attenuator and Power Meter

N7761A 1-Channel Optical Attenuator

N7762A 2-Channel Optical Attenuator

N7764A 4-Channel Optical Attenuator

N7766A 2-Channel Multimode Optical Attenuator

N7768A 4-Channel Multimode Optical Attenuator

**User's Guide** 



### **Notices**

© Agilent Technologies, Inc. 2010

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

### **Manual Part Number**

N7751-90B01

#### **Edition**

Fourth edition, July 2010

Agilent Technologies Deutschland GmbH Herrenberger Str. 130 71034 Böblingen, Germany

### **Warranty**

This Agilent Technologies instrument product is warranted against defects in material and workmanship for a period of one year from date of shipment. During the warranty period, Agilent will, at its option, either repair or replace products that prove to be defective.

For warranty service or repair, this product must be returned to a service facility designated by Agilent. Buyer shall prepay shipping charges to Agilent and Agilent shall pay shipping charges to return the product to Buyer. However, Buyer shall pay all shipping charges, duties, and taxes for products returned to Agilent from another country.

Agilent warrants that its software and firmware designated by Agilent for use with an instrument will execute its programming instructions when properly installed on that instrument. Agilent does not warrant that the operation of the instrument, software, or firmware will be uninterrupted or error free.

### **Limitation of Warranty**

The foregoing warranty shall not apply to defects resulting from improper or inadequate maintenance by Buyer, Buyer-supplied software or interfacing, unauthorized modification or misuse, operation outside of the environmental specifications for the

product, or improper site preparation or maintenance.

No other warranty is expressed or implied. Agilent Technologies specifically disclaims the implied warranties of Merchantability and Fitness for a Particular Purpose.

### **Technology Licenses**

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

### **Restricted Rights Legend**

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as "Commercial computer software" as defined in DFAR 252.227-7014 (June 1995), or as a "commercial item" as defined in FAR 2.101(a) or as "Restricted computer software" as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies' standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

### **Exclusive Remedies**

The remedies provided herein are Buyer's sole and exclusive remedies. Agilent Technologies shall not be liable for any direct, indirect, special, incidental, or consequential damages whether based on contract, tort, or any other legal theory.

### **Assistance**

Product maintenance agreements and other customer assistance agreements are available for Agilent Technologies products. For any assistance contact your nearest Agilent Technologies Sales and Service Office.

#### Certification

Agilent Technologies Inc. certifies that this product met its published specifications at the time of shipment from the factory. Agilent Technologies further certifies that its calibration measurements are traceable to the United States National Institute of Standards and Technology, NIST (formerly the United States National Bureau of Standards, NBS) to the extent allowed by the Institutes's calibration facility, and to the calibration facilities of other International Standards Organization members.

### **ISO 9001 Certification**

Produced to ISO 9001 international quality system standard as part of our objective of continually increasing customer satisfaction through improved process control.

### **Safety Notices**

### **CAUTION**

A **CAUTION** notice calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

### WARNING

A WARNING notice calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

# **Contents**

### 1 Getting Started

```
Overview
Agilent N7744A and N7745A Optical Power Meter 7
Agilent N775xA Optical Attenuator and Power Meter and N776xA Optical Attenuator
Agilent N7711/4A Single and 4-Port Tunable Laser System Source
                                                              10
Initial Inspection 10
Safety Considerations 11
  Safety Considerations - Overview 11
  Safety Symbols
                   11
  General 12
  Operating Environment 12
  Line Power Requirements
  Input/Output Signals
  Line Power Connectors
  Instrument Markings
AC Line Power Supply Requirements
  AC Line Power Supply Requirements - Overview
                                                14
  Line Power Requirements 15
  Line Power Cable 15
  Changing the Fuse
Operating and Storage Environment 16
  Operating and Storage Environment - Overview
                                                16
  Temperature 17
  Humidity 17
  Altitude 17
  Pollution Protection 17
  Instrument Cooling
```

```
Input and Output Connectors 18
  Input and Output Connectors - Overview
  Front Panels - Multiport Power Meters
  Front Panels - Optical Attenuators 19
  Front Panels - Tunable Laser Source
  Rear Panel 21
Making Optical Connections 22
  Making Optical Connections
  Making Quad-adapter Connections on the Agilent N7744/45A Multiport Power
  Meter 23
Optical Output
Electrical Connectors
LAN Interface 25
  LAN Interface
  Selecting a LAN Network 26
  Private LAN Considerations 26
  Security 26
  Performance
                26
  Reliability 27
  IP Address Availability 27
  Site LAN Considerations 27
  Connecting the LAN Cable 27
GPIB Interface 30
  GPIB Interface
  Cables and Adapters 30
  Connector 30
  GPIB Logic Levels 31
USB Interface 31
Powering Up the Instrument 32
  Powering Up the Instrument 32
  Front Panel Indicators 32
Connecting to the Instrument 33
  Connecting to the Instrument - Overview
                                         33
  Connecting over USB
  Finding the IP Address of an instrument
  Connecting over LAN
Claims and Repackaging
  Claims and Repackaging
  Return Shipments to Agilent Technologies
```

		Connecting to an Instrument 40 Viewing Measurements 42 Saving and Recalling Configurations 43 Controlling the Instrument Refresh 44
2	Attenuation	
		Attenuation 45
		Controlling the Optical Attenuator 45
		Setting Attenuation 47  Setting the Number of Digits displayed 47  Selecting the Power Units 47  Setting the Attenuation 48  How to Remove Electrical Offsets (Zero) 49  Enabling and Disabling the Output 51  How to Set a Power Level 51  How to Set the Wavelength 52  How to Set the Averaging Time 52  Compensating for Wavelength Dependencies in your test setup 52  Applying an Offset Table 52
3	Specifications	and Regulations Compliance
		Definition of Terms 55
		Literature 74
		Technical Specifications - Optical Variable Attenuator and Power Meter Optical Attenuator 75 Optical Power Meter 75 General Characteristics 76
		Declarations of Conformity - Optical Multiport Attenuator and Power Meter 78
		Regulations Information 79 EMC Canada 79 Acoustic Noise Information Germany 79
		Remove all doubt 79 Remove all doubt 79 Agilent E-mail Updates 79 Agilent Direct 80 Agilent Open 80 LXI 80 www.agilent.com 81

N77xx Viewer Common User Interface Functions 40

#### 4 Accessories

Accessories - Overview 83
Instrument and Options - Agilent N7751A, N7752A, N7761A, N7762A, N7764A, N7766A, N7768A 84
Cables 84

#### 5 Maintenance and Troubleshooting

```
Maintenance and Troubleshooting - Overview
                                            87
Cleaning Instructions 87
  Cleaning Instructions 87
  Cleaning Instructions for this Instrument 88
  Safety Precautions 88
  Why is it important to clean optical devices?
  What do I need for proper cleaning?
  Preserving Connectors 94
  Cleaning Instrument Housings
  Which Cleaning Procedure should I use?
                                          95
  How to clean connectors 95
  How to clean connector adapters
  How to clean bare fiber adapters
  How to clean optical glass plates 98
  Additional Cleaning Information
  Other Cleaning Hints 99
Firmware Upgrades 99
Error Messages 100
  SYST:ERR? 100
Data Sanitation 100
  Memory Inventory
  Clearing User Data
                     100
```

Index 103



### **Overview**

This chapter helps you start using your N77-Series instrument.

- "Agilent N7744A and N7745A Optical Power Meter" on page 7
- "Agilent N775xA Optical Attenuator and Power Meter and N776xA Optical Attenuator" on page 8
- "Initial Inspection" on page 10
- "Safety Considerations Overview" on page 11
- "AC Line Power Supply Requirements Overview" on page 14
- "Operating and Storage Environment Overview" on page 16
- "Input and Output Connectors Overview" on page 18
- " Making Optical Connections" on page 22
- "Electrical Connectors" on page 25
- "LAN Interface" on page 25
- "GPIB Interface" on page 30
- "USB Interface" on page 31
- "Powering Up the Instrument" on page 32
- "Connecting to the Instrument Overview" on page 33
- "Claims and Repackaging" on page 39

# Agilent N7744A and N7745A Optical Power Meter

The Agilent N7744A and Agilent N7745A are LXI compliant optical power meters with 4 or 8 power-sensor ports respectively. They address productivity challenges by increasing throughput and operational efficiency using industry-leading device connectivity, high-speed



### 1 Getting Started

acquisition of measurement data and fast data transfer for postprocessing.

The Agilent N7744A and Agilent N7745A are designed to characterize optical multiport components, such as multiplexers, PON splitters, wavelength selective switches (WSS) and ROADMs, in all product phases from R&D to manufacturing.

The one rack unit (1U half-rack high) instrument can be operated alone or with other units, by PC connected over LAN, USB, or GPIB.

The unprecedented N7740A fiber connectivity concept consists of a quadruple adapter (quad-adapter) with a quick-locking mechanism. Up to four fiber connectors from the device under test (DUT) can be connected in a comfortable working position to one quad-adapter. The push-on locking mechanism makes it easy and secure to attach the quad-adapter to and remove the quad-adapter from the instrument in seconds. With multiple quad-adapters, the measurement task can be fully separated from the DUT connection task. This means a further increase in test throughput. Where this is not an issue, traditional single connector connectivity is supported as well.

These power meters are highly compatible with setups and programs developed using the Agilent 816x modular power meters. They are also easy to synchronize with Agilent tunable laser sources, to make spectral measurements. The functionality for time-dependent measurement has been extended with sampling at up to 1 MHz and data logging of up to 1 million samples for each port.

# Agilent N775xA Optical Attenuator and Power Meter and N776xA Optical Attenuator

The Agilent N775xA and N776xA are compact multichannel attenuators and power meters. They are a new class of remote-controlled fiber optic instruments. All attenuators feature power control functionality to set the output power level of the attenuator.

The Agilent N7751A and Agilent N7752A are optical attenuators for single mode fiber with one and two single-mode fiber attenuator channels respectively and two extra optical power meter channels, for convenient power measurements and setup calibration. The Agilent N7761A, Agilent N7762A and Agilent N7764A are single-mode fiber optical attenuators for single mode fiber, with one, two and four attenuator channels respectively.

All N775xA and N776xA attenuators address productivity challenges by increasing throughput and operational efficiency.

The Agilent N7766A and the Agilent N7768A are optical attenuators for multimode fibers with two or four attenuator channels repectively. Both multimode attenuators are available with  $50\mu m$  and  $62.5\mu m$  multimode fiber.

The N7766A and N7768A are 2 and 4 channel multimode fiber attenuators that are optimized for quick power setting and preservation of the input modal power distribution.

The attenuators are designed to attenuate and control the power of light in optical fibers. They allow you to accurately and repeatably set the attenuation factor or power level manually using the Viewer software, or using the programming interface.

The attenuators include power control functionality that allow you to set the output power level of the attenuator. When power control mode is enabled, the module automatically corrects for power changes at the input so that the output power level you set is maintained.




Figure 1 Optical Attenuators with Power Control

The one rack unit (1U half-rack high) instrument can be operated alone or with other units, controlled over LAN, USB, or GPIB.

These attenuators and power meters are fully compatible with setups and programs developed using the Agilent 8157x modular attenuators and basic 8163x optical power meter commands, not including logging functionality.

#### **Modal Fidelity for Multimode Fiber Systems**

Signals in multimode fibers are distributed over a range of mode groups that can have different loss and delay in a link. For dependable multimode transceiver testing, the instrument used to set the power level should not change this modal distribution. The bulk-optic filter and collimated beam path of Agilent multimode attenuators are designed to assure homogeneous attenuation of all input modes.

# Agilent N7711/4A Single and 4-Port Tunable Laser System Source

The Agilent N7711A and N7714A tunable lasers are single-port and four-port sources, available with C-band or L-band wavelength coverage. The narrow linewidth and offset grid fine-tuning capability of the N7711A and N7714A make them ideal sources for realistic loading of the latest transmission systems.

To owners of Agilent's proven Lightwave Measurement System the 81950A tunable laser source module offers the same features as the N7711A. The 81950A plugs into the 8163B and 8164B mainframes. For additional information about the 81950A please refer to the Compact Tunable Laser data sheet, publication no. 5988-8518EN.

All models can reach any wavelength point within their specified wavelength range just like all other Agilent tunable lasers. In this mode, code compatibility with existing test setups based on Agilent's range of full-size and compact tunable lasers is a great asset. In system loading applications, it may be preferable to grid-tune the lasers like system transmitters, simply by changing the channel index. The channel grid is adjustable to standard ITU-T grid spacing like 50 GHz, and to arbitrary grids. Likewise, the zero frequency (base channel) of the chosen grid is adjustable. A 12 GHz fine-tuning range allows de-tuning the frequency.

Targeted for high test throughput, lowest cost-per-channel, and narrow footprint, all members of Agilent's 77-Series optical test instruments are built on a common platform and a common PC-based user interface. A complete set of control interfaces including LAN, USB2.0 and GPIB simplifies integration with manufacturing control systems. Code compatible to Agilent's Lightwave Measurement System modules, the new instrument generation can serve as plug-in replacements in existing test solutions. The N77 viewer software allows to control all features of the N7711/4A tunable laser sources.

# **Initial Inspection**

Inspect the shipping container for damage. If there is damage to the container or cushioning, keep them until you have checked the contents of the shipment for completeness and verified the instrument both mechanically and electrically. The Function Tests give a procedure for checking the operation of the instrument. If the contents are incomplete, mechanical damage or defect is apparent, or if an instrument does not pass the operator's checks, notify the nearest Agilent Technologies Sales/Service Office.

NOTE

To avoid hazardous electrical shock, do not perform electrical tests when there are signs of shipping damage to any portion of the outer enclosure (covers, panels, and

# **Safety Considerations**

## **Safety Considerations - Overview**

The following general safety precautions must be observed during all phases of operation, service, and repair of this instrument. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture, and intended use of the instrument. Agilent Technologies Inc. assumes no liability for the customer's failure to comply with these requirements. This product has been designed and tested in accordance with IEC Publication 61010-1, Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory, and has been supplied in a safe condition. The instruction documentation contains information and warnings that must be followed by the user to ensure safe operation and to maintain the product in a safe condition.

# **Safety Symbols**

CAUTION

The caution sign denotes a hazard. It calls attention to a procedure which, if not correctly performed or adhered to, could result in damage to or destruction of the product. Do not proceed beyond a caution sign until the indicated conditions are fully understood and met.

WARNING

The warning sign denotes a hazard. It calls attention to a procedure which, if not correctly performed or adhered to, could result in injury or loss of life. Do not proceed beyond a warning sign until the indicated conditions are fully understood and met.

### 1 Getting Started

### General

This is a Safety Class 1 instrument (provided with a protective earth terminal) and has been manufactured and tested according to international safety standards.

Before operation, you should review the instrument and manual for safety markings and instructions. You must follow these to ensure safe operation and to maintain the instrument in safe condition.

Some circuits are powered whenever the instrument is connected to the AC power source. To disconnect from the line power, disconnect the power cord either at the rear power inlet or at the AC line power source (receptacle). One of these must always be accessible. If the instrument is in a cabinet, it must be disconnected from the line power by the system's line power switch.

WARNING

To avoid hazardous electrical shock, do not perform electrical tests when there are signs of shipping damage to any portion of the outer enclosure (covers, panels, and so on).

### **Operating Environment**

WARNING

The instrument is not designed for outdoor use. To prevent potential fire or shock hazard, do not expose the instrument to rain or other excessive moisture.

## **Line Power Requirements**

CAUTION

The instrument complies with installation category II and can operate from the single-phase AC power source that supplies between 100 V and 240 V at a frequency in the range 50 to 60 Hz. The maximum voltage fluctuation is 10% of the nominal supply voltage. The maximum power consumption is 60 VA with all options installed.

### **Input/Output Signals**

### CAUTION

There is one input BNC connector: Trigger In. This is a TTL input. A maximum of 5 V can be applied as an external voltage to this input connector.

There is one output BNC connector: Trigger Out. This is a TTL output. Do not apply an external voltage to this connector.

### **Line Power Connectors**

In accordance with international safety standards, the instrument has a three-wire power cable. When connected to an appropriate AC power receptacle, this cable earths the instrument cabinet. The type of power cable shipped with each instrument depends on the country of destination. Please refer to "Accessories - Overview" on page 83" for the part numbers of available power cables.


### WARNING

To avoid the possibility of injury or death, you must observe the following precautions before switching on the instrument.

- Insert the power cable plug only into a socket outlet provided with a protective earth contact. Do not negate this protective action by the using an extension cord without a protective conductor.
- Do not interrupt the protective earth connection intentionally.
- Do not remove protective covers. Operating personnel must not remove instrument covers. Component replacement and internal adjustments must be made only by qualified service personnel.
- Instruments that appear damaged or defective should be made inoperative and secured against unintended operation until they can be repaired by qualified service personnel.
- Defective, damaged, or malfunctioning laser sources must be returned to an Agilent Technologies Service Center.

Do not operate the instrument in the presence of flammable gases or fumes. Operation of any electrical instrument in such an environment constitutes a definite safety hazard.

# **Instrument Markings**



# **AC Line Power Supply Requirements**

# **AC Line Power Supply Requirements - Overview**

This secton provides information on:

- Line Power Requirements
- Line Power Cable
- · Changing the Fuse

### **Line Power Requirements**

The instrument complies with installation category II and can operate from the single-phase AC power source that supplies between 100 V and 240 V at a frequency in the range 50 to 60 Hz. The maximum voltage fluctuation is 10% of the nominal supply voltage. The maximum power consumption is 60 VA with all options installed.

### Line Power Cable

In accordance with international safety standards, the instrument has a three-wire power cable. When connected to an appropriate AC power receptacle, this cable earths the instrument cabinet. For the part number of the power cable for your country and instrument see the "Accessories - Overview" on page 83 chapter.

CAUTION

Please note that the switch on the front panel of the instrument does not stop the flow of power to the instrument.

If you need to turn off the power, unplug the instrument at the mains or remove the power cable connector from the appliance coupler at the rear of the instrument. For this reason, the power cable connection should be easily accessible - allowing you to turn off the power quickly. If the instrument is in a cabinet, it must be disconnected from the line power by the system's line power switch.

The power switch allows you to switch between stand-by mode and power-on mode.




Figure 2 Power Switch

When the instrument is in stand-by mode, the Pwr LED is orange. When the instrument is powered-on, the Pwr LED is green.

### WARNING

To avoid the possibility of injury or death, you must observe the following precautions before switching on the instrument.

- Insert the power cable plug only into a socket outlet provided with a protective earth contact. Do not negate this protective action by using an extension cord without a protective conductor.
- Do not interrupt the protective earth connection intentionally.

The AC power requirements are summarized on the rear panel of the instrument.




Figure 3 AC Power Requirement Markings

# **Changing the Fuse**

CAUTION

There is no user replaceable fuse for the instrument. Changing the fuse should be carried out only by Agilent Technologies service personnel. If you need to get the fuse replaced refer to your nearest Agilent Technologies Sales/Service Office.

# **Operating and Storage Environment**

# **Operating and Storage Environment - Overview**

The following summarizes the operating environment ranges. In order for the instrument to meet specifications, the operating environment must be within these limits. WARNING

The instrument is not designed for outdoor use. To prevent potential fire or shock hazard, do not expose the instrument to rain or other excessive moisture.

## **Temperature**

The instrument should be protected from temperature extremes and changes in temperature that may cause condensation within it.

The operating temperature is from  $5^{\circ}$ C to  $+40^{\circ}$ C.

The storage temperature is from -40°C to +70°C.

# **Humidity**

The operating humidity is 15 to 95%, non-condensing.

### **Altitude**

The maximum operating altitude is 2000 m.

### **Pollution Protection**

The instruments are designed for pollution degree 2.

# **Instrument Cooling**

The instrument has a cooling fan mounted internally.

Mount or position your instrument upright and horizontally, as shown in Figure 4 so that air can circulate through it freely.

### **Operating Position**

When operating the instrument choose a location that provides at least 75 mm (3 inches) of clearance at the rear, and at least 25 mm (1 inch)

### 1 Getting Started

of clearance at each side. Failure to provide adequate air clearance may result in excessive internal temperature, reducing instrument reliability. The instrument should not be operated when resting on its rear or side panels.




Figure 4 Correct Operating Position

(shown here for the Agilent N7744/5A, valid also for the N7751/2A and Agilent N7761/2/4A, N7766/68A and the N7711/4A)


# **Input and Output Connectors**

# **Input and Output Connectors - Overview**


This section provides information on:

- Front Panels Multiport Power Meters
- Front Panels Optical Attenuators
- Front Panels Tunable Laser Source
- Rear Panel

# **Front Panels - Multiport Power Meters**



Front panel of the Agilent N7744A 4-Channel Optical Multiport Power Meter



Front panel of the Agilent N7745A 8-Channel Optical Multiport Power Meter

# **Front Panels - Optical Attenuators**





Figure 7 Front panel of the Agilent N7751A 1-Channel Optical Attenuator and 2-**Channel Power Meter** 



**Figure 8** Front panel of the Agilent N7752A 2-Channel Optical Attenuator and 2-Channel Power Meter






Figure 9 Front panel of the Agilent N7761A 1-Channel Optical Attenuator



**Figure 10** Front panel of the Agilent N7762A and N7766A 2-Channel Optical Attenuator



**Figure 11** Front panel of the Agilent N7764A and N7768A 4-Channel Optical Attenuator

### **Front Panels - Tunable Laser Source**




Figure 12 Front panel of the Agilent N7711A Single-Port Tunable Laser System Source




Figure 13 Front panel of the Agilent N7714A Single-Port Tunable Laser System Source

### **Rear Panel**




Figure 14 Rear panel of the Agilent N7744/5A Optical Multiport Power Meter, Agilent N7751/2A Optical Attenuator and Power Meter, and Agilent N7761/2/4/6/8A Optical Attenuator and Agilent N7711/4A Tunable Laser System Source

# **Making Optical Connections**

### **Making Optical Connections**

The optical attenuators are designed for use with FC connectors.

Attenuators ordered with option 021 have straight connector interfaces.

Attenuators ordered with option 022 have angled connector interfaces.

The tunable lasers are designed for use with FC connectors.

Lasers ordered with option 071 have straight connector interfaces.

Lasers ordered with option 072 have agled connector interfaces.

The N7751/2A power meters also use FC connectors but are non-contact. They can be used with both angled and straight connectors.

The N7766A and N7768A multimode attenuators have straight connector interfaces. The fiber option determines the dimensions of the internal fiber core, such as 50 or 62.5 $\mu$ m. Best performance is achieved by matching this to the fibers used in the setup.

CAUTION

If the connectors on your attenuator or lasers are angled, you can only use cables with angled connectors with the instrument.

If the connectors on your attenuator are straight, you can only use cables with straight connectors with the instrument.



Figure 15 Angled and Straight Connector Symbols

The above figure shows the symbols that tell you whether the optical connectors of your instrument are angled or straight. The angled contact connector symbol is colored green.

You should connect straight contact fiber end connectors with neutral sleeves to straight contact connectors, or connect angled contact fiber end connectors with green sleeves to angled contact connectors.

NOTE

You cannot connect angled non-contact fiber end connectors with orange sleeves directly to the attenuator or laser.

See "Accessories - Overview" on page 83 for further details on connector interfaces and accessories.

# Making Quad-adapter Connections on the Agilent N7744/45A Multiport **Power Meter**




Figure 16 Connecting Optical Fibers

NOTE

Before you attach a quad-adapter or fiber, you should make sure they are clean. This is important for accurate measurements.

The fibers do not contact the optical inputs (all contact with the optical inputs should be avoided). Usually the only cleaning necessary for the optical inputs is the gentle removal of dust with pressurized, dry air.

For further information see " Maintenance and Troubleshooting - Overview " on page 87.

To use a quad-adapter to connect optical fibers to the power meter inputs:

1 Select the quad-adapter that suits the optical connector that your optical fiber is terminated with.

### 1 Getting Started

- **2** Attach the optical connectors to the quad-adapter. You can attach up to 4 fibers to the quad-adapter at the same time.
- **3** Make sure the number labels on the quad-adapter are at the top.
- **4** Connect the quad-adapter to the optical inputs by pressing it on until it clicks securely into place.

# Removing Quad-adapter Interfaces from the Agilent N7744/5A Multiport Power Meter

NOTE

There are release levers on both sides of the quad-adapter for removing it from the optical inputs.




Figure 17 Removing Quad-adapter Interfaces

- 1 Press the release levers on both sides of the quad-adapter to release it from the front panel of the power meter.
- 2 You can now remove the optical connectors from the quad-adapter.

NOTE

Always place the dust cap on unused connector interfaces.

You can use the quad-adapter dust caps to protect the four power meter inputs and the rear of the quad-adapter, for example, when fibers are attached.

# **Optical Output**

### **Polarization Maintaining Fiber**

A Polarization maintaining fiber (PMF) output is standard for Agilent N7711/4A tunable laser system sources. PMF is aligned to maintain the state of polarization. A well defined state of polarization helps ensure

constant measurement conditions. The fiber is of Panda type, with TE mode in the slow axis in line with the connector key.




Figure 18 Polarization maintaining fiber

# **Electrical Connectors**

There are two BNC connectors on the rear panel of your instrument. These are the Trigger Out and the Trigger In connectors.

CAUTION

The Trigger In is a TTL input. A maximum of 5 V can be applied as an external voltage to this input connector.

The Trigger Output is the only output BNC connector. This is a TTL output. Do not apply an external voltage to this connector.

# **LAN Interface**

### **LAN Interface**

This section explains the concept of LAN in details.

### **Selecting a LAN Network**

For the purposes of this guide, a private (isolated) LAN network is defined as a network configuration in which instrument access is a direct connection between the computer and the instrument, or to multiple instruments connected via a dedicated router or switch. A site (company-wide) LAN is defined as a network in which instrument access is available to many users at on-site and remote locations. The instrument's application and/or your company's IT (Information Technology) department may have guidelines that help decide the type of network (private or site) used. If a network configuration has not been determined, refer to the following considerations concerning each type.

### **Private LAN Considerations**

Among the basic parameters of a private LAN network to consider are security, performance, reliability, and IP address availability.

## Security

A private network generally involves a direct connection between the computer and the instrument, or to multiple instruments using switches or routers. Access to the instrument is limited to users on the private network, as opposed to users on a site network that could locate and access the instrument from any location.

Private networks can reduce the possibility of tests being disrupted by unplanned or unauthorized access. Code generation for test systems on a private network is often simplified, as provisions against unauthorized users may not be required.

#### **Performance**

Test systems where large amounts of data are transferred usually have faster throughput on a private network. On a site network, heavy and unpredictable LAN traffic affects each instrument (node) on the network. The impact on a test system is that repeatability is difficult to achieve as latencies are difficult to account for.

### Reliability

Private networks are fundamentally more reliable than site networks as they host fewer users and are less complex than site networks.

Private networks are isolated from conditions that could bring down (crash) a site network.

# **IP Address Availability**

Every instrument (node) on a LAN (private or site) has an IP (Internet Protocol) address. Due to the expanding use of the internet, the number of site network IP addresses available is limited. By using a router with Dynamic Host Configuration Protocol (DHCP) capability on a private network, the router can assign an IP address to each instrument thus creating a sub-network (subnet) that does not consume site IP addresses.

### **Site LAN Considerations**

For applications requiring access by many users or by users at distributed sites, a site LAN network is required. In addition to supporting multiple users, site LANs often offer the advantage of being maintained by IT departments. When using a site LAN, consult your IT department regarding LAN configuration and security issues.

# **Connecting the LAN Cable**

LAN cables are connected between the LAN terminal on the instrument and the computer, or between the instrument and a router or switch if included in your network.

#### **Private Network Connections**

The figure below shows typical LAN cable connections for a private network.




Figure 19 Typical Private (isolated) LAN Network Connections

If the instrument is connected directly to the PC, use a LAN crossover cable, as provided with the instrument. If your computer supports Auto-MDIX or contains a LAN card with gigabit data transfer rates, the (blue) crossover cable is not required. A standard LAN cable can be used instead.

For private LAN networks that include a switch or router, use standard LAN cables for all network connections. Do not use a crossover cable.

### **Site Network Connections**

The following figure shows a typical LAN cable connection for a site network.




Figure 20 Typical Site LAN Network Connections.

On site networks, the instrument and the computer are connected directly to site LAN ports, or are connected to the site LAN through a switch. In each site network configuration, standard LAN cables are used.

#### The LAN LED

As the LAN connection is made, the DHCP server assigns an address and the LXI device identification proceeds you will see the following indicators:,




Figure 21 The LAN LED

There may be a delay between making the LAN link (yellow status) and getting the IP address (green status). This delay may be longer if there is no DHCP server, for example when the instrument is connected directly to a PC.

### The LAN Reset button

This recessed button has two functions.

Pressing the button briefly invokes a preset of the instrument and restores default measurement settings. This is equivalent to the programming command

:SYSTem:PRESet

Pressing and holding the button for 3 seconds will reset the LAN parameters to the factory default. This includes changing the password for modifying the configuration via LAN back to the default "agilent". This is equivalent to the following sequence of programming commands:

:SYSTem:COMMunicate:ETHernet:RESet

# **GPIB** Interface

### **GPIB** Interface

You can connect your GPIB interface into a star network, a linear network, or a combination star and linear network. The limitations imposed on this network are as follows:

- The total cable length cannot exceed 20 meters.
- The maximum cable length between devices is 4 meters, with an average separation of 2m between devices over the whole bus.
- No more than 15 devices may be interconnected on one bus.

# **Cables and Adapters**

See "Cables" on page 84 for details on cables and adapters.

### **Connector**

The following figure shows the connector and pin assignments.

Connector Part Number: 1251-0293




Figure 22 GPIB Connector

### CAUTION

Agilent Technologies products delivered now are equipped with connectors having ISO metric-threaded lock screws and stud mounts (ISO M3.5×0.6) that are black in color. Earlier connectors may have lock screws and stud mounts with Englishthreaded lock screws and stud mounts (6-32 UNC) that have a shiny nickel finish.

### CAUTION

It is recommended that you do not stack more than three connectors, one on top of the other.

Hand-tighten the connector lock screws. Do not use a screwdriver.

### **GPIB Logic Levels**

The instrument's GPIB lines use standard TTL logic, as follows:

True = Low = digital ground or 0 Vdc to 0.4 Vdc

False = High = open or 2.5 Vdc to 5 Vdc

All GPIB lines have LOW assertion states. High states are held at 3.0 Vdc by pull-ups within the instrument. When a line functions as an input, it requires approximately 3.2 mA to pull it low through a closure to digital ground. When a line functions as an output, it will sink up to 48 mA in the low state and approximately 0.6 mA in the high state.

NOTE

The GPIB line screens are not isolated from ground.

### **USB** Interface

NOTE

Before connecting the instrument over USB, install the Agilent I/O Libraries Suite, which is included on a CD with the instrument.

The instrument is a USB device, with a mini-USB connector.

# **Powering Up the Instrument**

# **Powering Up the Instrument**

When you switch on the instrument, the LEDs on the front panel show the various stages of booting.




Figure 23 Powering Up the Instrument

### **Front Panel Indicators**

During operation, you may notice the following LED indicators:




Figure 24 Front Panel Indicators

# **Connecting to the Instrument**

### **Connecting to the Instrument - Overview**

When you are first connecting the instrument, we recommend you use the Agilent Connection Expert included with the Agilent I/O libraries (available on the CD supplied with the instrument, or from www.agilent.com).

# **Connecting over USB**

NOTE

The screenshots show the Agilent N7745A Multiport Power Meter. The same procedure applies when connecting an Agilent N7744A Multiport Power Meter, an Agilent N7751/2A Optical Attenuator and Power Meter, or an Agilent N7761/2/4/6/8A Optical Attenuator or an Agilent N7711/4A Tunable Laser System Source.

1 If it is not already running, start the Agilent Connection Expert software.




Figure 25 Agilent Connection Expert software

2 Attach the instrument to the USB. The first time you attach the instrument, follow the instructions on the screen to select the driver for this instrument.

🏖 Refresh All 🌎 🥱 Undo The USB interface(s) connected directly to this computer Refresh this USB interface This item has been verified Change Properties... Change the label
Ignore
Delete Less << VISA interface ID: USBO SICL interface ID: Logical unit: Auto-discover General Tasks More Information How do I get drivers Where can I find programming sample:

When the instrument is connected, it is shown in the list.

Figure 26 Select the instrument in the list

**3** Select the instrument in the list to see the associated tasks (in the Task Guide on the left) and information (on the right).

# Finding the IP Address of an instrument

- 1 Connect to the instrument over USB, as described above.
- **2** Select the instrument in the list.




Figure 27 Select the instrument

Agilent VISA: The primary VISA library

**3** Click on "Send commands to this instrument". The default command is the \*IDN? query.





Figure 28 \*IDN? query

4 Enter the command

:SYSTem:COMMunicate:ETHernet:IPADdress:CURRent?

then click on [Send & Read].

The instrument returns its current IP address.



:SYSTem:COMMunicate:ETHernet:IPADdress:CURRent? Figure 29

# **Connecting over LAN**

NOTE

The screenshots show the Agilent N7745A Multiport Power Meter. The same procedure applies when connecting to other N77xx instruments like attenuators and laser sources.

- 1 Make sure the instrument is connected to the LAN, and that the LAN LED on the front panel is green. See, "LAN Interface" on page 25 for more information to connecting over LAN.
- **2** If it is not already running, start the Agilent Connection Expert software.
- **3** Select "LAN (TCPIP0)" in the list.




Figure 30

- 4 Click on "Add an instrument" in the Task Guide.
- **5** Wait for the scan to finish.
- **6** Select the instrument in the list.




Figure 31

7 To show the description of the N775xA or N776xA, check the box for \*IDN? query and click on "Find Again".




Figure 32

- **8** Select the instrument in the list to see the associated tasks (in the Task Guide on the left) and information (on the right).
- For the N774x Multiport power meters, click on the "Instrument Web Interface" button to control the instrument directly.




Figure 33

If it is not already installed, you will need to install a Java run-time engine to use this page.

Basic configuration and functionality are available from this web page.

- For the attenuators, and as an alternative for the power meters, you can control the instrument from the Agilent N77xx Viewer.
  - Connecting using the N77xx Viewer is described in .
  - Making power measurements with the N77xx Viewer is described in .
  - Controlling attenuation with the N77xx Viewer is described in "Controlling the Optical Attenuator" on page 45.

More advanced functionality is available through programming with SCPI commands (described in the Programming Guide available on the CD supplied with the instrument or from www.agilent.com), or using the 816x VXI Plug&Play driver, version 4.3 or higher. For the N774x Multiport power meters you can also use the IVI-COM driver.

# **Claims and Repackaging**

# **Claims and Repackaging**

If physical damage is evident or if the instrument does not meet specification when received, notify the carrier and the nearest Agilent Technologies Sales/Service Office. The Agilent Technologies Sales/ Service Office will arrange for repair or replacement of the unit without waiting for settlement of the claim against the carrier.

# **Return Shipments to Agilent Technologies**

If the instrument is to be shipped to an Agilent Technologies Sales/ Service Office, attach a tag showing owner, return address, model number and full serial number and the type of service required.

The original shipping carton and packing material may be reusable, but the Agilent Technologies Sales/Service Office will provide information and recommendations on materials to be used if the original packing is no longer available or reusable. General instructions for repackaging are as follows:

- Wrap instrument in heavy paper or plastic.
- Use strong shipping container. A double wall carton made of 350pound test material is adequate.
- Use enough shock absorbing material (3 to 4 inch layer) around all sides of the instrument to provide a firm cushion and prevent movement inside container. Protect control panel with cardboard.
- Seal shipping container securely.
- Mark shipping container FRAGILE to encourage careful handling.
- In any correspondence, refer to instrument by model number and serial number.

#### **Deleting user data**

If you need to delete all your logged data and user configurations, that is to reset the instrument completely:

- 1 Press and hold the Reset/standby button until the instrument starts to reboot.
- **2** When the instrument reboots, press the Reset/standby button once more briefly, to confirm you want to delete all the stored data.

# N77xx Viewer Common User Interface Functions

# **Connecting to an Instrument**

When you first launch the N77xx Viewer, it is not connected to any of your instruments.

- 1 Click on the "Connect" icon at the top right of the window.
- **2** Select the instrument to which you want to connect from the list.




Figure 34 Connecting to an Instrument

- 3 The list shows instruments that have been given a VISA address. If your instrument is not in the list,
  - 1 Use the Agilent Connection Expert to check your instrument has a VISA address
  - 2 Click on the "Scan Again" button. Once communication has been established to the instrument, its details are shown.
- 4 If you need to change the connection configuration for an instrument on the LAN, click on the "Configure Interfaces" button.

You can only configure the LAN parameters before the instrument is connected.




Figure 35 **Interface Parameters** 

- Make sure both DHCP and Auto IP are not selected.
- You can now edit the IP parameters.
- **5** For instruments you use regularly, you can add a link on your desktop by putting a check mark in the box to "Place a shortcut on the Desktop".
- **6** Click on the "OK" button to finish the connection.

# **Viewing Measurements**

There are two ways you can view results for a power meter or attenuator.

# Viewing the measurement on a single channel

1 Click on the "New" button to open a window with the current measurement and parameter settings.



**Figure 36** Viewing the measurement on a single channel attenuator/power meter (left) or tunable laser (right)

# Viewing all measurements

1 Click on Overview to open a window with the measurements on all channels.




Figure 37 Viewing all measurements

# **Saving and Recalling Configurations**

# Saving a configuration

You can save the current measurement configuration to either of two settings stored on the instrument.

- 1 In the File menu, select Save.
- 2 Enter or select the setting to which you want to save it.




Figure 38 Saving a configuration

**3** Click on the "OK" button.

# **Recalling a configuration**

- 1 In the File menu, select Load.
- 2 Enter or select the setting you want to retrieve.
- **3** Click on the "OK" button.

# **Controlling the Instrument Refresh**

1 In the File menu, select Update.
In the dialog, you can:





Figure 39 Controlling the Instrument Refresh

- Stop and start the display of all changes. Make sure to put a check in the box for "Update enabled", to enable the N77xx Viewer to display the status and results from the instrument. If there is no check in the box for "Update enabled", the display is not updated with the replies to queries sent to the instrument.
- Set the update interval for "Readings" to control how often the result is read from the instrument and the display updated.
- Set the update interval for "Polled Values" to control how often
  the measurement parameters are read from the instrument and
  the display updated. This is useful if, for example, more than one
  person can control the instrument.
- **2** When you have set all the values, click on the "OK" button.



# **Attenuation**

This chapter describes how to use the Agilent N7751A and N7752A Optical Attenuator and Power Meter, the Agilent N7761A, N7762A, and N7764A Single Mode Fiber Optical Attenuator and the Agilent N7766A and N7768A multimode fiber Optical Attenuators to attenuate and control optical power.

If you want to control the power meter channels, please refer to the previous chapter .

- "Controlling the Optical Attenuator" on page 45
- "Setting the Attenuation" on page 48

**Controlling the Optical Attenuator** 

There are different possibilities to control the instrument:

- LAN / USB / GPIB Connection: [N77xx Viewer]
- LAN / USB / GPIB: Programming with SCPI commands (described in the Programming Guide available on the CD supplied with the instrument or from www.agilent.com), or associated drivers, including the 816x VXI Plug&Play driver.

To connect to the instrument see "Input and Output Connectors - Overview" on page 18.

The most convenient way to control one or more Attenuators is using the Agilent N77xx Viewer, which is a control and display graphical user interface.

#### 2 Attenuation




Figure 40 The Viewer

The Viewer is available to download from the product webpage which you can access at

http://www.agilent.com/find/voa

Then select the tab "Technical Support" followed by "Drivers & Software".




Figure 41 Technical Support N7744A

On this webpage you will also find additional software tools like:

- Firmware update utility
- · Firmware update file
- VXI Plug&Play Drivers

# **Setting Attenuation**

# **Setting the Number of Digits displayed**

When the power is displayed with the Viewer program, you can set the maximum number of digits that are used in power measurement. This is the maximum number of digits after the decimal point.

# **Selecting the Power Units**

The Power Unit setting allows you to select either Watts (W) or dBm as the units in which the output power is displayed or controlled. Attenuation or offset values are in dB, and are the power relative to the input power.

# **Setting the Attenuation**

Before applying an attenuation factor, set the Attenuator module to the appropriate wavelength.

The instrument uses a filter technology with a low wavelength dependence. Nevertheless, for the best possible accuracy:

- If the source is of a particular wavelength, set the Attenuator module to that wavelength.
- If the source is a multi-wavelength DWDM signal, set the Attenuator module to the weighted wavelength of the spectrum.

You can set the attenuation for each channel. The power reading represents the average optical power over the selected averaging time. For optimal output power control the power meter continuously self-triggers new measurements.

# **Applying an Attenuation Offset**

Attenuating an input signal to a device under test (DUT) introduces additional connectors and other components into your optical test setup. Use the attenuation offset to compensate for the insertion loss of these components.

1 Measure the insertion loss of the attenuator and the connecting cables using a reference power meter.

If you are using an N775xA attenuator and power meter, you can use the internal power meter.

- Using all the components and connectors, without a DUT connected, connect the output of the attenuator to the input of the power meter.
- Calculate the difference between the power at the attenuator output, and the power measured by the power meter.
- **2** Enter the result into an <sub>Offset</sub> parameter.

The attenuator then compensates for this value.

## **Setting the Filter Transition Speed**

By default, the attenuation transition speed to a new setting is as fast as possible. If your test setup requires a slower attenuation transition speed, set a <Speed> between 0.1 and 1000 dB/s.

The N7766A and N7768A multimode attenuators accept controlled speeds between 0.1 and 80 dB/s and a full speed of about 1000 dB/s. Settings above 80 dB/s will result in full-speed operation.

# **How to Remove Electrical Offsets (Zero)**

The optical power meters in the attenuator measure optical power by converting optical power to electrical current, and then measuring electrical current. An electrical offset is electrical current that is always present, even if there is no optical power input. If electrical offsets are not removed, they affect the accuracy of attenuation and power control, especially at low optical power levels.

Performing a zero sets the zero power level to the average electrical offset level for the current environmental conditions.

NOTE

The environmental conditions and the temperature of the instrument affect electrical offset.

For the best results you must:

- Allow the instrument time to acclimatize (around 24 hours).
- Allow the instrument time to warm up (around 20 minutes).
- Make sure that the optical input and output is not receiving any light.
   If you are using multi-mode fiber-optic cable, you must disconnect the cable and cover the input and the output of the attenuator to perform a zero.

It is good practice to perform a zero before making any sensitive measurements.

NOTE

While zeroing an attenuator's power monitor or an integrated power meter, all attenuators in the unit are stopped and all shutters are closed to ensure the best possible zeroing performance.

1 In the File menu, select Zeroing.

#### 2 Attenuation




Figure 42 How to Remove Electrical Offsets (Zero)

2 Select the channel you want to zero from the drop down list and click on the "This" button.

OR

Click on the "All" button.

While the zero function is running, the traffic light flashes amber.

If you receive a zero error, the zeroing operation has failed because the power meter received light. The most common reason for zeroing to fail is if:

- a source is connected to the attenuator's input or output connector,
- the fiber connected to the attenuator's input or output connector is collecting light, or
- the attenuator receives ambient light because the input or output connector is uncovered.

When zeroing is finished, the traffic light is green. You can now proceed with your measurements.

When environmental conditions are stable, the zero setting will remain valid over long times and the same settings will be active after rebooting the instrument.

# **Enabling and Disabling the Output**

You can turn on and off the attenuator output both for protection purposes or to simulate channel drops.

Select or deselect the <Output>.

- When the "Output" box is checked, the output is on.
- When the "Output" box is not checked, the attenuation is set to the optical path blocking attenuation (>40 dB).
  - To turn on all the attenuator outputs, click on [All on].
  - To turn off all the attenuator outputs, click on [All off].

# **How to Set a Power Level**

Power Control automatically corrects for power changes at the attenuator input to maintain your desired output power. This feature is designed to compensate, for example, for unknown or drifting input power.

When Power Control is On:

- 1 After entering a new Pset, the attenuator is driven with the selected attenuation transition speed until the output power has reached Pset.
- **2** After settling, the power control uses the attenuation transition speed of 1000dB/s to achieve best input power correction.
- **3** Input power is measured using your set Averaging Time.
- 4 The attenuator can only maintain your desired  $P_{\rm set}$  for attenuations down to 0 dB, and within the specified attenuation and power setting range.
  - **a** Set the Wavelength to the wavelength of the source.
  - **b** Edit  $P_{set}$  to set the power output of the attenuator directly.

# **Applying a Power Offset**

Attenuation in your optical test setup can be compensated so the displayed output power value accurately reflects the power applied to the DUT.

NOTE

Before taking a power measurement, make sure the source is warmed up so that its power output is stable.

- 1 Measure the insertion loss of the components between the attenuator output and the DUT, using a reference power meter.
- **2** Enter the result into the P Offset parameter.

# How to Set the Wavelength

This is the wavelength value of the signal to be measured. The responsivity of the power meter varies with wavelength.

# **How to Set the Averaging Time**

This is the length of time over which a signal is averaged to produce a measured value. Longer averaging times increase the stability and improve the noise rejection.

The averaging time of the power meter and power monitor in the attenuator can be set over a wide range from 2ms to 1s.

Please note, even if you set a longer averaging time, internally the power control always works with an average time of 2ms to achieve best input power correction.

# **Compensating for Wavelength Dependencies in your test setup**

# **Applying an Offset Table**

For the  $P_{Offset}(l)$  parameter select <On> to enable or <Off> to disable the Wavelength Offset Table.

When  $P_{Offset}(l)$  is enabled, the offset applied depends on the operating wavelength you set for the attenuator.




Figure 43 Extrapolation and Interpolation of Offset Values

- If the operating wavelength matches a wavelength stored in the Offset Table, the exact offset value is used.
- If the operating wavelength is between two wavelengths stored in the Offset Table, the offset applied is calculated by linear interpolation.
- If the operating wavelength is greater than, or less than, the range of wavelengths stored in the Offset Table, the offset value for the closest wavelength is extrapolated.

#### **Creating or Editing the Offset Table**

You can create or edit the Offset Table if the P<sub>Offset</sub>(l) function is not on.

- **1** Make sure the the P<sub>Offset</sub>(1) function is NOT on.
- **2** Click the "Offset Table" button.
- **3** If you want to edit an existing offset table, click the "Load" button to load it from your PC, or click on the "Get from N77xx" button, to load the table being used by the attenuator.
  - To add a new pair, enter the wavelength:offset pair in the empty row at the bottom of the table.
  - The Offset Table is automatically ordered by wavelength; there is no need to enter wavelength: offset pairs in any particular order.
  - To edit the offset for an existing entry, click in the field and enter the new offset.
  - To remove an existing entry, click in the wavelength field and delete the wavelength.
- **4** To save a copy of the offset table to your PC, click the "Save" button.

# 2 Attenuation

To apply the offset table, click on "Set to N77xx" button, to save it to the attenuator for use.

**5** Click on the "OK" button.

# **Definition of Terms**

# (Applicable) fiber type

Connection type for which the specifications and characteristics apply (if not differently stated).

# Absolute wavelength (frequency) accuracy

The maximum difference between the displayed wavelength (frequency) and the actual wavelength (frequency) of the tunable laser source. Wavelength is defined as wavelength in vacuum.

#### **Attenuation**

Difference (in dB) between "Total loss" (with any attenuation setting, see "Total loss" on page 72) and "Insertion loss" (i.e. with attenuation set to zero, see "Insertion loss" on page 59) of the attenuator.




Figure 44 Attenuation



# Attenuation setting mode

Operating mode where the user sets the desired attenuation. The power control function is deactivated in this mode (see also "Power setting mode" on page 66).

# Attenuation accuracy (uncertainty)

The maximum possible difference (in dB) between the displayed "Attenuation" (see, "Attenuation" on page 55) and actual "Attenuation" (see, "Attenuation" on page 55).




Figure 45 Attenuation accuracy (uncertainty)

**Conditions:** Attenuator set to the wavelength of the source. Input polarization, temperature range, maximum power and input mode (multimode attenuators only) as specified. Other conditions as specified.

**Measurement:** With Fabry-Pérot laser source of applicable wavelength, polarization scrambler (if applicable), and optical power meter.

NOTE

When using a polarized source, the attenuation uncertainty increases by  $\pm$  half the specified polarization dependent loss of the attenuator (if specified).

# Averaging time

Specifies the period during which the power meter or monitor power meter takes readings from the photodetector. The displayed power result is the (arithmetical) average of these readings.

# **Constant operating conditions**

This generally includes constant values of temperature, humidity, wavelength, input power level, polarization state and mode distribution, if the quantity is not explicitly subject to variation.

# **Data logging capability**

Maximum number of power readings per port that can be captured during a logging operation.

# Drift (dark)

Specifies the maximum deviation from zero of the measured power with no optical input (dark) over time, excluding "Noise pp (dark)" on page 60, power expressed in Watts.

**Measurement:** From the dark measurement result  $P_i$ , the 3rd order polynomial fit  $F_i$  is calculated. The fit incorporates the least square-sum difference. The Drift D is calculated by

 $D = max_i\{F_i\}$  or  $mini\{F_i\}$ , whichever is worse.




Figure 46 Drift (dark)

**Conditions:** All ports dark. Observation time as specified. Other conditions as specified.

#### Dynamic range (logging mode)

For a fixed power meter range, dynamic range D specifies the difference between the maximum displayable value  $P_{max}$  in dBm in the corresponding power meter range and the minimum measurable power  $P_{min}$  in dBm.

Conditions: As specified.

**Measurement:**  $P_{max}$  is 3 dB above the nominal power meter range value (that is, in the +0 dBm power meter range  $P_{max}$  is +3 dBm), except for the +10 dBm range with  $P_{max}$  = +10 dBm.  $P_{min}$ , dBm is calculated as the root-sum-square (RSS) of the specified "Noise 2 sigma (dark)" (see, "Noise 2 sigma (dark)" on page 61 and "Drift (dark)" (see, "Drift (dark)" on page 57. RSS is calculated in linear units (Watts).

# **External analog modulation - modulation depth**

Specifies half the peak-to-peak optical power change divided by the average optical power for a sinusoidal input voltage at the analog modulation input. The average power is defined as half the sum of maximum and minimum power.

**Conditions:** Modulation input signal as specified, modulation frequency as specified.

**Measurement:** Using a photoreceiver (of sufficient bandwidth) and an oscilloscope.

# **Frequency Response**

Specifies the frequency  $f_{\rm resp}$  up to which the power meter's transfer function is above a defined cutoff level.

The transfer function is the measured optical amplitude on a modulated sinusoidal optical input divided by the actual amplitude over modulation frequency, expressed in dB.




Figure 47 Frequency Response

Conditions: Cutoff level as specified. Other conditions as specified.

#### **Insertion loss**

Specifies the "Total loss" at a displayed attenuation of 0 dB

**Conditions:** Temperature range, wavelength range, source polarization input polarization, input mode (multimode attenuators only) as specified.

**Measurement:** With Fabry-Pérot Laser Source, polarization controller (if applicable) and power meter, using high quality connectors in perfect optical condition.

NOTE

"Insertion loss excluding connectors" means connector loss is not part of the loss value.

"Insertion loss including connectors" includes the additional insertion loss from the connectors.

# Linearity

Specifies the relative difference between the measured (displayed) power ratio  $D_x/D_0$  and the actual (true) power ratio  $P_x/P_0$ , caused by changing the power level from the reference level  $P_0$  to an arbitrary level  $P_x$ .

Powers expressed in Watts.

Excluding "Noise pp (dark)" (see, "Noise pp (dark)" on page 60 and "Drift (dark)" (see, "Drift (dark)" on page 57).

$$L_{\%} = \left(\frac{D_x/D_0}{P_x/P_0} - 1\right)$$
 when expressed in %.

$$L_{\rm dB}$$
 =  $10 \log \left( \frac{D_x/D_0}{P_x/P_0} \right)$  when expressed in dB.

Conditions: Reference power level  $P_0$  = 0.01 mW (-20 dBm). Other conditions as specified.

#### Linewidth

The 3 dB width of the optical spectrum, expressed in Hertz.

**Conditions:** SBS suppression off.

#### **Measurement:**

*Using a self-heterodyning technique:* The output of the laser under test is sent through a Mach-Zehnder interferometer in which the length

difference of the two arms is longer than the coherence length of the laser. The electrical noise spectrum of the photodetector current is measured with Agilent lightwave signal analyzer, and the linewidth is calculated from the heterodyne spectrum.

Alternatively, *Using a heterodyning technique*: The output of the laser under test is mixed with another laser of the same type on a wide bandwidth photodetector. The electrical noise spectrum of the photodetector current is measured with Agilent Lightwave signal analyzer, and the linewidth is calculated from the heterodyne spectrum. (Lightwave signal analyzer settings: resolution bandwidth 1 MHz, video bandwidth 10 kHz, sweep time 20 ms, single scan).

# Maximum safe (input) power

The maximum input power that can be applied to any port of the power meter or attenuator without permanent change to its characteristics.

**Attention:** Applying more than the specified power might damage the power meter or attenuator!

NOTE

For input powers > 20 dBm, clean connectors in good condition are vital to avoid thermally induced fiber damage.

#### Minimum averaging time

The minimum selectable averaging time.

# **Maximum output power**

The maximum achievable output power of the tunable laser source and the maximum output power for which the tunable laser source specifications apply.

**Conditions:** As specified.

**Measurement:** Using a power meter at the output of the instrument.

# Noise pp (dark)

Specifies the peak-to-peak (pp) change of displayed power over time with zero input power (dark), excluding "Drift (dark)" (see, "Drift (dark)" on page 57). Powers expressed in Watts.

**Conditions:** All ports dark. Averaging time as specified. Other conditions as specified.

**Measurement:** As for "Drift (dark)". From the measured Pi and the fitted curve  $F_i$  the noise curve  $N_i = P_i - F_i$  is calculated.

The noise N is

 $N = max_i\{N_i\} - min_i\{N_i\}$ 

NOTE

Note: The noise pp specification (rather than noise 2 sigma) is given for legacy reasons.

# Noise 2 sigma (dark)

Specifies twice the standard deviation of the measured dark power (with zero input power), excluding "Drift (dark)". Power are expressed in Watts.

**Conditions:** All ports dark. Averaging time as specified. Other conditions as specified.

**Measurement:** As for "Drift (dark)" (see, "Drift (dark)" on page 57). From the measurement result Pi and the fitted curve Fi the noise N is calculated by

 $N = 2*StDev_i\{P_i - F_i\}$ , where  $StDev_i\{\}$  denotes the standard deviation over indices i.

# **Operating humidity**

Humidity range where the instrument is designed to be operated.

The instrument must not be operated outside this range. If previously stored beyond this range, wait for acclimation before turning on the instrument.

# **Operating temperature**

Temperature range for which the specifications apply if not differently stated.

The instrument must not be operated outside this range. If previously stored beyond this range, wait for acclimation before turning on the instrument.

# **Operating temperature and humidity**

The ambient temperature range and humidity range of the tunable laser source for which the specifications apply.

NOTE

If the tunable laser source module is rack- mounted the temperature and humidity within the rack apply.

# **Output** isolation

The insertion loss of the built-in isolator in the backward direction.

**Measurement:** This characteristic cannot be measured from outside the module. It is based on known isolator characteristics.

# **Optical path blocking**

Maximum possible attenuation.

# Polarization dependent loss (PDL)

The dependence of the "Total loss" (see, "Total loss" on page 72) on the input polarization state, expressed as the full difference (in dB) between the highest and the lowest total loss.

**Conditions:** Temperature range, input mode as specified.

**Measurement:** Using high quality connectors in perfect optical condition.

NOTE

PDL is a peak-to-peak value.

# Polarization dependent responsivity (PDR)

Specifies the dependence of the measured power on the input state of polarization (SOP) as ± half the difference between the highest and the lowest measured power over all SOPs, expressed in dB.

Powers expressed in dBm.


Conditions: As specified.

#### Polarization extinction ratio

Specifies the ratio of the optical power in the slow axis of a connected polarization-maintaining fiber to optical power in the fast axis, expressed in  $\ensuremath{\mathrm{d}} B$ 

**Conditions:** Applicable to tunable laser sources utilizing polarization maintaining fiber that has its TE mode in the slow axis and aligned with the connector key.

**Measurement:** Using a polarization analyzer at the end of a polarization-maintaining patchcord, by sweeping the wavelength to create circular traces on the Poincaré sphere. Calculate the polarization extinction ratio from the diameters of these circles.



# Port separation (static)

Specifies the difference between the power P applied to one port and the measured power D on any other port with no applied power (dark). Excluding "Noise pp (dark)" (see, "Noise pp (dark)" on page 60) and "Drift (dark)" (see, "Drift (dark)" on page 57).

Powers expressed in dBm, port separation S expressed in dB,

$$S = P - D$$

# Measured power (dBm) dark 1 2 3 4 Channel

Figure 48 Port separation (static)

**Conditions:** Applied power level as specified.

# Port separation, dynamic


Specifies the difference between step power input  $P_{step}$  (with positive or negative slope) applied to one port and the maximum measured power change  $D_{chng}$  (positive or negative) on any other port without applied power (dark). Excluding "Noise pp (dark)" (see, "Noise pp (dark)" on page 60) and "Drift (dark)" (see, "Drift (dark)" on page 57). Powers expressed in dBm, dynamic port separation SD expressed in dB.

$$SD = P_{step} - D_{chng}$$

Conditions: Power step height as specified (1 mW).

# Power flatness versus wavelength

Specifies  $\pm$  half the span (in dB) between the maximum and the minimum actual power levels of the tunable laser source when changing the wavelength.



**Conditions:** Uninterrupted tunable laser source output power, constant power setting, constant temperature.

# **Power repeatability**

The uncertainty in reproducing the power level after changing and resetting the power level. The power repeatability is  $\pm$  half the span between the highest and lowest actual power (in dBm).


**Conditions:** Uninterrupted tunable laser source output power, constant wavelength setting, constant temperature.

NOTE

The long-term power repeatability can be obtained by taking the power repeatability and power stability into account.

## **Power stability**

Specifies the change of the power level of the tunable laser source over time, expressed as  $\pm$  half the span (in dB) between the highest and lowest actual power.



**Conditions:** Time span as specified. Uninterrupted tunable laser source output power, constant wavelength and power level settings, constant temperature.

# **Power range**

Specifies the range from the smallest input power that causes a significant change of the measured (displayed) power to the highest power for which the specifications apply if not differently stated.

**Measurement:** The lower limit corresponds to half the "Noise pp (dark)" (see, "Noise pp (dark)" on page 60).

# Power setting mode

Operating mode where the user sets the desired output power. The power control function of the attenuator automatically adjusts the attenuation to get the set output power (see also, "Attenuation setting mode" on page 56).

## Range in attenuation setting mode

Attenuation range that can be set at the attenuator in the Attenuation Setting Mode and for which the specifications apply (if not differently stated).

#### Range in power setting mode

Range of optical output power that can be set at the attenuator in the Power Setting Mode and for which the specifications apply (if not differently stated).

NOTE

The actual minimum value is determined by the optical input power minus the maximum attenuation.

# Range settling time (one range step)

When changing from one power meter range to the next (one range higher or lower), the settling time specifies the time required for the power meter to settle into the new range. Measurements starting after that period are within specifications.

The time starts when the instrument receives the range change command (manual range) or when the instrument detects the requirement for power meter range change (autorange).

Conditions: As specified.

# Relative intensity noise (RIN)

Specifies the ratio between the mean-square of the optical power fluctuation amplitude Pf,B within a specified frequency range f and for bandwidth B, and the square of the average optical power  $P_{avg}$ .

$$RIN = \frac{\left\langle \Delta P_{f,B}^2 \right\rangle}{\left| P_{\text{avg}}^2 \cdot B \right|} \left[ \frac{1}{\text{Hz}} \right]$$

RIN, if expressed as "dB/Hz", is calculated by:

$$RIN_{\text{dB/Hz}} = 10 \log \left( \frac{\Delta P_{f,B}^2 \cdot 1 \text{Hz}}{P_{\text{avg}}^2 \cdot B} \right)$$

Conditions: As specified.

**Measurement:** Using an Agilent Lightwave signal analyzer and bandwidth set to 3 MHz.

# Relative port to port uncertainty

Specifies  $\pm$  half the maximum difference between the measured powers of different ports when applying the same stimulus. Measured powers expressed in dBm, the relative port to port uncertainty U expressed in  $\pm \mathrm{dB}$ .




Figure 49 Relative port to port uncertainty

**Conditions:** As specified.

## Relative uncertainty of monitor power meter

When changing the output power of the attenuator, the relative power meter uncertainty is the maximum error of the displayed output power ratio to the actual output power ratio. This uncertainty is caused by the internal power meter's nonlinearity and noise and by errors in the sensing hardware, expressed as  $\pm$  half the span of all possible errors with an offset due to the noise level of the power meter. Symbol RU.

**Conditions:** Reference power level 1mW, wavelength and polarization state constant; power limitations as in the Specifications; zeroing prior to measurement.

**Measurement:** The nonlinearity is calculated using the displayed power ratio D2 / D1 and the corresponding actual power ratio P2 / P1:

$$N_{dB}(1,2) = 10 \log \left( \frac{D_1/D_2}{P_2/P_1} \right)$$

The relative power meter uncertainty is then calculated using:

$$RU_{dB} = \pm \frac{max \left| N_{dB}(1,2) \right|}{2}$$


RU is then expressed in dB and the offset is expressed in pW.

NOTE

Absolute power accuracy attainable with the help of an external optical power meter.

# Relative wavelength (frequency) accuracy

When randomly changing the wavelength (frequency) of the tunable laser source and measuring the differences between the displayed and the actual wavelength (frequency), the relative wavelength (frequency) accuracy is  $\pm$  half the span between the maximum and the minimum value of all differences.



**Conditions:** Uninterrupted tunable laser source output power, constant power setting, constant temperature.

# Repeatability (of attenuation/total loss, or power level)

Specifies the uncertainty in reproducing the "Total loss" on page 72 (or "Attenuation" on page 55, or power level) after randomly changing

and resetting the attenuation or power level. The repeatability is defined as  $\pm$  half the difference between the highest and the lowest measured attenuation or power level over the repeated changing and resetting of the attenuation or power level.




Figure 50 Repeatability

**Conditions:** Constant operating conditions. Operating in power setting mode requires zeroing prior to measurement.

#### Resolution

The smallest possible increment or decrement in Attenuation Setting Mode or Power Setting Mode.

## **Return loss**

The ratio between the incident power and the power reflected back from the instrument into the input fiber, expressed in dB.

**Conditions:** As specified. (For attenuators: applicable to both ports with fiber cables with high quality connectors in perfect optical condition on both attenuator ports and the respective second port terminated (zero reflectance).)

**Measurement:** With return loss meter and non-coherent source of applicable wavelength.

NOTE

The measurement result includes internal reflections in the instrument, such as reflections from both attenuator ports.

# **SBS** suppression - effective linewidth

Specifies the peak-to-peak change of the periodically modulated wavelength resulting from the SBS suppression feature, expressed in Hertz.

# Side-mode suppression ratio

The ratio of optical power in the main mode to the optical power of the highest sidemode, expressed in dB:

$$SSR_{dB} = 10 \cdot \log \left( \frac{P_{signal}}{P_{highestsidemode}} \right)$$

**Conditions:** Within a distance from 0.1 to 6 GHz to the signal's optical frequency.


**Measurement:** Using the Agilent Lightwave signal analyzer, by analyzing the heterodyning between the main signal and the highest sidemode.

# **Settling Time**

Maximum time needed to change the attenuation or power level by a specified step, from the beginning until the end of the change.

#### Signal to source spontaneous emission ratio

Specifies the ratio between signal power and maximum spontaneous emission (SSE) power. The SSE power is determined in a specified bandwidth within a  $\pm 3$  nm window around the signal wavelength, where  $\pm 1$  nm around the signal wavelength are excluded, expressed in dB per nm.



**Conditions:** As specified.

**Measurement:** Using an optical spectrum analyzer at 0.5 nm resolution bandwidth (to address the possibility of higher SSE within a narrower bandwidth), then extrapolated to 1 nm bandwidth.

## Specification wavelength range

The wavelength range for which all specifications and characteristics apply, if not differently stated.

## Spectral ripple (power meter; due to interference)

A coherent input signal causes optical interference 1 between reflective interfaces within the power meter's optical assembly, including the end face of a connected fiber. Optical interference causes a wavelength-dependent periodic change of the power meter's responsivity. Spectral ripple is defined as  $\pm$  half the difference between the highest and the lowest power resulting from the periodic change.

**Conditions:** Constant input power, constant state of polarization (SOP). Connector as specified. Other conditions as specified.

**Measurement:** Using a laser source without spectral ripple or with known (monitored) spectral ripple and linewidth < 100 MHz and stepping the wavelength over the specified range. The periodic part in the measurement result is separated by subtracting the polynomial 2nd order fit (low frequency part).

# **Storage conditions**

Allowed temperature and humidity range for the non-operating instrument. Wait for acclimation to within the "Operating temperature" before turning on the instrument.

#### **Total loss**

The change of power level after inserting the attenuator between two connectorized patchcords, at an arbitrary attenuation setting, expressed in dB. Its value TL can be calculated from:

$$TL(dB) = 10\log\frac{P_a}{P_b} = P_a(dBm) - P_b(dBm)$$

where:

 $P_a$  = power measured at the end of the two patchcords.

 $P_b$  = power measured after the insertion of the attenuator.




Figure 51 Total loss

**Conditions:** Patchcord cables with high quality connectors in perfect optical condition on both attenuator ports.

NOTE

The total loss depends on the attenuation setting.

NOTE

Total loss includes the loss of one additional connector pair.

#### **Total uncertainty**

The maximum relative difference U between measured (displayed) power M and actual (true) power P for a specified set of operating conditions, expressed in percent. Powers expressed in Watts. Excluding "Noise pp (dark)" (see, "Noise pp (dark)" on page 60) and "Drift (dark)" (see, "Drift (dark)" on page 57).

```
U = 100\% * (M/P - 1)
```

**Conditions:** Operating conditions as specified. Other conditions as specified.

#### **Transition speed**

Specifies the settable attenuation change rate of the attenuator.

#### Tuning time

Specifies the time needed to tune to another wavelength.

**Conditions:** Tuning distance as specified. The time is measured from sending the command to the TLS until the tuning operation has finished.

#### **Uncertainty at reference conditions**

The maximum relative difference U between measured (displayed) power M and actual (true) power P for the specified set of reference conditions, expressed in percent. Powers expressed in Watts.

```
U = 100\% * (M/P - 1)
```

Conditions: Reference conditions as specified.

#### Wavelength range

The range of wavelengths for which the power meter is calibrated, or that can be set at the attenuator and for which the specifications apply (if not differently stated).

#### Wavelength (frequency) repeatability

The random uncertainty in reproducing a wavelength (frequency) of the tunable laser source after changing and re-setting the wavelength (frequency). The wavelength (frequency) repeatability is ± half the span

#### 3 Specifications and Regulations Compliance

between the maximum and the minimum of all actual values of this wavelength (frequency).

**Conditions:** Uninterrupted tunable laser source output power, constant power setting, constant temperature.

NOTE

The long-term wavelength repeatability can be obtained by taking the wavelength repeatability and wavelength stability into account.

#### **Wavelength resolution**

The smallest selectable wavelength (frequency) increment or decrement.

#### Wavelength (frequency) stability

Specifies the change of the actual wavelength (frequency) of the tunable laser source over time, expressed as  $\pm$  half the span between the maximum and minimum of all wavelengths (frequencies).

**Conditions:** Time span as specified, uninterrupted tunable laser source output power, constant wavelength and power level settings, constant temperature.

#### Warm-up time

Time after power-up of the acclimated instrument after which the specifications and characteristics apply.

#### Literature

 $1\ {\rm Fiber}$  Optic Test and Measurement, Hewlett Packard Professional Books, edited by Prentice Hall, ISBN 0-13-534330-5

### **Technical Specifications - Optical Variable Attenuator and Power Meter**

### **Optical Attenuator**

N7751A, N	17752A, N7761A, N7762A, N7764A	
Connectivity	FC/APC angled (Option -022) or	
	FC/PC straight (Option -021) contact co	nnector interface
Fiber type	9/125 µm SMF 28	
Wavelength range	1260 nm to 1640 nm	
Attenuation range	0 dB to 40 dB (45 dB typ.)	
Attenuation resolution	0.01 dB	
	Attenuation setting mode	Power setting mode
Range	0 dB to 40 dB	-50 dBm to +20 dBm
Resolution	0.01 dB	0.01 dB
Repeatability <sup>1</sup>	typ. +/- 0.05 dB for attenuation 0 dB to 30 dB	± 0.025 dB
	typ. +/- 0.10 dB for attenuation 30 dB to 40 dB	
Accuracy (uncertainty) <sup>1, 2, 3</sup>	typ. $\pm$ 0.10 dB for attenuation 0 dB to 10 dB	
	typ. $\pm$ 0.15 dB for attenuation 10 dB to 20 dB	
	typ. $\pm$ 0.40 dB for attenuation 20 dB to 40 dB	
Relative accuracy (uncertainty) 1,2,4	± 0.05 dB ± 300	
Polarization dependent loss 5	typ. $\leq$ 0.15 dB <sub>PP</sub> for attenuation 0 dB to 10 dB	
	typ. $\leq$ 0.25 dB <sub>PP</sub> for attenuation 10 dB to 20 dB	≤ 0.15 dB <sub>PP</sub> <sup>9</sup>
	typ. $\leq 0.5 \text{ dB}_{pp}$ for attenuation 20 dB to 40 dB	
Settling time	typ. 20 ms <sup>10</sup>	typ. 100 ms <sup>9</sup>
Insertion loss <sup>2, 6</sup>	typ. ≤ 1.2 dB (excluding connectors)	
	≤ 2.2 dB (including connectors) <sup>7</sup>	
Attenuation transition speed	selectable from 0.1 to 1000 dB/s	
Relative uncertainty of monitor power meter <sup>2, 8</sup>	±0.05 dB ± 300 pW	
Averaging time of monitor power meter	2 ms to 1 s	
Return loss	typ. 45 dB	
Maximum safe input power	+23 dBm	
Optical path blocking	typ. 45 dB	

- 1. At constant operating conditions.
- 2. For unpolarized light. Temperature constant and between (23  $\pm$  5) °C.
- 3. For (1310  $\pm$ 15) nm, (1490  $\pm$ 10) nm and (1550  $\pm$ 15) nm.
- 4. Output power > -40 dBm, input power < +10 dBm. For input power > +10 dBm add typ. ± 0.02 dB.
- 5. For (1550  $\pm$ 15) nm. Add typ. 0.05 dB for (1310  $\pm$ 15) nm and (1490  $\pm$ 10) nm.
- 6. For (1550  $\pm$ 15) nm. Add typ. 0.20 dB for (1310  $\pm$ 15) nm and (1490  $\pm$ 10) nm.
- 7. Measured with Agilent reference connectors.
- 8. For (1550 ±15) nm.
- 9. Output power > -30 dBm, input power < +10 dBm.
- 10. For 20 dB step.

### **Optical Power Meter**

N7751A and N7752A		
Sensor element	InGaAs	
Wavelength range	1260 nm to 1640 nm	
Specification wavelength range	(1310 ±15) nm, (1490 ±10) nm, (1550 ±15) nm	
Power range	-80 dBm to +10 dBm	
Maximum safe power	+16 dBm	
Averaging time	2 ms to 1s	
Applicable fiber type	Standard SM and MM $\leq$ 62.5 $\mu$ m core size, NA $\leq$ 0.24	
Uncertainty at reference conditions <sup>11,13</sup>	±2.5%	
Total uncertainty 12, 15, 16	±4.5%	
Linearity 17		
at (23 ± 5)°C	±0.02 dB	
over operating temperature	±0.04 dB	
Polarization dependent responsivity (PDR) <sup>13,18</sup>	typ. < ±0.01 dB (1260 nm to 1580 nm)	
Spectral ripple (due to interference) <sup>20</sup>	typ. < ±0.01 dB	
Drift (dark) <sup>14</sup>	±9 pW	
Noise pp (dark) <sup>13</sup> , (1s averaging time, 300s observation time)	< 7 pW <sub>pp</sub>	
Return loss <sup>19</sup>	typ. > 57 dB	

11. Reference conditions:

reterence conditions:
Single mode fiber SMF 9 µm. Power Level: –20 dBm to 0 dBm.
On day of calibration (add ±0.3% for aging over one year; add ±0.6% for aging over two years).
Spectral width of source < 10 nm FWHM.
Wavelength setting of power sensor corresponds to source wavelength ± 0.4 nm.

Operating conditions:

Single mode fiber SMF. For multimode fiber, typical.

Within one year of calibration; add  $\pm$  0.3% for second year.

Spectral width of source < 10 nm FWHM.

Wavelength setting of power sensor corresponds to source wavelength  $\pm$  0.4 nm.

- Ambient temperature (23 ±5) °C
- Temperature constant within ±1 K after zeroing
- 15. Excluding noise and offset drift
- 16. Power range -60 dBm to +10 dBm
- Power range -40 dBm to +10 dBm. For power <-40 dBm add 0.05 dB.
- Straight connector, SMF
- Angled connector 8°, ceramic ferrule, SMF
- For constant state of polarization, source linewidth < 100 MHz, angled connector 8°, wavelength range 1260 nm to 1625 nm. Add ±0.01 dB typ. within specification wavelength range for straight connector with ceramic ferrule.

### **General Characteristics**

#### Dimensions (H x W x D):

1U half-rack

372 mm x 212 mm x 43 mm (excluding front and back rubber cushions)

#### Weight:

Approx. 3 kg

#### Recommended recalibration period

24 month

#### Operating temperature

+5 °C to +40 °C

#### Operating humidity

15% to 95%, non-condensing

The maximum operating altitude is 2000 m.

#### Pollution protection

The Agilent N775xA/6xA is designed for pollution degree 2.

#### Warm-up time

20 minutes

#### Interfaces

The instruments can be controlled via LAN, USB or GPIB interfaces

#### Power consumption

Line power: AC 100 - 240 V  $\pm$  10%, 50/60 Hz, 60 VA max.

### **Declarations of Conformity - Optical Multiport Attenuator and Power** Meter



Manufacturer's Name: Agilent Technologies Singapore (International) Pte. Ltd.

No. 1 Yishun Ave 7 SINGAPORE 768923 Manufacturer's Address:

Singapore

Declares under sole responsibility that the product as originally delivered

Product Name: Optical Multiport Attenuator and Power Meter 1-Channel Optical Attenuator and Power Meter 2-Channel Optical Attenuator and Power Meter Product Number:

N7752A 1-Channel Optical Attenuator N7761A N7762A 2-Channel Optical Attenuator

4-Channel Optical Attenuator 2-Channel Multimode Optical Attenuator N7764A N7766A N7768A 4-Channel Multimode Optical Attenuator This declaration covers all options of the above system.

complies with the essential requirements of the following applicable European Directives, and carries the CE marking accordingly:

The Low Voltage Directive 2006/95/EC

The EMC Directive 2004/108/EC

Product Options:

#### and conforms with the following product standards:

Standard IEC 61326-1:2005 / EN 61326-1:2006 EMC

IEC 61326-1:2005 / EN 61326-1:2006 Reference Standards CISPR 11-2003 / EN 55011:1998+A1:1999+A2:2002 IEC 61000-4-2: 2001 / EN 61000-4-2:1995+A1:1998+A2:2001 IEC 61000-4-3:2002+A1:2002 / EN 61000-4-3:2002+A1:2002 IEC 61000-4-5: 2001 / EN 61000-4-5:1995+A1:2001 IEC 61000-4-6:2003 / EN 61000-4-6:1996+A1:2001 IEC 61000-4-1:2004 / EN 61000-4-11:2004

Canada: ICES/NMB-001:2004 Australia/New Zealand: AS/NZS CISPR 11:2004

IEC 61010-1:2001 / EN 61010-1:2001 Safety

Canada: CAN/CSA-C22.2 No. 61010-1:2004 USA: UL 61010-1:2004

#### Supplementary Information:

The products were tested in a typical configuration with Agilent Technologies test systems.

This DoC applies to above-listed products placed on the EU market after:

2010-May-02 Date

Product Quality & Compliance

triber

For further information, please contact your local Agilent Technologies sales office, agent or distributor or Agilent Technologies Deutschland GmbH, Herrenberger Straße 130, 71034 Boeblingen, Germany.

Template:A5971-5302-2, Rev. F.00

Year of first marking: 2009

Revision: B

Limit

Group 1 Class A 4 kV CD, 4 kV AD 3 V/m / 80 MHz-2 GHz, 1V/m / 2-2.7 GHz 0.5 kV signal lines, 1 kV power lines 0.5 kV line-line, 1 kV line-ground 3 V, 0.15-80 MHz 0% for 1/0.5 (0°, 180°) cycle 0% for 250/300 cycles 70% for 25/30 cycles

### **Regulations Information**

#### **EMC Canada**

These ISM devices comply with Canadian ICES-001. Ces appareils ISM sont conformes à la norme NMB-001 du Canada.

### **Acoustic Noise Information Germany**

ISO 7779  $L_{pA}$  < 70 dB(A), normal operation, operator position.

### Remove all doubt

#### Remove all doubt

Our repair and calibration services will get your equipment back to you, performing like new, when promised. You will get full value out of your Agilent equipment throughout its lifetime. Your equipment will be serviced by Agilent-trained technicians using the latest factory calibration procedures, automated repair diagnostics and genuine parts. You will always have the utmost confidence in your measurements. Agilent offers a wide range of additional expert test and measurement services for your equipment, including initial start-up assistance, onsite education and training, as well as design, system integration, and project management. For more information on repair and calibration services, go to

www.agilent.com/find/removealldoubt

### **Agilent E-mail Updates**

Get the latest information on the products and applications you select.



www.agilent.com/find/emailupdates

### **Agilent Direct**

Quickly choose and use your test equipment solutions with confidence.



www.agilent.com/find/agilentdirect

### **Agilent Open**

Agilent Open simplifies the process of connecting and programming test systems to help engineers design, validate and manufacture electronic products. Agilent offers open connectivity for a broad range of systemready instruments, open industry software, PC-standard I/O and global support, which are combined to more easily integrate test system development.



www.agilent.com/find/open

#### LXI

LXI is the LAN-based successor to GPIB, providing faster, more efficient connectivity. Agilent is a founding member of the LXI consortium.




www.lxistandard.org

### www.agilent.com

For more information on Agilent Technologies products, applications and services, please contact your local Agilent office. The complete list  $% \left( 1\right) =\left( 1\right) \left( 1\right)$ is available at:

www.agilent.com/find/contactus

**3** Specifications and Regulations Compliance



### **Accessories - Overview**

The instruments are available in various configurations for the best possible match to the most common applications.

This chapter provides information on the available options and accessories.

- "Instrument and Options Agilent N7751A, N7752A, N7761A, N7762A, N7764A, N7766A, N7768A" on page 84
- "Cables" on page 84

# Instrument and Options - Agilent N7751A, N7752A, N7761A, N7762A, N7764A, N7766A, N7768A

All systems have 1 year warranty

Model number	
N7751A	Optical Attenuator (1 Channel) with Optical Power Meter (2 Channels)
N7752A	Optical Attenuator (2 Channels) with Optical Power Meter (2 Channels)
N7761A	Optical Attenuator (1 Channel)
N7762A	Optical Attenuator (2 Channels)
N7764A	Optical Attenuator (4 Channels)
N7766A	Optical Attenuator (2 Channels Multimode)
N7768A	Optical Attenuator (4 Channels Multimode)

Connector interface option		
-021	Straight contact connector (N7751A, N7752A, N7761A, N7762A, N7764A)	
-022	Angled contact connector (N7751A, N7752A, N7761A, N7762A, N7764A)	
-050	50/125 μm MMF Interface (N7766A, N7768A)	
-062	62.5/125 µm MMF Interface (N7766A, N7768A)	

Accessories		
N7744-100	Rack Mount Kit for 1 or 2 Units	

Warranty	
All systems have 1	year warranty
R-51B-001-3C	1 year Return-to-Agilent warranty extended to 3 years
R-51B-001-5C	1 year Return-to-Agilent warranty extended to 5 years

Calibration	
R-50C-011-3	Agilent calibration upfront support plan 3 year coverage
R-50C-011-5	Agilent calibration upfront support plan 5 year coverage

### **Cables**

- GPIB Cable, 10833A, 1 m (3.3 feet)
- GPIB Cable, 10833B, 2 m (6.6 feet)
- GPIB Cable, 10833C, 4 m (13.2 feet)
- GPIB Cable, 10833D, 0.5 m (1.6 feet)

- GPIB Adapter, 10834A, 2.3 cm extender. Use this adapter if there is no space to connect your GPIB cable directly to a GPIB interface.
- USB Cable, 8121-1243, 2 m (6.6 feet)

#### 4 Accessories



### **Maintenance and Troubleshooting - Overview**

This section provides detailed information on:

- · Cleaning Instructions
- Firmware Upgrades
- · Error Messages
- Data Sanitation

### **Cleaning Instructions**

### **Cleaning Instructions**

The following Cleaning Instructions contain some general safety precautions, which must be observed during all phases of cleaning. Consult your specific optical device manuals or guides for full information on safety matters.

Please try, whenever possible, to use physically contacting connectors, and dry connections. Clean the connectors, interfaces, and bushings carefully after use.

If you are unsure of the correct cleaning procedure for your optical device, we recommend that you first try cleaning a dummy or test device.

Agilent Technologies assume no liability for the customer's failure to comply with these requirements.



### **Cleaning Instructions for this Instrument**

The Cleaning Instructions apply to a number of different types of Optical Equipment. Most relevant for this instrument is the section "How to clean connector adapters" on page 97 for cleaning the quadadapters, and "How to clean optical glass plates" on page 98 for cleaning the front of the instrument.

### **Safety Precautions**

Please follow the following safety rules:

- Do not remove instrument covers when operating.
- Ensure that the instrument is switched off throughout the cleaning procedures.
- Use of controls or adjustments or performance of procedures other than those specified may result in hazardous radiation exposure.
- Make sure that you disable all sources when you are cleaning any optical interfaces.
- Under no circumstances look into the end of an optical device attached to optical outputs when the device is operational. The laser radiation is not visible to the human eye, but it can seriously damage your eyesight.
- To prevent electrical shock, disconnect the instrument from the mains before cleaning. Use a dry cloth, or one slightly dampened with water, to clean the external case parts. Do not attempt to clean internally.
- Do not install parts or perform any unauthorized modification to optical devices.
- Refer servicing only to qualified and authorized personnel.

### Why is it important to clean optical devices?

In transmission links optical fiber cores are about 9m~(0.00035") in diameter. Dust and other particles, however, can range from tenths to hundredths of microns in diameter. Their comparative size means that they can cover a part of the end of a fiber core, and as a result will reduce the performance of your system.

Furthermore, the power density may burn dust into the fiber and cause additional damage (for example, 0 dBm optical power in a single mode fiber causes a power density of approximately 16 million W/m2). If this happens, measurements become inaccurate and non-repeatable.

Cleaning is, therefore, an essential yet difficult task. Unfortunately, when comparing most published cleaning recommendations, you will discover that they contain several inconsistencies. In this section, we want to suggest ways to help you clean your various optical devices, and thus significantly improve the accuracy and repeatability of your lightwave measurements.

### What do I need for proper cleaning?

Some Standard Cleaning Equipment is necessary for cleaning your instrument. For certain cleaning procedures, you may also require certain Additional Cleaning Equipment.

#### **Standard Cleaning Equipment**

Before you can start your cleaning procedure you need the following standard equipment:

- · Dust and shutter caps
- Isopropyl alcohol
- Cotton swabs
- Soft tissues
- Pipe cleaner
- Compressed air

#### **Dust and shutter caps**

All of Agilent Technologies' lightwave instruments are delivered with either laser shutter caps or dust caps on the lightwave adapter. Any cables come with covers to protect the cable ends from damage or contamination.

We suggest these protected coverings should be kept on the equipment at all times, except when your optical device is in use. Be careful when replacing dust caps after use. Do not press the bottom of the cap onto the fiber too hard, as any dust in the cap can scratch or pollute your fiber surface.

If you need further dust caps, please contact your nearest Agilent Technologies sales office.

#### Isopropyl alcohol

This solvent is usually available from any local pharmaceutical supplier or chemist's shop.

If you use isopropyl alcohol to clean your optical device, do not immediately dry the surface with compressed air (except when you are cleaning very sensitive optical devices). This is because the dust and the dirt is solved and will leave behind filmy deposits after the alcohol is evaporated. You should therefore first remove the alcohol and the dust with a soft tissue, and then use compressed air to blow away any remaining filaments.

If possible avoid using denatured alcohol containing additives. Instead, apply alcohol used for medical purposes. Never try to drink this alcohol, as it may seriously damage to your health.

Do not use any other solvents, as some may damage plastic materials and claddings. Acetone, for example, will dissolve the epoxy used with fiber optic connectors. To avoid damage, only use isopropyl alcohol.

#### **Cotton swabs**

We recommend that you use swabs such as Q-tips or other cotton swabs normally available from local distributors of medical and hygiene products (for example, a supermarket or a chemist's shop). You may be able to obtain various sizes of swab. If this is the case, select the smallest size for your smallest devices.

Ensure that you use natural cotton swabs. Foam swabs will often leave behind filmy deposits after cleaning.

Use care when cleaning, and avoid pressing too hard onto your optical device with the swab. Too much pressure may scratch the surface, and could cause your device to become misaligned. It is advisable to rub gently over the surface using only a small circular movement.

Swabs should be used straight out of the packet, and never used twice. This is because dust and dirt in the atmosphere, or from a first cleaning, may collect on your swab and scratch the surface of your optical device.

#### **Soft tissues**

These are available from most stores and distributors of medical and hygiene products such as supermarkets or chemists' shops.

We recommend that you do not use normal cotton tissues, but multilayered soft tissues made from non-recycled cellulose. Cellulose tissues are very absorbent and softer. Consequently, they will not scratch the surface of your device over time.

Use care when cleaning, and avoid pressing on your optical device with the tissue. Pressing too hard may lead to scratches on the surface or misalignment of your device. Just rub gently over the surface using a small circular movement.

Use only clean, fresh soft tissues and never apply them twice. Any dust and dirt from the air which collects on your tissue, or which has gathered after initial cleaning, may scratch and pollute your optical device.

#### Pipe cleaner

Pipe cleaners can be purchased from tobacconists, and come in various shapes and sizes. The most suitable one to select for cleaning purposes has soft bristles, which will not produces scratches.

There are many different kinds of pipe cleaner available from tobacco shops.

The best way to use a pipe cleaner is to push it in and out of the device opening (for example, when cleaning an interface). While you are cleaning, you should slowly rotate the pipe cleaner.

Only use pipe cleaners on connector interfaces or on feed through adapters. Do not use them on optical head adapters, as the center of a pipe cleaner is hard metal and can damage the bottom of the adapter.

Your pipe cleaner should be new when you use it. If it has collected any dust or dirt, this can scratch or contaminate your device.

The tip and center of the pipe cleaner are made of metal. Avoid accidentally pressing these metal parts against the inside of the device, as this can cause scratches.

#### **Compressed air**

Compressed air can be purchased from any laboratory supplier.

It is essential that your compressed air is free of dust, water and oil. Only use clean, dry air. If not, this can lead to filmy deposits or scratches on the surface of your connector. This will reduce the performance of your transmission system.

When spraying compressed air, hold the can upright. If the can is held at a slant, propellant could escape and dirty your optical device. First spray into the air, as the initial stream of compressed air could contain

some condensation or propellant. Such condensation leaves behind a filmy deposit.

Please be friendly to your environment and use a CFC-free aerosol.

#### **Additional Cleaning Equipment**

Some Cleaning Procedures need the following equipment, which is not required to clean each instrument:

- Microscope with a magnification range about 50X up to 300X
- Ultrasonic bath
- · Warm water and liquid soap
- · Premoistened cleaning wipes
- · Polymer film
- · Infrared Sensor Card

#### Microscope with a magnification range about 50X up to 300X

A microscope can be found in most photography stores, or can be obtained through or specialist mail order companies. Special fiberscopes are available from suppliers of splicing equipment.

Ideally, the light source on your microscope should be very flexible. This will allow you to examine your device closely and from different angles.

A microscope helps you to estimate the type and degree of dirt on your device. You can use a microscope to choose an appropriate cleaning method, and then to examine the results. You can also use your microscope to judge whether your optical device (such as a connector) is severely scratched and is, therefore, causing inaccurate measurements.

#### **Ultrasonic bath**

Ultrasonic baths are also available from photography or laboratory suppliers or specialist mail order companies.

An ultrasonic bath will gently remove fat and other stubborn dirt from your optical devices. This helps increase the life span of the optical devices.

Only use isopropyl alcohol in your ultrasonic bath, as other solvents may damage.

#### Warm water and liquid soap

Only use water if you are sure that there is no other way of cleaning your optical device without corrosion or damage. Do not use hot water, as this may cause mechanical stress, which can damage your optical device.

Ensure that your liquid soap has no abrasive properties or perfume in it. You should also avoid normal washing-up liquid, as it can cover your device in an iridescent film after it has been air-dried.

Some lenses and mirrors also have a special coating, which may be sensitive to mechanical stress, or to fat and liquids. For this reason we recommend you do not touch them.

If you are not sure how sensitive your device is to cleaning, please contact the manufacturer or your sales distributor.

#### Premoistened cleaning wipes

Use pre-moistened cleaning wipes as described in each individual cleaning procedure. Cleaning wipes may be used in every instance where a moistened soft tissue or cotton swab is applied.

#### Polymer film

Polymer film is available from laboratory suppliers or specialist mail order companies.

Using polymer film is a gentle method of cleaning extremely sensitive devices, such as reference reflectors and mirrors.

#### Infrared Sensor Card

Infrared sensor cards are available from laboratory suppliers or specialist mail order companies.

With this card you are able to control the shape of laser light emitted. The invisible laser beam is projected onto the sensor card, then becomes visible to the normal eye as a round spot.

Take care never to look into the end of a fiber or any other optical component, when they are in use. This is because the laser can seriously damage your eyes.

### **Preserving Connectors**

Listed below are some hints on how best to keep your connectors in the best possible condition.

#### **Making Connections**

Before you make any connection you must ensure that all cables and connectors are clean. If they are dirty, use the appropriate cleaning procedure.

When inserting the ferrule of a patchcord into a connector or an adapter, make sure that the fiber end does not touch the outside of the mating connector or adapter. Otherwise you will rub the fiber end against an unsuitable surface, producing scratches and dirt deposits on the surface of your fiber.

#### **Dust Caps and Shutter Caps**

Be careful when replacing dust caps after use. Do not press the bottom of the cap onto the fiber as any dust in the cap can scratch or dirty your fiber surface.

When you have finished cleaning, put the dust cap back on, or close the shutter cap if the equipment is not going to be used immediately.

Keep the caps on the equipment always when it is not in use.

All of Agilent Technologies' lightwave instruments and accessories are shipped with either laser shutter caps or dust caps. If you need additional or replacement dust caps, contact your nearest Agilent Technologies Sales/Service Office.

#### **Immersion Oil and Other Index Matching Compounds**

Where it is possible, do not use immersion oil or other index matching compounds with your device. They are liable to impair and dirty the surface of the device. In addition, the characteristics of your device can be changed and your measurement results affected.

### **Cleaning Instrument Housings**

Use a dry and very soft cotton tissue to clean the instrument housing and the keypad. Do not open the instruments as there is a danger of electric shock, or electrostatic discharge. Opening the instrument can cause damage to sensitive components, and in addition your warranty will be voided.

Never open the instruments as they can be damaged. Opening the instruments puts you in danger of receiving an electrical shock from your device, and renders your warranty void.

### Which Cleaning Procedure should I use?

#### Light dirt

If you just want to clean away light dirt, observe the following procedure for all devices:

- 1 Use compressed air to blow away large particles.
- **2** Clean the device with a dry cotton swab.
- 3 Use compressed air to blow away any remaining filament left by the swab.

#### **Heavy dirt**

If the above procedure is not enough to clean your instrument, follow one of the procedures below. Please consult "Cleaning Instructions for this Instrument" on page 88 for the procedure relevant for this instrument.

If you are unsure of how sensitive your device is to cleaning, please contact the manufacturer or your sales distributor

#### How to clean connectors

Cleaning connectors is difficult as the core diameter of a single-mode fiber is only about 9m. This generally means you cannot see streaks or scratches on the surface. To be certain of the condition of the surface of your connector and to check it after cleaning, you need a microscope.

In the case of scratches, or of dust that has been burnt onto the surface of the connector, you may have no option but to polish the connector. This depends on the degree of dirtiness, or the depth of the scratches. This is a difficult procedure and should only be performed by skilled personal, and as a last resort as it wears out your connector.

#### WARNING

Never look into the end of an optical cable that is connected to an active source.

To assess the projection of the emitted light beam you can use an infrared sensor card. Hold the card approximately 5 cm from the output of the connector. The invisible emitted light is project onto the card and becomes visible as a small circular spot.

#### **Preferred Procedure**

Use the following procedure on most occasions.

- 1 Clean the connector by rubbing a new, dry cotton-swab over the surface using a small circular movement.
- **2** Blow away any remaining lint with compressed air.

#### **Procedure for Stubborn Dirt**

Use this procedure particularly when there is greasy dirt on the connector:

- 1 Moisten a new cotton-swab with isopropyl alcohol.
- **2** Clean the connector by rubbing the cotton-swab over the surface using a small circular movement.
- **3** Take a new, dry soft-tissue and remove the alcohol, dissolved sediment and dust, by rubbing gently over the surface using a small circular movement.
- 4 Blow away any remaining lint with compressed air.

#### **An Alternative Procedure**

A better, more gentle, but more expensive cleaning procedure is to use an ultrasonic bath with isopropyl alcohol.

- 1 Hold the tip of the connector in the bath for at least three minutes.
- **2** Take a new, dry soft-tissue and remove the alcohol, dissolved sediment and dust, by rubbing gently over the surface using a small circular movement.
- **3** Blow away any remaining lint with compressed air.

### How to clean connector adapters

#### **Preferred Procedure**

Use the following procedure on most occasions.

- 1 Clean the adapter by rubbing a new, dry cotton-swab over the surface using a small circular movement.
- 2 Blow away any remaining lint with compressed air.

#### **Procedure for Stubborn Dirt**

Use this procedure particularly when there is greasy dirt on the adapter:

- **1** Moisten a new cotton-swab with isopropyl alcohol.
- **2** Clean the adapter by rubbing the cotton-swab over the surface using a small circular movement.
- **3** Take a new, dry soft-tissue and remove the alcohol, dissolved sediment and dust, by rubbing gently over the surface using a small circular movement.
- 4 Blow away any remaining lint with compressed air.

### How to clean bare fiber adapters

Bare fiber adapters are difficult to clean.

Protect from dust unless they are in use.

CAUTION

Never use any kind of solvent when cleaning a bare fiber adapter as solvents can damage the foam inside some adapters.

They can deposit dissolved dirt in the groove, which can then dirty the surface of an inserted fiber.

#### **Preferred Procedure**

Use the following procedure on most occasions.

1 Blow away any dust or dirt with compressed air.

### How to clean optical glass plates

The multiport power meter has an optical glass plate to protect the sensor. This has a special coating that is sensitive to solvents, grease, liquid and mechanical abrasion. Take extra care when cleaning this plate.

#### **Preferred Procedure**

Use the following procedure on most occasions.

- 1 Clean the plate by rubbing a new, dry cotton-swab over the surface using a small circular movement.
- **2** Blow away any remaining lint with compressed air.

#### **Procedure for Stubborn Dirt**

Use this procedure particularly when there is greasy dirt on the plate:

- 1 Moisten a new cotton-swab with isopropyl alcohol.
- **2** Clean the plate by rubbing the cotton-swab over the surface using a small circular movement.
- **3** Using a new, dry cotton-swab remove the alcohol, any dissolved sediment and dust.
- 4 Blow away any remaining lint with compressed air.

### **Additional Cleaning Information**

The following cleaning procedures may be used with other optical equipment:

· How to clean bare fiber ends

#### How to clean bare fiber ends

Bare fiber ends are often used for splices or, together with other optical components, to create a parallel beam. The end of a fiber can often be scratched. You make a new cleave. To do this:

- 1 Strip off the cladding.
- **2** Take a new soft-tissue and moisten it with isopropyl alcohol.
- **3** Carefully clean the bare fiber with this tissue.

**4** Make your cleave and immediately insert the fiber into your bare fiber adapter in order to protect the surface from dirt.

### **Other Cleaning Hints**

Selecting the correct cleaning method is an important element in maintaining your equipment and saving you time and money. This Appendix highlights the main cleaning methods, but cannot address every individual circumstance.

This section contain some additional hints which we hope will help you further. For further information, please contact your local Agilent Technologies representative.

#### Making the connection

Before you make any connection you must ensure that all lightwave cables and connectors are clean. If not, then use appropriate the cleaning methods.

When you insert the ferrule of a patchcord into a connector or an adapter, ensure that the fiber end does not touch the outside of the mating connector or adapter. Otherwise, the fiber end will rub up against something which could scratch it and leave deposits.

#### Lens cleaning papers

Note that some special lens cleaning papers are not suitable for cleaning optical devices like connectors, interfaces, lenses, mirrors and so on. To be absolutely certain that a cleaning paper is applicable, please ask the salesperson or the manufacturer.

### Firmware Upgrades

Updated versions of the firmware for your instrument may be released from time to time. These will then be available from the Agilent web page (www.agilent.com/find/octfirmware) and will include the update application and instructions.

### **Error Messages**

#### SYST:ERR?

The syst:err? command returns the next error from the error queue (refer to the Programming Guide for details).

Each error has an error code and a short description of the error, separated by comma. Please note that the error codes are returned as signed integer numbers in a range from -32768 to +32767 (INT16). Negative error numbers are defined by the SCPI standard, positive error numbers are device dependent.

#### **Data Sanitation**

### **Memory Inventory**

The instruments have the following memory size:

**EEPROM** 1Mb + 1Mb for each block of ports That is a N7745A, 2-port power meter has 3Mb. This memory does not contain any customer data.

**Flash ROM** 16MB This memory can contain calibration data, and so on, as well as user parameter settings, zeroing data and logged information (the log contains limited operational information).

**DRAM** volatile working memory. This contains data for current and recent calibrating and measuring.

### **Clearing User Data**

#### Clearing the DRAM

All data in the DRAM is cleared by powering the instrument off.

#### Clearing the Flash ROM

The user data can be cleared from the Flash ROM by:

- 1 Reset the instrument to factory default.
  - **a** Send the command:

SYSTem:PRESet

OR

- **b** Press the LAN Reset button on the front panel briefly (<1 s).
- **2** Clear the instrument log.
  - a Send the special command

DIAG:SYST:FLAS:LOGF:ERAS

The contents of the log can be checked with the query READ:LOGF?.

- **3** Clear the Ethernet settings
  - **a** Connect to the instrument by some means other than the LAN,  $using\,VISA\,Assistant\,(for\,example\,over\,USB).\,Send\,the\,command:$

SYSTem:COMMunicate:ETHernet:RESet

OR

**b** With the instrument connected over the LAN (the LAN LED is green), press and hold the LAN Reset button for more than 3s.

In either case, the instrument reboots.

**5** Maintenance and Troubleshooting

## Index

Α	Front Panel, 32	Polymer film, 93 Power Meters, 8
accessories, 83 adapters, 30	G	pre-moistened cleaning, 93
Additional Cleaning Equipment, 89 Agilent N7744A, 7	GPIB Interface, 7 GPIB interface, 30	S
Agilent N7751A, 8 Agilent N7752A, 8	T.	safety, 11 Safety Considerations, 7 Standard Cleaning Equipment, 89
Agilent N7766A, 9 Agilent N7768A, 9 Averaging Time, 52	Infrared sensor cards, 93 Initial Inspection, 7	Т
C	IP address, 26 isopropyl alcohol, 90	Temperature for operation and storage, 17
cables, 30 Cleaning Instructions, 87	L	U
Cleaning Procedures, 92 Connector, 30 connectors, 94	LAN Interface, 7 Line Power Cable, 15	Ultrasonic baths, 92 USB Interface, 7
D	M	W
DRAM, 100	microscope, 92	wavelength, 52
E	0	
EEPROM, 100 Electrical Connectors, 7	Optical Equipment, 88 options, 83	
F	P	
Flash ROM, 100	Pipe cleaners, 91	

Index

### www.agilent.com

© Agilent Technologies 2010 Fourth edition, July 2010



N7751-90B01

