

# Microfunnels

**Microfunnels** facilitate the addition of dry reagents or slurries to the Quest 210 reaction vessels. The Microfunnels can be used individually or in a group.



**Instructions.** To add a reagent to the Quest 210, remove the reaction vessel upper manifold Luer Plug and insert a Microfunnel. Add reagent slowly to prevent blockage of the Microfunnel. If the reagent sticks to the Microfunnel, rinse the Microfunnel with an appropriate solvent to transfer the remaining reagent into the reaction vessel. If the funnel becomes plugged, gently tap the Microfunnel against the upper manifold to dislodge the obstruction.

**Part Number:** 900158  
**Quantity:** 20

## Quest Accessories

P/N 102587