



BARIATRIC SOLUTIONS

B-Series




*Shower Trolley
(XXX lb. capacity)*



*Ready Lift
(700 lb. capacity)*



*I-1000 Convertible™ Chair
(I-660 also available)*



*B-1000 TransBed™ with
Positioning and Transfer System (PTS™)*

BARTON™
MEDICAL CORPORATION
www.bartonmedical.com

“Your Safest Move.”

(512) 476-7199
fax (512) 476-7190
email info@bartonmedical.com

Barton® Medical Corporation has developed *Bariatric Solutions* that include a bed, convertible transfer chair and the *Positioning and Transfer System (PTS™)*, high-capacity lifts, and bariatric shower beds. They are all designed to provide a safe, no-lift solution to the increasing number of non- and limited-ambulatory bariatric patients.

The Barton® *TransBed™* includes unique features such as the swing-away trapeze and “Helping Hand™” side-egress bar for patients up to 1000 pounds. The *I-1000* and *I-660 Convertible™* chairs, with 1000-lb. and 660-lb. capacities and the Barton® *PTS™* integrate with the *TransBed™* to allow a single caregiver to move and transfer bariatric patients safely, carefully, and efficiently. Floor Lift and Shower Bed make day-to-day tasks easy and dignified.



BARIATRIC SOLUTIONS

B-Series



Barton's® *Bariatric Solutions* address the increased risk of physical injury to both caregiver and patient when caring for bariatric patients. The Barton® *B-1000 TransBed™* with *Positioning and Transfer System (PTS™)* and *I-1000/I-660* Barton chair, allows a single caregiver to transfer a non- or limited-ambulatory bariatric patient weighing up to 1000/660 lbs. safely, easily and efficiently.

Barton® B-1000 TransBed™

- The *B-1000 TransBed™*, with a 1000-lb. capacity and 5-function system, will quietly and smoothly adjust your patient into a variety of positions. Trendelenberg and reverse Trendelenberg can be accomplished with the push of one button.
- Expandable width (39 to 48 inches) and length (86 to 92 inches). Fits through standard hospital doorways, and accepts oversized mattresses.
- Swing-away trapeze allows patients to reposition themselves yet swings away to eliminate confinement or allow the easy access of lifting equipment.
- "Helping Hand"™ side-egress bars allow patients to help themselves into and out of bed, promoting general sense of independence.
- Solid polypropylene headboard & footboard designed for strength and easy integration with the Barton® *PTS™*.
- Built-in digital scale allows easy patient weight measurement without moving patient out of bed.

Barton® I-1000 Chair

- The *I-1000* chair is designed to recline into a stretcher position to easily facilitate a lateral transfer, using the innovative *PTS* transfer system or overhead lift.
- The power-assisted tilt and recline allow for easy and frequent repositioning of patients, allowing them to tolerate sitting for a longer time.
- The 1000-lb. weight capacity and 31" seat width allow the chair to fit a variety of patients, while still being able to fit through standard doorways.

Barton® Ready Lift 700

- The 700-lb. weight capacity and long boom and 4-point carry bar allow for the lifting or positioning of large patients.
- The low profile base is quickly adjusted and designed to access the most difficult areas.
- Barton® slings are available in multi or single patient use models and are uniquely designed to maximize fit, comfort, and ease of application.

Barton® Ceiling Track 1000

- The 1000-lb. weight capacity is designed into a single, easily used motor.
- The 4-point carry bar uses the same slings as our line of floor lifts.

Barton® Shower Trolley

- Height adjustable and portable, the XXX-capacity shower trolley allows for easy bedside bathing, helping maintain patient dignity.
- Flip down rails and tilt accommodate lateral transfers. Tilt allows for easy drainage into toilet.

The Barton® Positioning and Transfer System



I-1000 with accompanying Barton® *PTS™*.




1 Place chair at bed side in horizontal position.



2 Attach sheet to the Barton® *PTS™*.




3 Transfer patient by turning handle.




4 Incline chair and you're ready to go!

Process showing patient transfer. To transfer from the chair to the bed, reverse the process.



\ upright position
 — prone position
 / Trendelenburg position

B-1000 TransBed™ features



A Integrated digital scale for no-stress weight measurement.



B Solid polypropylene headboard and footboard has easy-grip holds, and includes built-in *PTS™* holder.



C "Helping Hand"™ side egress bar for patient safety and convenience.

D Swing-away overhead trapeze can be put to side when not in use.

Call Barton® or your local dealer for information.



(877) 8-BARTON

(877) 822-7866

"Your Safest Move"

BARTON
MEDICAL CORPORATION
www.bartonmedical.com

Technology that puts people first.

