



## Merlin Wand Version 2.00 Release Notes

This is CATC Merlin's Wand Software Version 2.0. You can find a summary of previous releases in Previous Releases section.

The application was updated to include new functionalities, first introduced in the BTTrainer exerciser, like the Device List archive.

New profiles are now supported by the Merlin's Wand: Serial port profile, File Transfer, Headset, Headset Gateway, LAN Access profile, Dial-up profile, FAX profile, PAN profile (GN and NAP), HCRP profile, Intercom profile (in Script Manager and CmdGen).

Support for selecting which analyzer to control (the Merlin or the Merlin Mobile) was added.

### **Automation:**

- The Merlin's Wand application now includes automation API (DCOM) allowing external client applications to control it. Please refer to the new "Automation API for CATC Bluetooth Analyzers Reference Manual" for further information.
- The Merlin's Wand can be run in minimized mode, so it would not visually interfere in case an automation client controls it.
- Sample files Installation: To install the sample files, please select the 'Custom' installation option and select the 'Automation' group. The samples are suitable for the Merlin's Wand analyzer. Please refer to the Automation manual for additional information.
- To run the sample files, you should make sure that the Merlin's Wand software version 2.0 is already installed.
- For editing the HTML sample code you can use a text editor or other HTML-editing tools.
- The library file name to be used for exporting the API is BTEexerciser.tlb.