

Release Notes PETracer Protocol Analyzers – 5.00

What's New

- PETracer Edge support
- PCI Express specification 2.0 support at Gen1 speeds.
- New Trace File format that provides smaller files and faster operation. Old traces can be automatically converted to the new format on open.
- "Simple Mode" is added to Recording Options for quick triggering and filtering setup
- More Errors are tracked in the Error summary
- Improved display of Ordered Sets
- New BusEngine for PETracer ML fixes the problem with the analyzer not always being able to lock to the data after power cycling of the DUT
- Bug fixes