

“PMM Emission Suite” Software History 2.16

Here are reported only the last 3 releases, the complete release history is available on request

VER.	DATE	APPLIED MODIFICATIONS
2.16	17/02/14	<ul style="list-style-type: none"> – Manual mode: when this modality is chosen the Pre selector is set as the last session – Modified conversion of TXT file to create a Correction Factor, now the software is able to manage more TXT file in a different format – Sweep Archive mode: added, in the save image option, the visualization of the selected sub measure – Modified report creation, now the "margin" parameter is able to change from -999 to +999 instead from -20 to +20
2.15	11/12/13	<ul style="list-style-type: none"> – Added the management of 7010/00, 7010/01, 7010/02 – Analyzer node: when is read a hold time > respect of what set up, the software memorizes the original hold time and visualized the hold time read – Modified management of Reference Level – End of Scale value so that they remain the same in every mode (Sweep Archive, Sweep, Analyzer, Manual) – Modified Reference level set up to have no differences between software and instrument.
2.14	17/10/13	<ul style="list-style-type: none"> – Software modification: now it is possible to interrupt a measure with a high number points made with USB port – Improved the procedure to use the COMM port to avoid communication problems. – Minor graphic improvements – Bugs correction