


- SA Shaper Amplifier
- PF Thyatron Pulse Former
- PA Pulse Amplifier
- T Trigger Tube
- D Delay
- DC A DC Amplifier


SINYDER

Note: Leads labelled with Command Symbols are connected to TTS or CTS output as per the particular Command.  
Nos. in paren. identify misc. F.F.'s


**start clock**  
**4100** Start Clock Unit




**WAIT-RESUME**  
**4200** WAIT-RESUME UNIT




**CLOCK**  
**4600** Clock generator




**STOP CLOCK PULSED RELAY**  
**4400**



**RELEASE LOAD-RUN SELECTION**



**STOP RELAY CONTROL**  
**4500**



MASTER CLEAR NEG  
RAW CLOCK PULSES  
STOP CLOCK

Stop after Command exec. 2-E 07 08 09

INCH BY COMM STOP 2-F 07 08 09

PD 760