

Please make the following alterations to the User's Manual IM 04L41B01-01E (see underlined text).

## ■ Page 12-1 “Measurement range and measurable range”

The measurable range for the type N thermocouple;

0.0 to 1300.0°C, on models with release number 2 or earlier.

-270.0 to 1300.0°C, on models with release number 3 and later.

## ■ Pages App-1 and App-2 “File Size”

### Size of Information Other Than the Sampled Data

Item	Size [Bytes]
Release number 3 expansion information	<u>1,696 + 80xN</u> <u>Expansion information is always stored.</u>

#### Example 1:

$216 + (88 \times 36 + 32) + 992 + (104 \times 50 + 32) + 832 + (80 + 32 + 36 \times 8 + 16 + 2) + (24 + 8) + \underline{4,576}$

= 15,498 bytes

#### Example 4:

From examples one and two, we obtain 15,498 + 357,120 = 372,618 bytes = 0.355MB

#### Example 5:

From examples one and three, we obtain 15,498 + 921,600 = 937,098 bytes = 0.893MB