

well below thresholds for physiological RyR1 activation. Contributions of MsC to the depolarisation were minor as judged from only slightly more negative resting potentials in the presence of GsMTx4 (-64 mV). We suggest that DHPR inhibition on RyR1 is modulated by MsC in mammalian skeletal muscle and is partly relieved in mdx muscle probably either due to some Ca^{2+} influx through aberrant MsC or direct interactions with the DHPR. A direct DHPR-RyR activation by depolarised membrane potentials in mdx fibers is more unlikely to be a consequence of osmotic membrane stress.

113-Plat

Fluvastatin Alters Both The Calcium Homeostasis And Cell Proliferation In Cultured Myotubes And The Calcium Release Events In Adult Muscle Fibers Of The Rat

Marta Fuzi¹, Zoltan Palitz¹, Laszlo Szabo², Janos Vincze¹, Peter Szentesi¹, Gyorgy Paragh¹, Pal Kertai¹, **Laszlo Csernoch¹**.

¹University of Debrecen, Debrecen, Hungary, ²Hungarian University of Transylvania, Targu Mures, Romania.

Skeletal muscle cells of newborn rats were cultured in the absence and presence of fluvastatin (0.1 μM). Resting levels of intracellular calcium concentration ($[\text{Ca}^{2+}]_i$) were determined from Fura-2 fluorescence while proliferation was assessed by counting the number of nuclei in myogenic cells. The presence of the drug did not alter $[\text{Ca}^{2+}]_i$ (130 ± 11 vs 121 ± 7 nM; $n=12$ vs 10 rats; mean $\pm$ SE; control and treated, respectively) but reduced the number of myogenic nuclei by 50% after 24 hours of treatment. To assess the chronic effects of fluvastatin on skeletal muscles of adult animals, female rats were kept on diets that either contained 62.5 mg/kg (daily intake 6 mg/kg body weight) fluvastatin or not. Animals were either fed with an otherwise normal chow or with a chow that induced an increase in blood cholesterol. Similarly to cultured cells, $[\text{Ca}^{2+}]_i$ of adult fibers was unaltered by the drug, however, a clear reduction of muscle mass was observed. Single fibers were enzymatically isolated from the *m. extensor digitorum communis*, permeabilized with Saponin and loaded with Fluo-4. Calcium release events (CRE) were captured using laser scanning confocal microscopy and analyzed with an automated computer program. Fluvastatin increased the frequency of CRE on both normo- and hypercholesterolaemic animals from 0.028 to 0.042 and 0.034 to 0.047 sarc-1s-1 (normo- and hypercholesterolaemic, respectively; $n=14$ vs 14 and 7 vs 15). While leaving the full width at half maximum unchanged the drug significantly increased the amplitude of sparks in both groups (0.405 ± 0.005 vs 0.436 ± 0.004 and 0.354 ± 0.004 vs 0.422 ± 0.005 ; $n=605$ vs 1429 and 741 vs 1052). This gave rise to an increased amount of released calcium in statin treated animals as reflected in the elevated signal mass.

Platform K: Ion Channels, Other

114-Plat

A Microfluidic Approach Enables Ligand Gated Ion Channel Recording from Cell Ensembles

Cristian Ionescu-Zanetti, Ian Spencer, Nianzhen Li, Jeff Jensen. Fluxion Biosciences, South San Francisco, CA, USA.

Electrophysiology remains the preferred method for characterizing ion channel function and kinetics. For screening, it is the most functionally relevant assay, and supersedes flux and fluorescent assays in terms of information content. Many of the medium to high throughput pharmaceutical screens are performed using the 'population patch' approach, which measures current from as many as 64 cells in parallel, eliminating a good deal of the cell-to-cell variability of single cell recordings. A major drawback of this method is the inability to exchange solutions during voltage clamp or apply multiple compounds to the same ensemble of cells.

Here we present the first data showing that, by integrating a microfluidic network design, a large number of cells under voltage clamp can be exposed to a compound within short time scale (50ms) in parallel. A comparison between competing compound injection designs will be presented along with validation data for a number of important ion channel targets. The fluid dynamics of the microfluidic networks were characterized by measuring the time domain response of channel activation and block. Another important aspect that we addressed is the rate of adsorption and desorption of compounds from the channel surfaces. Small fluorescent molecules were used to measure the surface properties as a function of molecular LogP values.

Continuous recording coupled with fast compound additions (50ms for a 20-cell ensemble) opens the way to ensemble recording for ligand gated ion channels, including fast desensitizing channels. We will present ion channel ensemble recording from cell lines expressing GABA-A, P2X3 and TRP-V1 that were obtained using this microfluidic approach.


115-Plat

Discovery Of Photochromic Ligands That Block Voltage-gated K⁺ Channels At The Internal TEA Binding Site

Matthew R. Banghart^{1,2}, Alexandre Mourot¹, Doris L. Fortin¹, Richard H. Kramer¹, Dirk Trauner^{1,3}.

¹University of California, Berkeley, Berkeley, CA, USA, ²Harvard Medical School, Boston, MA, USA, ³University of Munich, Munich, Germany.

We recently introduced an azobenzene derivative (AAQ) that affords photochemical control of endogenous potassium channels and cellular excitability. This electrophilic molecule was designed to function as a covalently attached photoswitchable affinity label for the external tetraethylammonium (TEA) binding site of potassium channels and was shown to have pronounced and long-lasting effects in various neural tissue preparations. Investigation into the mechanism of action revealed that AAQ is an open-channel blocker of the internal TEA binding site and that covalent modification of the channel is not required to achieve persistent photosensitivity. When applied externally to Shaker IR, AAQ causes fast inactivation in a voltage-dependent manner while exhibiting frequency-dependent block as well as reduced potency when external potassium is increased. Structure-activity studies of analogues lacking electrophilic groups identified blockers with increased potency and revealed structural features that allow the channel to discriminate between isomers. Although most analogues block in the trans form, one analogue (PrAQ) was found to be a cis-blocker. Importantly, the most potent trans-blocker BzAQ was established as an effective substitute for AAQ, affording photocontrol of endogenous K⁺ channels and action potential firing in dissociated hippocampal neurons.



116-Plat

Gaining control over membrane potential by light using Channelrhodopsin

Anton Rösler, Peter Hegemann, Matthias Prigge.

Experimentelle Biophysik, Humboldt Universität zu Berlin, Berlin, Germany.

Ion channels play a critical role in many pathological diseases. They can react on very different stimuli like for example ligands, heat, light and voltage. Especially voltage-gated channels are of enormous interest for pharmaceutical companies. Studying those channels in a high-throughput manner is a critical step towards developing specific drugs.

The light driven cation selective ionchannel ChR2 from *Chlamydomonas reinhardtii* was already successfully used to depolarize a variety of different neuron cells upon blue light stimulation.

We combined the light guided membrane potential change with calcium-sensitive dyes to analyze the activity of a voltage gated calcium channel. This is also possible because of a very low conductance of calcium ions by light activated ChR2.

Therefore, a stable HEK293 cell line expressing a leak potassium channel (mTrek) and voltage gated calcium channel (CaV3.2) was used. This cell line allows opening of the calcium channel by adding a low KCl concentration in order to activate the CaV3.2 (hyperpolarisation) followed by exchange with a high KCl concentration in order to open the CaV3.2 (depolarisation).

Upon transiently transfecting those cells with ChR2 we made them susceptible to light. In this combination we were able to depolarize the cells with blue light causing an opening of the calcium channel and hence an increased intracellular calcium concentration, which was monitored by fura2. This assay allows us to repeatedly open CaV3.2 simply by light.

Also fraction of open channels could be controlled by different light intensity. In a proof of principle experiment the dose-dependent inhibition of known inhibitor could be reproduced. This assay can be easily implemented into a multi-well reader assay and can be used for screening different voltage-gated calcium channels.