

of normetanephrine (10^{-4} M) had no effect on the neuronal accumulation but completely prevented the accumulation in the small non-neuronal cells.

These findings strongly indicate that not all the noradrenaline or α -methyl-noradrenaline present in the heart immediately after the perfusion with high concentrations of the amine (that is, Uptake₂) is located in the adrenergic nerves, but that it is also accumulated extraneuronally, most evidently in the muscle cells and in small non-neuronal cells.

REFERENCES

- IVERSEN, L. L. (1965). The uptake of catecholamines at high perfusion concentrations in the isolated rat heart: a novel catecholamine uptake process. *Br. J. Pharmac. Chemother.*, **25**, 18-33.
MALMFORS, T. (1965). Studies on adrenergic nerves. *Acta physiol. scand.*, **64**, Suppl. 248, 1-93.

A new programmed interval timer for use with automatic assay apparatus

W. T. S. AUSTIN and G. HOOPER (introduced by R. B. BARLOW), *E.P.I. Limited, Edinburgh, and Department of Pharmacology, University of Edinburgh*

The apparatus is designed to give accurate programmed intervals from 1 to 1,000 sec. An electronic pulse generator provides the basic frequency, which is divided by dekatron tubes into tens and units to give the required time intervals. The cathodes of the dekatrons are brought to switches, two for each channel, and the circuit (called a coincidence circuit) is so arranged that a signal from the pulse generator is made available to both of the switches when the required ten and unit selected is reached.

Each channel has a bistable circuit which controls its own output relay. If a channel is live, the coincidence pulse from the generator will switch it off, and automatically the circuit will then return the dekatron to the zero position and the new count will start.

The front panel of the apparatus has plugs so that a programme of any sequence can be selected. The "off" pulse of any one bistable circuit is used to trigger the "on" position of the next bistable circuit, which is switched off by its selected time pulse, and so on. "Hold" and "cancel" controls are fitted on each channel. An important feature is the time multiplier which allows the apparatus to be infinitely variable up to a thousand seconds in each individual channel and therefore for any number of operations the machine is required to perform. Spare bistable circuits are available so that sequences of operations may be repeated within a programme.

The machine can be used to produce an infinitely variable programme. It can be used to control, say, a 90 sec cycle for use with acetylcholine-like compounds on the guinea-pig ileum or a frog rectus preparation for the assay of nicotine-like compounds in which the time cycle is 30 min.

The apparatus will drive standard G.P.O. relays and can also be used to control solenoid valves.

The use of intracerebral pyretogenins in testing for antipyretic activity in conscious mice

C. H. CASHIN and CHRISTINE E. HEADING, *Lilly Research Centre Ltd., Erl Wood Manor, Windlesham, Surrey*

The injection of pyretogenins into the cerebral ventricles of conscious mice produces a pyresis of short duration, which can be antagonized by prior administration of anti-pyretic drugs at doses approximating those used clinically.

Temperatures are measured at 15 min intervals for 1 hr with a thermistor (Standard Telephones and Cables Ltd., Type F15) inserted 1 cm into the rectum of each mouse.