High Sensitivity Optical Resolution of DL-Amino-acids by Gas Chromatography

By B. HALPERN and J. W. WESTLEY

(Department of Genetics, Stanford Medical School, Palo Alto, California)

RECENT advances in synthetic peptide chemistry impel attention to optical purity of the aminoacids, which should have a steric impurity of less than 0.1% for use in peptide synthesis.¹ Optical rotation does not usually reveal this amount of contamination. Furthermore, there is no theoretical value for an optical rotation, but only the value on which most observers agree. Optically specific amino-acid oxidases and decarboxylases are more sensitive,² but are usually not chemically specific.

The advantages of gas-liquid chromatography (g.l.c.) for the separation of diastereoisomeric TFA-dipeptide esters have already been demonstrated by Weygand,³ and we have recently

¹ J. P. Greenstein and M. Winitz, "Chemistry of the Amino-Acids," Wiley, New York, 1961, p. 946.

² A. Meister, L. Levintow, R. B. Kingsley, and J. P. Greenstein, *J. Biol. Chem.*, 1951, **192**, 535. ³ F. Weygand, Angew. Chem. Internat. Edn., 1963, **2**, 183.

NUMBER 12, 1965

applied this technique to show the asymmetry of some neutral amino-acids.⁴ While this procedure gives excellent results with some neutral aminoacids, many TFA-dipeptide esters have such low volatility that g.l.c. analysis results in very long retention times and subsequently poor separation of the diastereoisomeric pairs.

In typical assays (Table 1), the amino-acid was esterified with thionyl chloride-methanol⁶ and the excess of reagent and solvent removed. An excess of $L-\alpha$ -chloroisovaleryl chloride in an inert solvent was added to the residue and the suspension cooled and neutralized with triethylamine. After being washed (H₂O) and dried (Na₂SO₄), the

TABLE 1. Gas-chromatographic separation of diastereoisomeric N-a-chloroisovalerylamino-acid methyl esters*

TUDDD 1			<i>c cop a</i> . <i>a</i> . <i>cop a</i> .			<i></i>	
			Separation	N_2 Flow	R.t. [‡] of diasteroisomer		Ratio of R.t.‡ LL/LD
Amino-acid		Column†	temp. (°C)	(ml./min.)	LD	LL	
Alanine		в	1 61	28	8.05	9.15	1.13
Valine		ſΑ	161	28	5.9	6.6	1.12
	••	ίв	161	28	8.1	9.9	$1 \cdot 22$
Leucine		ſВ	161	28	11.85	13.65	1.12
	••	٦c	175	56	11-1	12.7	1.14
Isoleucine		Ъ	161	28	10.8	12.8	1.18
Proline		ſB	185	28	16-1	18.9	1.17
	•••	íc	200	56	18.0	$21 \cdot 2$	1.17
Threonine		A	145	28	18.0	19.9	1.11
Phenylalanin	.e	Α	185	28	15.4	16.5	1.08

* G.l.c. analyses were carried out on a Wilkens 600C aerograph equipped with a flame ionization detector.

A 5-foot steel column 5% SE 30 on Chromosorb W.
B 5-foot steel column 5% FFAP on Chromosorb W.

С 5-foot steel column 20% DEGS on Chromosorb W.

 \ddagger R.t. = Retention time.

TABLE 2. Gas-chromatographic separation of diastereoisomeric N-chloroacyl-valine methyl esters

			Separation	N ₂ Flow	R.t. of diastereoisomer		Ratio of R.t.
Chloro-acid from:		Column	temp. (°C.)	(ml./min.)	DL	LL	ll/dl
Alanine	••	Α	140	28	5.9	6.25	1.06
		в	161	28	5.75	6.5	1.13
Valine	••	Α	140	28	11.85	13.5	1.14
		в	161	28	8.1	9.9	1.22
Isovaline	••	Α	140	28	9.3	9.8	1.05
		в	161	28	4.6	4.95	1.08
Leucine	••	Α	140	28	17.05	18.8	1.10
		в	161	28	10.15	12.0	1.18
Isoleucine		Α	140	28	17.9	20.6	1.12
		в	161	28	10.35	12.9	1.25
Alloisoleucine	••	Α	140	28	18.9	21.7	1.12
		в	161	28	10.5	12.85	1.22
Ornithine	••	Α	175	28	10.4	11.4	1.10
		в	214	28	8.4	9.9	1.17
Lysine		Α	175	28	15.7	17.5	1.10
•		в	214	28	11.45	12.8	1.12

We have now found that α -halogeno-acyl chlorides are excellent reagents for quantitative g.l.c. analysis of amino-acids, both with respect to chemical structure and resolution of stereoisomers. In particular, $L-\alpha$ -chloroisovaleryl chloride can be used to analyze many DL-amino-acids and in most cases quantitative separation of the diastereoisomeric pairs is readily achieved on packed columns only 5 feet long. The reagent can be prepared from L-valine by Renard's procedure⁵ followed by treatment with thionyl chloride; the product is 100% sterically pure.

solution was injected into the gas chromatograph.

The amino-acid may also be converted into its α -chloro-derivative and then resolved with L-valine methyl ester (Table 2).

Finally, α -halogeno-acids are of special interest in the studies concerned with Walden inversions in the amino-acid series. Reliable and sensitive methods of optical analysis therefore suggest a re-investigation of the work of Fischer and others.7,8

(Received, May 3rd, 1965.)

⁴ B. Halpern and J. W. Westley, Biochem. Biophys. Res. Comm., 1965, 19, 361.

⁵ M. Renard, Bull. Soc. Chim. biol., 1946, 28, 497.

- ⁶ M. Brenner and W. Huber, Helv. Chim. Acta, 1953, 36, 1109.
- ⁷ E. Fischer et al., Annalen, 1905, 340, 123; 1907, 357, 1; 1911, 303, 337; Ber., 1908, 41, 889, 2891; 1909, 42, 1219.
- ⁸ A. Neuberger, Adv. Protein Chem., 1948, 4, 297.