## The Reaction of Imidazole with Dichlorocarbene in the Vapour Phase: a New Route to the Pyrimidine Ring System

By R. E. Busby, M. IQBAL, J. PARRICK, and C. J. G. SHAW\*

(Department of Chemistry, Brunel University, Woodlands Avenue, London, W.3)

Summary Reaction of imidazole with dichlorocarbene at 550° gives a mixture of 5-chloropyrimidine and chloropyrazine (10:1).

RECENTLY we reported the preparation of 2- and 3-chloropyridine in high yield (86%) by the reaction of pyrrole with dichlorocarbene generated from CHCl<sub>3</sub> in the vapour phase.<sup>1</sup> Using essentially the same conditions, we have now studied the reaction between imidazole and dichlorocarbene. The best yields were obtained by introducing a solution of imidazole in CHCl<sub>3</sub> (5 M) into a stream of N<sub>2</sub>, pre-heating the reactants at 350°, and then passing them through a vertically mounted Pyrex glass tube maintained at 550°. Unchanged imidazole (30%) was recovered and re-cycled, and a 33% yield† of ring-expanded products obtained. These were separated by preparative g.l.c. [7 ft.  $\times \frac{3}{8}$  in. glass column, 20% w/w Carbowax 20 M on Diatomite CAW (60-72 mesh), 100°] into 5-chloropyrimidine (30%) and chloropyrazine (3%).‡

Addition of dichlorocarbene to the C=C or C=N double bonds of imidazole followed by ring expansion would give 5-chloropyrimidine and chloropyrazine, respectively. It appears that the C=C double bond is more susceptible to attack by dichlorocarbene than the C=N double bond. This trend is in keeping with the work of Jones and Rees,<sup>2</sup> who treated 2,4,5-trimethylimidazole with dichlorocarbene, generated either under neutral conditions from the thermal decarboxylation of sodium trichloroacetate, or under basic conditions from CHCl<sub>3</sub> and EtONa. In each case only a single ring-expanded product, 5-chloro-2,4,6-trimethylpyrimidine (best yield: 12.4%), was obtained. This must arise from attack of dichlorocarbene on the 4,5-C=C double bond. No product arising from addition to the 2,3-C=N double bond could be detected.

We are investigating this type of selective reactivity of carbenes in other systems (benzimidazole, methylimidazoles, etc.), and preliminary results confirm our findings with imidazole.

5-Chloropyrimidine was reduced at atm. press. using a Pd/BaSO<sub>4</sub> catalyst in the presence of MgO.<sup>3</sup> The dehalogenated compound was separated (85%) and shown to be pyrimidine. It is thus possible to convert imidazole into pyrimidine in 28% overall yield by this method.

We thank the Science Research Council for financial support, Dr. G. W. Cheeseman for a sample of chloropyrazine, and Mr. R. J. Langston for technical assistance.

(Received, September 29th, 1969; Com. 1467.)

- Based on imidazole changed.
- f Structures confirmed by elemental analysis, i.r., u.v., n.m.r., and mass spectra.
- <sup>1</sup> F. S. Baker, R. E. Busby, M. Iqbal, J. Parrick, and C. J. G. Shaw, *Chem. and Ind.*, 1969, 1344.

  <sup>2</sup> R. L. Jones and C. W. Rees, *J. Chem. Soc.* (C), 1969, in the press. We thank Professor Rees for permission to quote his results prior to publication. <sup>3</sup> B. Lythgoe and L. S. Raynor, J. Chem. Soc., 1951, 2323.