

Metastable Phase Equilibrium in the Quaternary System NaCl + KCl + Na₂CO₃ + K₂CO₃ + H₂O at 273 K

Ying Zeng,* Junmei Liu, and Zhiyuan Zheng

College of Materials and Chemistry & Chemical Engineering, Chengdu University of Technology, Chengdu, 610059, P. R. China

The metastable phase equilibrium in the quaternary system NaCl + KCl + Na₂CO₃ + K₂CO₃ + H₂O was studied at 273 K using an isothermal evaporation method. The solubility and the density of the equilibrated solution were determined. On the basis of the solubility data, the metastable phase diagram and the water diagram were plotted. The phase diagram consists of three invariant points, seven univariant curves, and five metastable crystallization fields. Among the three invariant points, points E₁ and E₃ are a commensurate invariant point type, while point E₂ is an incommensurate invariant point type. The three cosaturated salts for points E₁, E₂, and E₃ are NaCl + KCl + Na₂CO₃·10H₂O, KCl + Na₂CO₃·10H₂O + KNaCO₃·6H₂O, and KCl + K₂CO₃·3/2H₂O + KNaCO₃·6H₂O, respectively. The five metastable crystallization fields correspond to the single salts potassium chloride (KCl), sodium chloride (NaCl), sodium carbonate decahydrate (Na₂CO₃·10H₂O), potassium carbonate sesquihydrate (K₂CO₃·3/2H₂O), and a hexahydrate double salt of potassium and sodium carbonate (KNaCO₃·6H₂O). Sodium carbonate decahydrate (Na₂CO₃·10H₂O) is the only crystallization form of sodium carbonate at 273 K in this system; the monohydrate (Na₂CO₃·H₂O) and heptahydrate (Na₂CO₃·7H₂O) are not formed.

Introduction

Salt lakes are widely distributed in the west of China, especially in the area of the Qinghai-Xizang (Tibet) Plateau. Salt lakes of the Qinghai-Xizang Plateau are of a carbonate brine type or sulfate brine type. Zabuye Salt Lake, located in Tibet, is a typical carbonate brine type famous for its abundance of lithium, potassium, and boron resources and also with a scarce magnesium resource.¹ The lowest concentration ratio of magnesium to lithium makes the comprehensive utilization technology of Zabuye Lake quite different from other salt lakes in the Qinghai-Xizang Plateau, such as Dongtai Lake, Xitai Lake, and so forth.

The composition of Zabuye Salt Lake brine falls under the complex system which includes lithium, sodium, potassium, carbonate, sulfate, chloride, and borate. The climate conditions in the region of Zabuye Lake are generally windy and arid, with an average temperature of about 273 K.¹ To economically exploit the brines and mineral resources, it is important to adequately adopt local wind and sun energy; thus a technique like solar ponds is widely used. Earlier researchers found that in the evaporation process of brine, the equilibria among the salts are always metastable;^{2,3} therefore, metastable phase equilibrium research is essential to predict the crystallization path of evaporation of the salt lake brine. Although a number of papers describing the metastable equilibrium of the seawater system Na + K + Mg + Cl + SO₄ + H₂O have been reported,^{4–7} little work has been done on the carbonate brine type, and data for the quaternary system NaCl + KCl + Na₂CO₃ + K₂CO₃ + H₂O at 273 K are not reported in the literature yet.

The quaternary system NaCl + KCl + Na₂CO₃ + K₂CO₃ + H₂O is one of the most important and basic subsystems of carbonate type brines. In this paper, the metastable equilibria

of the quaternary system NaCl + KCl + Na₂CO₃ + K₂CO₃ + H₂O at 273 K are presented. The solubilities and densities of the equilibrium solution of this system were measured.

Experimental Section

Apparatus and Reagents. A SHH-250 type thermostatic evaporator made by the Chongqing Inborn Instrument Corporation, China, was used for the metastable phase equilibrium experiments. The equipment has a temperature-controlling apparatus with a blower accessory to control the temperature and the evaporation quantity of the evaporating system. The evaporation rate is (4.0 to 5.5) mm·d⁻¹, and the temperature precision is ± 0.1 K.

A Siemens D500 X-ray diffractometer with Ni-filtered Cu K α radiation was used to analyze the crystalloid form of the solid phases. The operating conditions are 35 kV and 25 mA. An atomic absorption spectrometer (type WYD-YII, Beijing Geology Instrument Corporation) was employed for the determination of the sodium ion concentration in solution.

The chemicals used in this work were of analytical purity grade and obtained from the Chengdu Chemical Reagent Ltd. or the Shanghai Chemical Reagent Ltd. They were sodium chloride (NaCl, 99.5 % (by mass)), potassium chloride (KCl, 99.5 % (by mass)), potassium carbonate (K₂CO₃, 99.0 % (by mass)), and sodium carbonate (Na₂CO₃, 99.0 % (by mass)). Potassium carbonate and sodium carbonate were heated at 553 K to eliminate bicarbonate before use. Doubly deionized water was obtained from a Millipore water system with an electrical conductivity less than 1·10⁻⁴ S·m⁻¹ and pH = 6.6. For each experiment, the required amounts of reagents were dissolved into enough deionized water to produce the experimental solutions.

Experimental Method. The isothermal evaporation method was employed in this study. According to the phase equilibrium composition, the appropriate quantity of salts were mixed

* Corresponding author. Tel.: +86-28-84079016. Fax: +86-28-84079074. E-mail address: zengy@gmail.com.

Table 1. Experimental Solubility and Density Values of Equilibrated Solutions in the Quaternary System NaCl + KCl + Na₂CO₃ + K₂CO₃ + H₂O at 273 K

no.	composition of liquid phase, $w(B)^a \cdot 100$				Janěcke index of the dry salt, $J(B)/[\text{mol}/100 \text{ mol dry salt}]$			density $\text{g} \cdot \text{cm}^{-3}$	equilibria solid phase
	$w(\text{CO}_3^{2-})$	$w(\text{Cl}^-)$	$w(\text{K}^+)$	$w(\text{Na}^+)$	$J(\text{CO}_3^{2-})$	$J(\text{K}^+)$	$J(\text{H}_2\text{O})$		
1, C	19.17	0.00	17.05	4.63	100.0	68.47	812.3	1.4223	Na ₂ CO ₃ ·10H ₂ O + KNaCO ₃ ·6H ₂ O
2	21.72	0.86	20.07	5.40	96.76	68.62	770.9	1.4388	Na ₂ CO ₃ ·10H ₂ O + KNaCO ₃ ·6H ₂ O
3	21.94	0.86	20.74	5.18	96.77	70.19	753.4	1.4018	Na ₂ CO ₃ ·10H ₂ O + KNaCO ₃ ·6H ₂ O
4	19.77	0.89	19.98	3.98	96.31	74.69	898.5	1.4192	Na ₂ CO ₃ ·10H ₂ O + KNaCO ₃ ·6H ₂ O
5, E ₂	11.36	0.44	12.99	1.33	96.83	85.15	2097	1.3393	KCl + Na ₂ CO ₃ ·10H ₂ O + KNaCO ₃ ·6H ₂ O
6	55.36	13.00	12.12	8.30	87.36	46.26	2897	1.3480	KCl + Na ₂ CO ₃ ·10H ₂ O
7	57.57	13.80	7.40	8.80	83.16	33.16	2789	1.3500	KCl + Na ₂ CO ₃ ·10H ₂ O
8	21.03	13.30	1.28	10.9	65.16	6.45	2795	1.2330	KCl + Na ₂ CO ₃ ·10H ₂ O
9	7.59	14.50	0.66	12.10	38.24	3.10	2210	1.2301	KCl + Na ₂ CO ₃ ·10H ₂ O
10	4.36	15.10	0.62	12.80	25.47	2.79	1857	1.2248	KCl + Na ₂ CO ₃ ·10H ₂ O
11	3.46	15.80	0.64	12.30	20.57	2.96	1785	1.2241	KCl + Na ₂ CO ₃ ·10H ₂ O
14	2.44	14.83	0.42	11.23	16.29	2.18	1582	1.2253	KCl + Na ₂ CO ₃ ·10H ₂ O
13	2.26	15.02	0.40	11.23	15.11	2.06	1585	1.2153	KCl + Na ₂ CO ₃ ·10H ₂ O
12	2.00	15.10	0.39	11.08	13.55	2.04	1613	1.2186	KCl + Na ₂ CO ₃ ·10H ₂ O
15, A	1.94	14.95	0.00	11.17	13.32	0.00	1646	1.2153	NaCl + Na ₂ CO ₃ ·10H ₂ O
16	1.79	14.08	0.32	10.31	13.08	1.82	1791	1.2053	NaCl + Na ₂ CO ₃ ·10H ₂ O
17	1.95	15.18	0.37	11.11	13.20	1.91	1611	1.2184	NaCl + Na ₂ CO ₃ ·10H ₂ O
18, E ₁	1.93	15.21	0.44	11.08	13.06	2.27	1609	1.2231	NaCl + Na ₂ CO ₃ ·10H ₂ O + KCl
19	12.04	14.10	0.74	11.83	50.27	3.56	2523	1.2323	NaCl + Na ₂ CO ₃ ·10H ₂ O
20, B	0.00	15.64	3.68	7.96	0.00	21.42	1758	1.2173	NaCl + KCl
21	0.16	17.22	4.19	8.81	1.09	21.90	1577	1.2159	NaCl + KCl
22	0.17	17.43	4.10	9.00	1.14	21.17	1550	1.2233	NaCl + KCl
23	0.73	23.5	2.47	8.10	3.56	15.25	1555	1.2261	NaCl + KCl
24	1.11	21.12	2.15	8.00	5.87	13.68	1563	1.2302	NaCl + KCl
25	1.41	19.55	1.51	7.81	7.89	10.21	1576	1.2354	NaCl + KCl
26	1.66	18.56	0.74	7.55	9.59	5.47	1586	1.2360	NaCl + KCl
27, D	22.47	0.00	26.34	1.39	100.0	91.79	592.4	1.5361	K ₂ CO ₃ ·3/2H ₂ O + KNaCO ₃ ·6H ₂ O
28	24.87	0.32	28.73	1.98	98.93	89.72	575.6	1.4804	K ₂ CO ₃ ·3/2H ₂ O + KNaCO ₃ ·6H ₂ O
29	24.50	0.30	28.91	2.08	98.96	89.06	597.4	1.5015	K ₂ CO ₃ ·3/2H ₂ O + KNaCO ₃ ·6H ₂ O
30, E ₃	24.96	0.30	28.58	2.11	99.00	89.09	573.1	1.5072	K ₂ CO ₃ ·3/2H ₂ O + KNaCO ₃ ·6H ₂ O + KCl
31	24.69	0.31	28.18	2.32	98.94	87.88	588.6	1.4962	KCl + KNaCO ₃ ·6H ₂ O
32	24.59	0.33	26.54	2.48	98.86	86.97	594.7	1.5040	KCl + KNaCO ₃ ·6H ₂ O
33, F	20.60	0.39	27.21	0.00	98.32	100.0	712.1	1.5199	KCl + K ₂ CO ₃ ·3/2H ₂ O
34	22.18	0.40	27.54	1.02	98.80	94.09	723.1	1.4941	KCl + K ₂ CO ₃ ·3/2H ₂ O
35	22.40	0.38	26.82	1.60	98.59	90.82	716.2	1.5198	KCl + K ₂ CO ₃ ·3/2H ₂ O

^a $w(B)$ is the mass fraction of B.

together as a series of artificial synthesized brines and loaded into clean polyethylene containers (24 cm long, 14 cm wide, and 6 cm high), and then the containers were placed in a thermostatic evaporator for isothermal evaporation. The temperature of the solution inside the evaporator was controlled to (273 ± 0.1) K and measured by a J thermocouple with an operating range of $(258.15 \text{ to } 373.15)$ K. The crystal behavior of the solid phase was observed periodically. When enough new solids appeared in the evaporating bottles, the liquid and solid phases were separated by filtration. The obtained wet solids were divided into two parts. One was analyzed by chemical methods to obtain the composition of the salts. The other was dried at 273 K, pestled into a powder, and then identified with X-ray diffraction to determine the crystalloid form of the solid phase. Meanwhile, a 5.0 mL sample of the clarified solution was taken from the liquid phase and diluted to a 50 mL final volume in a volumetric flask filled with deionized water to analyze the liquid-phase components. Another 5.0 mL sample of the clarified solution was taken to measure the density. The remainder of the solution continued to be evaporated to reach the next measuring point. The same procedure was repeated until the solution was fully evaporated.

The densities of solution determined in this study were used for the mass fraction calculation of the liquid components. The specific gravity bottle method with a correction for the air floating force was used. The precision of the density measurements was $0.0001 \text{ g} \cdot \text{cm}^{-3}$. The sample was obtained at 273 K and measured at room temperature. The excursion caused by

the temperature difference was calibrated by using a multipoint temperature-revision method.⁸

Analytical Method. The composition of K^+ in the liquids and their corresponding wet solid phases were analyzed by using a sodium tetraphenylborate (STPB)/hexadecyl trimethyl ammonium bromide (CTAB) back-titration method with the existence of propanetriol, with a precision within 0.5 % (by mass).⁹ The sodium ion concentration was determined by atomic absorption spectrometry, with a precision of 0.06 % (by mass). The chlorine ion concentration was determined by titration with a mercury nitrate standard solution in the presence of a mixed indicator of diphenylcarbazone and bromophenol blue with a precision of 0.3 % (by mass). The carbonate ion concentration was determined by an acid–base neutralization titration with a precision within 0.5 % (by mass).⁹

Results and Discussion

The determined values of solubilities and densities of the equilibrated solution in the quaternary system NaCl + KCl + Na₂CO₃ + K₂CO₃ + H₂O at 273 K are presented in Table 1. The ion concentration values in the equilibrated solution are expressed in mass fraction $w(B)$. The composition of dry salts are expressed in Janěcke index $J(B)$, with $J(\text{Na}_2^{2+}) + J(\text{K}_2^{2+}) = 100 \text{ mol}$ or $J(\text{CO}_3^{2-}) + J(\text{SO}_4^{2-}) = 100 \text{ mol}$. On the basis of the experimental data in Table 1, the metastable equilibrium phase diagram of the system at 273 K is plotted in Figure 1.




Figure 1. Metastable equilibrium phase diagram of the quaternary system $\text{NaCl} + \text{KCl} + \text{Na}_2\text{CO}_3 + \text{K}_2\text{CO}_3 + \text{H}_2\text{O}$ at 273 K.

The metastable phase diagram of the system in Figure 1 consists of five crystallization fields corresponding to potassium chloride (KCl), sodium chloride (NaCl), sodium carbonate decahydrate ($\text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O}$), potassium carbonate sesquihydrate ($\text{K}_2\text{CO}_3 \cdot 3/2\text{H}_2\text{O}$), and a hexahydrate double salt of potassium and sodium carbonate ($\text{KNaCO}_3 \cdot 6\text{H}_2\text{O}$). There are three invariant points in this system, noted as E_1 , E_2 , and E_3 , respectively. Point E_1 is a commensurate invariant point, which is saturated with three salts $\text{NaCl} + \text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O} + \text{KCl}$, and the mass fraction composition of the corresponding liquid phase is $w(\text{CO}_3^{2-}) = 1.93\%$, $w(\text{Cl}^-) = 15.21\%$, $w(\text{K}^+) = 0.44\%$, and $w(\text{Na}^+) = 11.08\%$. Point E_2 is an incommensurate invariant point, which is saturated with three salts, $\text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O} + \text{KCl} + \text{KNaCO}_3 \cdot 6\text{H}_2\text{O}$, and the mass fraction composition of the corresponding liquid phase is $w(\text{CO}_3^{2-}) = 11.36\%$, $w(\text{Cl}^-) = 0.44\%$, $w(\text{K}^+) = 12.99\%$, and $w(\text{Na}^+) = 1.33\%$. Point E_3 is also a commensurate invariant point. The cosaturated three salts for invariant point E_3 are $\text{KCl} + \text{K}_2\text{CO}_3 \cdot 3/2\text{H}_2\text{O} + \text{KNaCO}_3 \cdot 6\text{H}_2\text{O}$, and the mass fraction composition of the corresponding liquid phase is $w(\text{CO}_3^{2-}) = 24.96\%$, $w(\text{Cl}^-) = 0.30\%$, $w(\text{K}^+) = 28.58\%$, and $w(\text{Na}^+) = 2.11\%$. There are seven isothermal evaporation curves corresponding to curves E_1A , BE_1 , E_1E_2 , E_2C , E_2E_3 , E_3D , and E_3F , indicating the cosaturation of two salts. Points A, B, C, D, and F, cosaturated with two salts, are the invariant points of the ternary subsystems of the quaternary system $\text{NaCl} + \text{KCl} + \text{Na}_2\text{CO}_3 + \text{K}_2\text{CO}_3 + \text{H}_2\text{O}$, respectively. The double salt $\text{KNaCO}_3 \cdot 6\text{H}_2\text{O}$ is formed in the ternary subsystem $\text{Na}_2\text{CO}_3 + \text{K}_2\text{CO}_3 + \text{H}_2\text{O}$ and also formed in the quaternary system $\text{NaCl} + \text{KCl} + \text{Na}_2\text{CO}_3 + \text{K}_2\text{CO}_3 + \text{H}_2\text{O}$ at 273 K. The crystallization field of salt $\text{K}_2\text{CO}_3 \cdot 3/2\text{H}_2\text{O}$ is the smallest, while the crystallizing field of potassium chloride is the largest. These results indicate that potassium chloride is easy to saturate and crystallize from solution at 273 K.

Figure 2 is the relevant water-phase diagram of the system at 273 K. Figure 2 shows that the Janěcke indexes of water (mol/100 mol of dry salt) change with the Janěcke indexes of carbonate and reach the smallest value at the commensurate invariant point E_3 .

There are three hydrated salts of sodium carbonate, that is, $\text{Na}_2\text{CO}_3 \cdot \text{H}_2\text{O}$, $\text{Na}_2\text{CO}_3 \cdot 7\text{H}_2\text{O}$, and $\text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O}$. A lot of work has been done on the metastable or stable equilibria




Figure 2. Water-phase diagram of the quaternary system $\text{NaCl} + \text{KCl} + \text{Na}_2\text{CO}_3 + \text{K}_2\text{CO}_3 + \text{H}_2\text{O}$ at 273 K.

of the systems containing sodium carbonate from (273 to 318) K.^{10–13} According to previous studies on the stable equilibria,^{9–12} the crystallization forms of sodium carbonate are $\text{Na}_2\text{CO}_3 \cdot \text{H}_2\text{O}$ at 318 K, $\text{Na}_2\text{CO}_3 \cdot 7\text{H}_2\text{O}$ at 308 K, and $\text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O}$ at (273 to 298) K. However, Fang's study on the metastable equilibria at 298 K showed that the crystallization form of sodium carbonate is $\text{Na}_2\text{CO}_3 \cdot 7\text{H}_2\text{O}$ at 298 K.¹⁴ The transition temperature of $\text{Na}_2\text{CO}_3 \cdot 7\text{H}_2\text{O}$ and $\text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O}$ is 298 K. The crystal water in the sodium carbonate decreases along with an increase of temperature. No matter whether in the stable equilibria or in the metastable equilibria, the crystallization form of sodium carbonate is $\text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O}$ at 273 K. Our research result also supports this conclusion.

Conclusions

Metastable equilibria of the quaternary system $\text{NaCl} + \text{KCl} + \text{Na}_2\text{CO}_3 + \text{K}_2\text{CO}_3 + \text{H}_2\text{O}$ at 273 K were studied using an isothermal evaporation method. The solubilities and densities of the equilibrated solution were determined experimentally. According to the experimental data, the metastable equilibrium phase diagram and the water diagram of the system were plotted. The results show that there is a double salt of potassium and sodium carbonate ($\text{KNaCO}_3 \cdot 6\text{H}_2\text{O}$) formed in this quaternary system. $\text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O}$ is the only crystallization form of sodium carbonate at 273 K in this system.

Literature Cited

- Zheng, M. P.; Liu, W. G.; Xiang, J. *Qinghai-Xizhang Plateau Saline Lakes*; Beijing Science and Technology Press: Beijing, China, 1989; pp 192–270.
- Lepeshko, I. N.; Romaschova, N. N. Solubility in system $\text{Li}_2\text{SO}_4\text{-Na}_2\text{SO}_4\text{-MgSO}_4\text{-H}_2\text{O}$ at 298 K. *J. Inorg. Chem.* **1959**, *12*, 2812–2815.
- Gao, S. Y.; Song, P. S.; Xia, S. P.; Zheng, M. P. *Chemistry of Salt Lakes - New Salt Lake of Boron and Lithium*; Beijing Science and Technology Press: Beijing, China, 2007; pp 61–73.
- Shu, Y. G.; Li, J.; Jiang, C. F. Studies on the Metastable Phase Equilibrium of Na^+ , K^+ , $\text{Mg}^{2+}/\text{Cl}^-$, $\text{SO}_4^{2-}\text{-H}_2\text{O}$ Quinary System at 288 K. *J. Chem. Ind. Eng. (China, Chin. Ed.)* **1992**, *43*, 549–555.
- Jin, Z. M.; Zhou, H. N.; Wang, L. S. Studies on the Metastable Phase Equilibrium of Na^+ , K^+ , $\text{Mg}^{2+}/\text{Cl}^-$, $\text{SO}_4^{2-}\text{-H}_2\text{O}$ Quinary System at 288 K. *Chem. J. Chin. Univ.* **2002**, *23*, 690–694.
- Jin, Z. M.; Xiao, X. Z.; Liang, S. M. Study of the Metastable Equilibrium for Pentanary System of (Na^+ , K^+ , Mg^{2+}), (Cl^- , SO_4^{2-}), H_2O . *Acta Chim. Sin.* **1980**, *38*, 313–321.
- Jin, Z. M.; Zhou, H. N.; Wang, L. S. Studies on the Metastable Phase Equilibrium of (Na^+ , K^+ , Mg^{2+}), (Cl^- , SO_4^{2-}), H_2O Quinary System at 308.15 K. *Chem. J. Chin. Univ.* **2001**, *22*, 634–638.
- Xu, G. L. The Determination and Application of the Temperature Correction Coefficient for the Density of Vinyl Acetate. *Fujian Anal. Test.* **1999**, *3*, 1104–1108.
- Institute of Qinghai Salt-Lake of Chinese Academy of Sciences. *Analytical Methods of Brines and Salts*, 2nd ed.; Chinese Science Press: Beijing, 1984; pp 75–80.

- (10) Zeng, Y.; Yin, H. A.; Tang, M. L.; Feng, X. X. Studies on the Phase Equilibrium of Li^+ , Na^+ , $\text{K}^+/\text{CO}_3^{2-}$, Cl^- - H_2O Quinary System at 298.15 K. *Chem. J. Chin. Univ.* **2003**, *24*, 968–972.
- (11) Yin, H. A.; Hao, L. F.; Zeng, Y.; Tang, M. L. Studies on the Phase Equilibrium and Physico-chemical Properties of the Equilibrated Solution in the Quinary System Li^+ , $\text{Na}^+/\text{CO}_3^{2-}$, $\text{B}_4\text{O}_7^{2-}$, Cl^- - H_2O at 298 K. *J. Chem. Eng. Chin. Univ.* **2003**, *17*, 1–5.
- (12) Li, J.; Yang, Q.; Zhang, F. S.; Li, H. X.; Wang, Z. J. Isothermal Solubility of $\text{Na}_2\text{B}_4\text{O}_7$ - Na_2CO_3 - NaHCO_3 - NaBO_2 - H_2O System at 25, 35 and 45°C. *Acta Phys. Chim. Sin.* **2001**, *17*, 1045–1048.
- (13) Li, J.; Yang, Q.; Zhang, F. S.; Li, H. X.; Wang, Z. J. Study on the Isothermal Solubility of $\text{Na}_2\text{B}_4\text{O}_7$ - Na_2CO_3 - NaHCO_3 - NaBO_2 - H_2O System at 0 and 15 °C. *Chin. J. Inorg. Chem* **2001**, *17* (6), 849–852.
- (14) Fang, C. H.; Niu, Z. D.; Liu, Z. Q.; Cheng, J. Q. Studies on the Metastable Phase Diagram in the Quinary System Na^+ , K^+/Cl^- , CO_3^{2-} , SO_4^{2-} - H_2O at 25°C. *Acta Chim. Sin.* **1991**, *49*, 1062–1070.

Received for review July 27, 2009. Accepted January 12, 2010. Financial support for this work was provided by the National Nature Science Foundation of China (No. 40673050), the Research Fund for the Doctoral Program of Higher Education from the Ministry of Education of China (20070616008), and the Scholarship Leaders Training Fund from Sichuan Province (2008-140).

JE900640Z