Facile Activation of Carbon-Fluorine Bonds in Saturated Fluoroalkyl Ligands by Coordinated Water in Fluoroalkyl Aqua Complexes of Rhodium

Russell P. Hughes* and Danielle C. Lindner

Department of Chemistry, 6128 Burke Laboratory Dartmouth College, Hanover, New Hampshire 03755

Arnold L. Rheingold and Louise M. Liable-Sands

Department of Chemistry, University of Delaware Newark, Delaware 19716

Received July 14, 1997

Activation of normally inert carbon-fluorine bonds in saturated fluorocarbons is a topic of considerable recent interest.¹ The methodology of choice usually involves strongly reducing conditions; the fate of the fluorine is usually fluoride ion, stabilized as an alkali metal salt. An alternative method for activation of CF₂ groups is hydrolysis, with strong H-F and C=O bonds providing driving force for the reaction. However, conditions necessary to hydrolyze a CF₂ group in saturated fluorocarbons often require a combination of strong bases or acids, heat, and long reaction times, depending on the position of the hydrolyzable group in the molecule.² Hydrolysis of CF₂ or CF₃ groups bound to transition metal centers is sometimes more facile, but initial activation of an α -C-F bond still requires a strong Lewis acid, such as BF₃,³ SbF₅,⁴ or SiMe₃⁺,⁵ although there are some indications that a proton can act as the initial fluoride acceptor.⁶ We have reported remarkably facile hydrolysis of a Rh-CF2 group in a coordinatively unsaturated perfluorometallacyclobutene complex by traces of moisture present on glassware surfaces; a coordinatively saturated analogue was inert, leading to the hypothesis that the vacant coordination site on the metal could bind and activate water in the hydrolysis mechanism.⁷ Here, we describe the synthesis and structures of two cationic complexes containing adjacent fluoroalkyl and water ligands, both of which undergo hydrolysis of a CF₂ group, the facility of which depends strongly on the hydrogen-bonding ability of the counterion.

Addition of perfluorobenzyl complex 1^8 or perfluoropropyl analogue 2^9 to AgBF₄ in moist CH₂Cl₂ affords BF₄⁻ salts of the cationic aqua complexes 3 and 4,⁹ the structures of which were confirmed by X-ray crystallographic studies.¹⁰ ORTEP diagrams for the two ion pairs, with selected bond distances and angles, are presented in Figures 1 and 2. Each molecule exhibits close hydrogen-bonding interactions of coordinated

(2) See: Chemistry of Organic Fluorine Compounds II; Hudlicky, M., Pavlath, A. E., Eds.; ACS Monograph 187; American Chemical Society: Washington, DC, 1995; Chapter 4. Chemistry of Organic Fluorine Compounds, 2nd ed.; Hudlicky, M., Ed.; John Wiley & Sons: New York, NY, 1976; Chapter 5. Chemistry of Organic Fluorine Compounds; Hudlicky, M., Ed.; Macmillan; New York, 1962

(3) Richmond, T. G.; Crespi, A. M.; Shriver, D. F. Organometallics 1984, 3, 314.

(4) Reger, D. L.; Dukes, M. D. J. Organomet. Chem. 1978, 153, 67.

(4) Reger, D. L.; Dukes, M. D. J. Organomet. Chem. 1978, 155, 67.
(5) Koola, J. D.; Roddick, D. M. Organometallics 1981, 10, 591.
(6) Clark, B. R.; Hoskins, S. V.; Roper, W. R. J. Organomet. Chem.
1982, 234, C9. Appleton, T. G.; Berry, R. D.; Hall, J. R.; Neale, D. W. J. Organomet. Chem. 1989, 364, 249. Burrell, A. K.; Clark, G. R.; Rickard, C. E. F.; Roper, W. R. J. Organomet. Chem. 1994, 482, 261.

(7) Hughes, R. P.; Rose, P. R.; Rheingold, A. L. Organometallics **1993**, *12*, 3109.

(8) Hughes, R. P.; Lindner, D. C.; Rheingold, A. L.; Yap, G. P. A. Organometallics **1996**, *15*, 5678.

(9) Experimental details and spectroscopic data for this compound are provided in the Supporting Information.



Figure 1. Molecular structure of 3; fluorines are shown as spheres for clarity. A full ORTEP appears in the Supporting Information. Selected bond distances (Å) and angles (deg): Rh-O, 2.164(7); Rh-P, 2.338(2); Rh-C(17), 2.113(9); C(17)-F(7), 1.388(10); C(17)-F(6), 1.412(10); F(7)-C(17)-F(6), 101.8(6); C(17)-Rh-O, 89.5(3); O-Rh-P, 90.4(2); C(17)-Rh-P, 88.1(2).



Figure 2. Molecular structure of 4; fluorines are shown as spheres for clarity. A full ORTEP appears in the Supporting Information. Selected bond distances (Å) and angles (deg): Rh-O, 2.219(5); Rh-P, 2.319(2); Rh-C(11), 2.086(7); C(11)-F(1), 1.396(8); C(11)-F(2), 1.389(8); F(1)-C(11)-F(2), 102.1(6); C(11)-Rh-O, 85.2(3); O-Rh-P, 85.8(2); C(11)-Rh-P, 92.8(2).

water molecule with the tetrafluoroborate counterion, with O-F(9A) and O-F(11) distances of 2.640 and 2.855 Å in 3, and 2.723 and 2.636 Å in 4; such interactions are well precedented.¹¹ Close approaches of O–F(1) [2.749 Å] in ${\bf 3}$ and O-F(2) [2.767 Å] in 4 are also well within potential hydrogenbonding distance.

Solution spectra of 3 are consistent with the solid state structure. Asymmetric and symmetric stretches of coordinated

⁽¹⁾ For leading recent references, see: Kiplinger, J. L.; Richmond, T. G. J. Am. Chem. Soc. 1996, 118, 1805. Burdeniuc, J.; Chupka, W.; Crabtree, R. H. J. Am. Chem. Soc. 1996, 118, 2525. Burdeniuc, J.; Jedlicka, B.; Crabtree, R. H. Chem. Ber. (Recueil) 1997, 130, 145. Saunders, G. C. Angew. Chem., Int. Ed. Engl. 1996, 35, 2615.

⁽¹⁰⁾ Crystal data for **3**: C₂₀H₂₆BF₁₁OPRh, red block, triclinic, $P\overline{1}$, a =8.682(3), b = 11.033(2), and c = 14.376(3) Å, $\alpha = 74.50(1)$, $\beta = 76.11-(3)$, and $\gamma = 73.02(3)^\circ$, V = 1249.4(7) Å³, Z = 2, $D_x = 1.691$ g cm⁻³, T (3), and $\gamma = 75.02(3)$, $\gamma = 1249.4(7)$ A, Z = 2, $D_X = 1.091$ g cm⁻³, T = 2.98 K, R(F) = 6.74%, $R(wF^2) = 19.49\%$. Crystal data for 4: C₁₆H₂₆-BF₁₁OPRh, orange block, monoclinic, $P2_1/n$, a = 10.115(1), b = 16.270-(2), and c = 14.070(2) Å, $\beta = 92.74(1)^\circ$, V = 2312.7(7) Å³, Z = 4, $D_X = 1.689$ g cm⁻³, T = 298 K, R(F) = 4.81%, $R(wF^2) = 8.41\%$. Details of the crystallographic determinations are provided in the Supporting Information. (11) For a compilation of some crystallographically characterized organometallic complexes containing water ligands and a discussion of hydrogenbonding to counterions, see: Kubas, G. J.; Burns, C. J.; Khalsa, G. R. K.; Van Der Sluys, L. S.; Kess, G.; Hoff, C. D. *Organometallics* **1992**, *11*, 3390. For other more recent references, see: Veltheer, J. E.; Burger, P.; Bergman, R. G. *J. Am. Chem. Soc.* **1995**, *117*, 12478. Canty, A. J.; Jin, H.; Skelton, B. W.; White, A. H. J. Organomet. Chem. 1995, 487, 1. Dadci, L.; Elias, H.;
 Frey, U.; Hörnig, A.; Koelle, U.; Merbach, A. E.; Paulus, H.; Schneider, J. Frey, U., Holnig, A., Koene, U., Merbach, A. E., Patus, H., Schneider, J. S. Inorg. Chem. 1995, 34, 306. Carmona, D.; Cativiela, C.; García-Correas, R.; Lahoz, F. J.; Lamata, M. P.; López, J. A.; López-Ram de Víu, M. P.; Oro, L. A.; San José, E.; Viguri, F. J. Chem. Soc., Chem. Commun. 1996, 1247. Fries, A.; Green, M.; Mahon, M. F.; McGrath, T. D.; Nation, C. B. M.; Walker, A. P.; Woolhouse, C. M. J. Chem. Soc., Dalton Trans. 1996, 4517

Scheme 1



water are observed in the IR spectrum at 3630 and 3370 cm⁻¹.12 ¹⁹F NMR resonances for the diastereotopic fluorines of the CF₂ group are broad at room temperature, but sharpen when cooled to -50 °C to give a typical AX pattern. This temperature dependence may reflect dissociation of water¹³ with inversion at rhodium, thereby equilibrating the CF₂ environments; ¹⁰³Rh coupling to ³¹P is retained indicating that phosphine dissociation is not responsible for this phenomenon. In contrast, perfluoropropyl analogue 4^9 shows a sharp AX pattern for the α -CF₂ resonances at room temperature indicating that the molecule is not fluxional on the NMR time scale.

While perfluoropropyl complex 4 appears to be indefinitely stable on standing in solution at room temperature, perfluorobenzyl analogue 3 is cleanly transformed on standing overnight in CDCl₃ solution into pentafluorophenyl carbonyl complex 5;9 the CO stretch is observed in the IR spectrum at 2074 cm⁻¹, and the pentafluorophenyl ligand exhibits five resonances in the ¹⁹F NMR spectrum due to restricted rotation, as previously observed for other pentafluorophenyl rhodium complexes.14 HF can be observed in the volatiles, after vacuum transfer, as previously noted.8 A suggested mechanism is shown in Scheme 1; hydrolysis of the CF₂ group affords a coordinatively unsaturated pentafluorobenzoyl complex, which undergoes pentafluorophenyl migration into the vacant coordination site to afford 5. The initial step is suggested to involve a proton from coordinated water acting as a fluoride acceptor with loss of fluoride enhanced by resonance stabilization from Rh: presumably the difference in reactivity of the CF₂ groups toward hydrolysis in perfluorobenzyl complex 3 compared to perfluoropropyl analogue 4 rests in additional enhancement of the leaving group ability of the benzylic fluoride.

This transformation does not take place in acetone- d_6 , which displaces water from the coordination sphere of **3**,¹³ suggesting that it is a reaction of coordinated water with the adjacent CF₂ group which is important; addition of more water to an acetone solution failed to accomplish hydrolysis. In addition, reaction of 1 with silver triflate affords the triflate complex 6;⁹ the sharp AX pattern observed for the CF₂ group at room temperature indicates the triflate ligand is not labile. The triflate ligand is not displaced by water alone, and no hydrolysis of the CF2 group is observed. Bergman has shown that addition of the sodium salt of weakly coordinating anion tetrakis{3,5-bis(trifluoromethyl)phenyl}borate (B{Ar_F}₄) will displace triflate from iridium as sodium triflate;¹⁵ however, attempts to prepare the $B{Ar_F}_4$

salt of cation 3 from 6 using this method in the presence of water afforded very rapid hydrolysis to give the $B{Ar_F}_4^-$ salt 7^9 with no signs of any intermediate water complex analogous to 3.

The inability to isolate the perfluorobenzyl aqua cation using the B{Ar_F}₄⁻ counterion suggests that an anion incapable of either displacing water, or of hydrogen-bonding to it, might strongly enhance the ability of coordinated water to hydrolyze the adjacent CF₂ group. Loss of hydrogen bonding to the counterion will allow the proton of coordinated water a stronger interaction with the adjacent CF₂ group. In agreement with this hypothesis, the perfluoropropyl complex 2 reacts with silver triflate to afford 8^{9} ; treatment of 8 with NaB{Ar_F}₄ in the presence of water results in rapid hydrolysis to give the perfluoroethyl carbonyl complex 9.9 The BF₄⁻ analogue 10 was prepared unambiguously by the reaction of $[Rh(\eta^5-C_5Me_5) (C_2F_5)I(PMe_3)$] with AgBF₄ in the presence of CO.



In summary, we provide evidence that the hydrolysis reaction of CF₂ groups in perfluoroalkyl aqua complexes is a reaction of coordinated water, which can be extremely facile if the counterion is incapable of hydrogen bonding. This represents an unusually mild activation of saturated C-F bonds and suggests that other reactions of water in the coordination sphere of transition metals, particularly those in which it acts as a protic acid, might also benefit from the presence of a non-hydrogenbonding counterion.

Acknowledgment. R.P.H. acknowledges support of this research by the National Science Foundation and by the Donors of the Petroleum Research Fund, administered by the American Chemical Society.

```
JA9723388
```

⁽¹²⁾ Nakamoto, K. Infrared and Raman Spectra of Inorganic and Coordination Compounds, 4th ed.; Wiley: New York, NY, 1986.

⁽¹³⁾ The water ligands in [Rh(C₅Me₅)(OH₂)₃²⁺][OTf⁻]₂ are labile to exchange and substitution: Eisen, M. S.; Haskel, A.; Chen, H.; Olmstead, M. M.; Smith, D. P.; Maestre, M. F.; Fish, R. H. Organometallics 1995, 14, 2806.

⁽¹⁴⁾ Jones, W. D.; Partridge, M. G.; Perutz, R. N. J. Chem. Soc., Chem. Commun. 1991, 264. Belt, S. T.; Helliwell, M.; Jones, W. D.; Partridge, (15) Arndtsen, B. A.; Bergman, R. G. Science 1995, 270, 1970.

Supporting Information Available: Details of syntheses and characterizations of complexes 2-9, and the X-ray structure determinations, tables of atomic coordinates and anisotropic thermal parameters, bond lengths and bond angles, anisotropic displacement parameters, and hydrogen atom coordinates for 3 and 4 (24 pages). See any current masthead page for ordering and Internet access instructions.