


Communication

Solid-State O NMR as a Probe for Structural Studies of Proteins in Biomembranes

Vincent Lematre, Maurits R. R. de Planque, Andy P. Howes, Mark E. Smith, Ray Dupree, and Anthony Watts *J. Am. Chem. Soc.*, **2004**, 126 (47), 15320-15321• DOI: 10.1021/ja0473283 • Publication Date (Web): 06 November 2004 **Downloaded from http://pubs.acs.org on April 5, 2009**



More About This Article

Additional resources and features associated with this article are available within the HTML version:

- Supporting Information
- Links to the 3 articles that cite this article, as of the time of this article download
- Access to high resolution figures
- Links to articles and content related to this article
- Copyright permission to reproduce figures and/or text from this article

View the Full Text HTML





Published on Web 11/06/2004

Solid-State ¹⁷O NMR as a Probe for Structural Studies of Proteins in Biomembranes

Vincent Lemaître, †,‡ Maurits R. R. de Planque,‡ Andy P. Howes,§ Mark E. Smith,§ Ray Dupree,*,\$ and Anthony Watts*,\$

BioAnalytical Department Vers-Chez-Les-Blanc, Nestlé Research Center, CH-1000 Lausanne 26, Switzerland, Department of Biochemistry, University of Oxford, South Parks Road, Oxford OX1 3QU, U.K., and Department of Physics, University of Warwick, Coventry CV4 7AL, U.K.

Received May 7, 2004; E-mail: anthony.watts@bioch.ox.ac.uk (A.W.); R.Dupree@warwick.ac.uk (R.D.)

A key experimental challenge to understand conformational details of membrane proteins is to provide unambiguous atomic-scale information about the molecular bonding arrangement and any changes that occur upon receptor activation. This demands the development of experimental probe techniques to deliver this information of biological and pharmaceutical importance. Solid-state NMR is a nonperturbing approach which can be used to study ligandprotein interactions where molecular size is not limiting and crystallinity is not a requirement. As a first step in addressing this challenge by exploiting ¹⁷O NMR in biomembrane applications, we report here the ¹⁷O solid-state NMR spectra at high field of an ¹⁷O selectively labeled transmembrane peptide in a biomimetic

Oxygen plays a key role in the molecular conformation of biopolymers. Among the important information that can be obtained from ¹⁷O, the only NMR-active oxygen isotope, is the isotropic chemical shift (δ_{iso}), the quadrupolar coupling constant (C_0), and the asymmetry parameter (η) . These parameters are extremely sensitive to the electronic distribution around the nucleus; more specifically, they are sensitive to its protonation state² and its involvement in hydrogen bonds. Furthermore, they contain structural information,²⁻⁵ and several methods for determining internuclear distances between oxygen and other nuclei have been developed.^{2,6,7} This suggests that oxygen could play a central role in biological NMR studies. However, ¹⁷O has a low natural abundance (0.037%) and a spin $(I = \frac{5}{2})$ with a corresponding quadrupole interaction that is manifested as significantly broadened signals in NMR spectra. Consequently, ¹⁷O solid-state NMR studies are still relatively uncommon, and selective labeling is essential. Despite these difficulties, in recent years, with the advent of higher magnetic fields and techniques for improving resolution beyond magic angle spinning experiments (MAS), there has been a significant increase in ¹⁷O NMR reports from inorganic materials, such as glasses and zeolites.^{8,9} There has been much less ¹⁷O NMR reported from organic materials since ¹⁷O presents even more of a challenge due to the larger quadrupole interaction and, hence, larger line widths. 19 Recent reports of ¹⁷O NMR from biologically relevant materials have included inorganic molecules interacting with hemeproteins, 10 polypeptides,^{11,12} amino acids,² and nucleic acid bases.¹³

Here, we report the first example of ¹⁷O NMR spectra from a selectively labeled transmembrane peptide, ¹⁷O-[Ala12]-WALP23, as a lyophilized powder and incorporated in hydrated vesicles, opening up new possibilities for applications of ¹⁷O solid-state NMR on real biological systems. WALP23 is a synthetic peptide which represents a consensus sequence for transmembrane protein segments.¹⁴ This hydrophobic peptide forms well-defined and wellcharacterized transmembrane α-helices¹⁴ and has special relevance




Figure 1. ¹⁷O NMR spectra of lyophilized ¹⁷O-[Ala12]-WALP23 at room temperature: (A) at 14.1 T, 4 mm rotor spinning at ~15 kHz, with simulation; (B) at 18.8 T, 2.5 mm rotor spinning at ~22 kHz, with simulation. Spinning sidebands are marked with an asterisk (*).

for the packing of α-helices in polytopic proteins, such as G-proteincoupled receptors.

Fmoc-L-alanine was prepared starting from ¹⁷O uniformly enriched L-alanine following a previously described procedure¹⁵ using H₂¹⁷O (70% enriched at ¹⁷O). WALP23 was synthesized and transmembranously incorporated into multilamellar vesicles (MLVs) of 1,2-distearoyl-sn-glycerol-3-phosphatidylcholine (DSPC) lipids. 14 The ¹⁷O NMR experiments were carried out on Chemagnetics Infinity 600 and Varian/Magnex 800 wide-bore 89 mm spectrometers at frequencies of 81.345 and 108.419 MHz for ¹⁷O, respectively. MAS was performed at a spin rate of 11-22 kHz. Spectral simulations were carried out using DMFit software.¹⁶

The ¹⁷O MAS NMR spectra of lyophilized ¹⁷O-[Ala12]-WALP23, acquired at magnetic field strengths of 14.1 and 18.8 T, show a single resonance centered 17 at $\sim 280-300$ ppm from the single label at alanine-12 (Figure 1). Spinning sidebands of the central transition are also visible. The decrease of the second-order quadrupole broadening and the associated increase in spectral resolution with increasing applied magnetic field is clearly evident. The simulation parameters for ¹⁷O in this alanine environment are summarized in Table 1.

Since WALP23 is lyophilized from a solution of 2,2,2-trifluoroethanol, the α -helical conformation should be retained. The isotropic chemical shift for the label (317.5 ppm) is very close to that of another polymer, poly-L-alanine, which also forms ideal α -helices (319 ppm¹¹). The difference between these two samples

Nestlé Research Center.

University of Oxford.
University of Warwick.

Table 1. 17O NMR Interaction Parameters for L-Alanine in Two Protonation States, [Ala12]-WALP23 Lyophilized and in Dry and Hydrated MLVs (1:10 peptide/lipid molar ratio), and Poly-L-alanine

site	$C_{\rm Q}$ (MHz)	η	δ_{iso} (ppm)	assignment
L-Alanine ^a				
1	7.86 ± 0.05	0.28 ± 0.03	284.0 ± 0.5	O1
2	6.53 ± 0.05	0.70 ± 0.03	260.5 ± 0.5	O2
L-Alanine•HCla				
1	8.31 ± 0.02	0.00 ± 0.02	327.8 ± 0.5	C=O
2	7.29 ± 0.02	0.25 ± 0.02	176.7 ± 0.5	C-OH
Fmoc Alanine ^a				
1	7.89 ± 0.05	0.16 ± 0.02	-	C=O
2	6.95 ± 0.05	0.12 ± 0.02	175.7 ± 0.5	C-OH
Lyophilized [Ala12]-WALP23				
1	8.45 ± 0.10	0.21 ± 0.03		C=O
[Ala12]-WALP23 in DSPC Vesicles				
1	8.42 ± 0.10	0.21 ± 0.03	311 ± 1	C=O
[Ala12]-WALP23 in Hydrated DSPC Vesicles				
1	8.55 ± 0.15	0.24 ± 0.03	315 ± 1	C=O
Poly-L-alanine ^b				
1	8.59	0.28	319	C=O

^a From ref 2. ^b From ref 11.

is only \sim 1.5 ppm. This suggests only small differences in local environment for the two peptides, although poly-L-alanine is crystalline and 17 O-[Ala12]-WALP23 is not.

The ¹⁷O MAS NMR spectrum of ¹⁷O-[Ala12]-WALP23 in hydrated DSPC vesicles (Figure 2) also shows a single resonance from the label inserted in alanine-12 centered at ~280 ppm. There is a partial overlap between the spinning sidebands and the central transition because of the lower spinning speed, which was necessary for this sample to maintain vesicle integrity. The hydration water is also visible (0 ppm) but does not interfere with the main signal from the ¹⁷O label. The lower spinning speed does not prevent determination of the ¹⁷O NMR parameters, which are given in Table 1. This experiment, which was designed as a test to evaluate the sensitivity of the method, shows that it is possible to study ¹⁷O selectively labeled peptides in hydrated biomembranes.

We recently obtained a correlation between the isotropic chemical shift and the C=O bond lengths determined by diffraction.² In addition, *ab initio* calculations on glutamic acid polymorphs also showed a correlation between these parameters.¹⁸ These two correlations were very similar. On the basis of this data, the C=O bond length at the labeled site for lyophilized WALP is 1.217 Å, which is very close to the 1.220 Å bond length for an ideal α -helix. For the WALP peptide incorporated in hydrated vesicles, the corresponding value is 1.223 Å, an increase in C=O bond length of \sim 0.006 Å. To further test this correlation, it was applied on other polymers studied previously by ¹⁷O solid-state NMR (polyglycine II and poly-L-alanine I and II; data not shown) and was shown to give bond lengths consistent with those of the known structures.

The data presented here are the first applications of ¹⁷O solid-state NMR to characterize and estimate bond distances within a selectively labeled peptide of biological interest. The availability of high-field NMR spectrometers (here, up to 18.8 T, with 89 mm bore magnet) allows useable S/N ratios to be readily achieved. It should be emphasized that this included a sample with <5% of the amino acid residues labeled and which was further diluted by a factor of 10–20 because of the presence of lipids and water. The ability to detect and characterize a signal from a selectively labeled peptide suggests that ¹⁷O NMR will be a fruitful experimental approach for studying hydrogen bonding in macromolecules. It could find widespread application in the study of selectively labeled organic molecules interacting with crystallized proteins or with membrane proteins in biological membranes. The high sensitivity




Figure 2. ¹⁷O NMR spectra of ¹⁷O-[Ala12]-WALP23 in hydrated DSPC vesicles (1:10 peptide/lipid molar ratio) at 14.1 T, 4 mm rotor spinning at \sim 11 kHz, with simulation. The lyophilized MLV sample was hydrated with one weight equivalent of water, and the spectrum was acquired at room temperature, with the lipids in the liquid—crystalline phase.

of solid-state ¹⁷O NMR to the local bonding environment also allows very accurate determination of small changes in C=O bond lengths.

Acknowledgment. This work was supported by EPSRC, BBSRC, and EMBO. We thank Dr. Dirk Rijkers and Dr. Dahlia Fischer for their advice on the peptide synthesis. Magnex and Varian are acknowledged for their effort in building the first commercially available 800 MHz ¹H WB 89 mm NMR spectrometer.

Supporting Information Available: Details for the sample preparation and the NMR experiments. This material is available free of charge via the Internet at http://pubs.acs.org.

References

- (1) Watts, A. Curr. Opin. Biotechnol. 1999, 10, 48-53.
- (2) Pike, K. J.; Lemaître, V.; Kukol, A.; Anupõld, T.; Samoson, A.; Howes, A. P.; Watts, A.; Smith, M. E.; Dupree, R. J. Phys. Chem. B 2004, 108, 9256–9263.
- (3) Seliger, J. Chem. Phys. 1998, 231, 81-86.
- (4) Asakawa, N.; Kameda, T.; Kuroki, S.; Kurosu, H.; Ando, S.; Ando, I.; Shoji, A. Annu. Rep. NMR Spectrosc. 1998, 35, 55-137.
- (5) Kuroki, S.; Ando, S.; Ando, I. Chem. Phys. 1995, 195, 107-116.
- (6) Chopin, L.; Vega, S.; Gullion, T. J. Am. Chem. Soc. 1998, 120, 4406–4409.
- Hu, K. N.; Iuga, D.; Griffin, R. G. Dynamic Nuclear Polarization Enhanced ¹⁷O Solid State NMR Spectroscopy on Biological Solids. Presented at the 44th Experimental NMR Conference, Savannah, Georgia, Mar 30–Apr 4, 2003; p 198.
 MacKenzie, K. J. D.; Smith, M. E. Multinuclear Solid State NMR of
- (8) MacKenzie, K. J. D.; Smith, M. E. Multinuclear Solid State NMR of Inorganic Materials; Pergamon Press: Oxford, 2002.
 (9) Bull, L. M.; Bussemer, B.; Anupöld, T.; Reinhold, A.; Samoson, A.; Sauer,
- (9) Bull, L. M.; Bussemer, B.; Anupöld, T.; Reinhold, A.; Samoson, A.; Sauer J.; Cheetham, A. K.; Dupree, R. J. Am. Chem. Soc. 2000, 122, 4948–4958.
- (10) McMahon, M. T.; deBios, A. C.; Godbout, N.; Salzmann, R.; Laws, D. D.; Le, H.; Havlin, R. H.; Oldfield, E. J. Am. Chem. Soc. 1998, 120, 4784–4797
- (11) Yamauchi, K.; Kuroki, S.; Ando, I.; Ozaki, T.; Shoji, A. Chem. Phys. Lett. **1999**, 302, 331–336.
- (12) Yamauchi, K.; Kuroki, S.; Ando, I. J. Mol. Struct. 2002, 602–603, 171–175.
- (13) Wu, G.; Dong, S.; Ida, R.; Reen, N. J. Am. Chem. Soc. 2002, 124, 1768–1777.
- (14) de Planque, M. R.; Kruijtzer, J. A.; Liskamp, R. M.; Marsh, D.; Greathouse, D. V.; Koeppe, R. E., II.; de Kruijff, B.; Killian, J. A. J. Biol. Chem. 1999, 274, 20839–20846.
- (15) Steinschneider, A.; Burgar, M. I.; Buku, A.; Fiat, D. Int. J. Pept. Protein Res. 1981, 18, 324–333.
- (16) Massiot, D.; Fayon, F.; Capron, M.; King, I.; Le Calvi, S.; Alonso, B.; Durand, J.-O.; Bujoli, B.; Gan, Z.; Hoatson, G. Magn. Reson. Chem. 2002, 40, 70
- (17) For a quadrupolar nucleus, the second-order quadrupolar broadening means that the isotropic shift is near the left-hand edge of the line (see ref 8 above).
- (18) Yates, J. R.; Pickard, C. J.; Paynes, M. C.; Dupree, R.; Profeta, M.; Mauri, M. J. Phys. Chem. A **2004**, 108, 6032–6037.
- (19) Lemaître, V.; Smith, M. E.; Watts, A. Solid State Nucl. Magn. Reson. 2004, 26, 215–235.

JA0473283