This article was downloaded by:


On: 24 January 2011

Access details: Access Details: Free Access

Publisher Taylor & Francis

Informa Ltd Registered in England and Wales Registered Number: 1072954 Registered office: Mortimer House, 37-

41 Mortimer Street, London W1T 3JH, UK



# Journal of Liquid Chromatography & Related Technologies

Publication details, including instructions for authors and subscription information: http://www.informaworld.com/smpp/title~content=t713597273

# The Determination of Methylene Bis-Thiocyanate by High Performance Liquid Chromatography

R. Shustina<sup>a</sup>; J. H. Lesser<sup>a</sup>

<sup>a</sup> Department of Analytical Research Research and Development Division, Makhteshim Chemical Works, Ltd., Beer Sheva, Israel

**To cite this Article** Shustina, R. and Lesser, J. H.(1984) 'The Determination of Methylene Bis-Thiocyanate by High Performance Liquid Chromatography', Journal of Liquid Chromatography & Related Technologies, 7: 13, 2653 — 2658 **To link to this Article: DOI:** 10.1080/01483918408067032

URL: http://dx.doi.org/10.1080/01483918408067032

## PLEASE SCROLL DOWN FOR ARTICLE

Full terms and conditions of use: http://www.informaworld.com/terms-and-conditions-of-access.pdf

This article may be used for research, teaching and private study purposes. Any substantial or systematic reproduction, re-distribution, re-selling, loan or sub-licensing, systematic supply or distribution in any form to anyone is expressly forbidden.

The publisher does not give any warranty express or implied or make any representation that the contents will be complete or accurate or up to date. The accuracy of any instructions, formulae and drug doses should be independently verified with primary sources. The publisher shall not be liable for any loss, actions, claims, proceedings, demand or costs or damages whatsoever or howsoever caused arising directly or indirectly in connection with or arising out of the use of this material.

# THE DETERMINATION OF METHYLENE BIS-THIOCYANATE BY HIGH PERFORMANCE LIQUID CHROMATOGRAPHY

R. Shustina and J. H. Lesser Department of Analytical Research Research and Development Division, Makhteshim Chemical Works, Ltd. Beer Sheva, 84100, Israel

#### **ABSTRACT**

A fast, straightforward and specific method for the determination of methylene bis-thiocyanate (MBT) in technical materials and aqueous samples is described. Reversed phase high precision liquid chromatography using acetonitrile-water as eluent at a flow rate of 2 ml/min. gave a capacity factor of 2.73 for MBT. Using the method for aqueous samples allowed the determination of 2 ppm MBT by direct injection.

#### INTRODUCTION

Methylene bis-thiocyanate (MBT) is a fungicide and slimicide used extensively in the wood-pulp industry. It can be synthesized by mixing dibromomethane with sodium thiocyanate in an appropriate solvent (1). Methods of analysis for thiocyanates in general include those involving conversion to heavy metal derivatives (2,3), or to cvanide ion (4). Additional methods include reaction of alkyl and aryl thiocyanates with sulfide (5) and a micromethod using polysulfide (6). These methods are general, work intensive and require several hours for completion. An infrared (7,8) as well as an NMR (9) method represent instrumental analysis reported. however neither is applicable to low concentrations or to aqueous thiocyanate samples.

A rapid and specific analytical method for the determination of MBT applicable to technical, aqueous, and low concentration samples was required. This paper presents an HPLC method for the determination of MBT which meet the above requirements.

#### **EXPERIMENTAL**

#### Materials and Reagents

Acetonitrile Spectroscopic grade (Aldrich Chemical Co., U.S.A.), hexane, Analytical grade (Merck, Dramstadt, FGR) and dichloromethane (Analytical grade, Frutarom, Haifa, Israel) were used. Water was distilled and passed through a Norganic cartridge (Millipore Corporation, Massachussets, U.S.A.). Technical MBT samples were supplied by Makhteshim Chemical Works, Ltd. Aqueous solutions of MBT were prepared by dissolving a quantity of technical material in water and diluting to the desired concentration. Non-aqueous technical samples were prepared by dissolving MBT in acetonitrile prior to analysis.

### High Performance Liquid Chromatography

A Lichrosorb RP-18 octadecylsilane reverse phase column (Merck, Darmstadt, FGR), 12.5 cm 4.6 mm i.d. was used in a Varian model 5020 HPLC equipped with a Rheodyne model 7125 injector and UV-1 ultraviolet detector (254nm). MBT was eluted using acetonitrile: water (12:88) and a flow rate of 2 ml/min. Peak areas were integrated electronically using a Varian model CDS-111L electronic integrator and recorded using a Pantos Unicorder model U-228 (Pantos, Japan).

# Linearity Study

Into a 100 ml volumetric flask was weighed 0.222 g. of MBT and acetonitrile was added to the mark. From this stock solution a series of dilutions were made ranging from 0.148 mg/ml to 2.2 mg/ml. 10  $\mu$ l of these solutions were injected onto the column in duplicate. The resultant peak areas were plotted against the concentration and are represented in Figure 1.




Figure 1. Curve obtained for MBT when peak area versus concentration injected was plotted.

#### Relative Standard Deviation

A solution of 0.139 g of MBT in 100 ml of acetonitrile was prepared and 13 successive 10  $\mu l$  injections made. From the resultant peak areas the relative standard deviation was calculated.

#### RESULTS AND DISCUSSION

A rapid specific method for the analysis of MBT in aqueous samples and technical materials has been developed using HPLC and




Figure 2. Typical Chromatogram of MBT

ultraviolet detection. An octadecylsilane reversed phase column was used with an acetonitrile; water mixture as eluant. A typical chromatogram obtained using this method is shown in Figure 2. Both thiocyanate and bromide ion elute rapidly from the system and their elution time was used as the to value in the capacity factor (k') calculation. The k' was calculated and found to be 2.73. Unreacted dibromomethane did not interfere with the analysis. For technical samples the total analysis time was about 10 min. When peak areas were plotted against concentration the graph shown in Figure 1 was obtained showing a linear response curve from 0.148 mg/ml to 2.2 mg/ml sample concentration. Calcul-




Figure 3. Chromatogram of 2 ppm MBT Aqueous Sample

ation of the relative standard deviation yielded a value of  $\pm$  0.6% based on peak areas.

In order to confirm that the chromatographic peak represented MBT alone several samples were chromatographed on a silica column using hexane; dichloromethane (10:90) as eluent. Only a single sharp peak was observed in these instances.

Aqueous samples were filtered and injected without further workup. As little as 2 ppm (S/N>2) could be detected in this manner as shown in Figure 3. The fungicides captan, foltan and captafol did not elute under the above conditions and therefore are not potential interferents.

In conclusion a method for the analysis of MBT, applicable to aqueous and non-aqueous samples, has been presented. Method

reproducibility, linear range, and limit of detection have been established under given experimental conditions.

#### REFERENCES

- Matt, J., "Sequesting solvents in ethylene bis (Thiocyanate) production", U. S. Pat. 3,524,871 (Nalco Chemical Co.) 18 Aug. 1970.
- M. Schuessler, and R. Pohloudek Fabini, "Identification of organically bound thiocyanate III. Organic thiocyanates", Pharm. Zentralhalle 100, 452-8, 1961.
- Mieczyslaw Wronski, "Mercurimetric determination of thiocyanates and isothiocyanates", Talanta 13 (8), 1145-9, 1966.
- Kemp W. E. "The Analysis of Aliphatic and Aromatic Thiocyanates", Analyst (London), 64, 648-53, 1939.
- Maija Saarivirta and Artturi T. Virtanen, "A method for estimating benzyl isothiocyanate, benzyl thiocyanate, and benzyl nitrile in the crushed, moistened seeds of Lepidium sativum", Acta Chem. Scand. 17, Suppl. 1, S74-S78, 1963.
- R. Pohloudeic Fabini, D. Goeckeritz, and H. Brueckner, "Organic thiocyanate compounds XXXIII. Quantitative determination of organic thiocyanates", Mikrochim. Acta 3, 588-98, 1967.
- R. B. Hannan, "Analysis of insecticides for thiocyanates by infrared spectroscopy", Chem. Specialties Mfrs. Assoc. Proc. Dec., 151-3, 1955.
- N. G. McTaggart, E. Thornton, and A. D. Harford, "The determination of pyrethrins and other insecticidal compounds by infrared spectrometry", Pyrethrum Post 4, No. 4, 12-15, 1958.
- A. Mathias, "Analysis of alkyl thiocyanates and isothiocyanates by N. M. R.", Tetrahedron <u>21</u> (5), 1073-5, 1965.