



Figure S18. Autoradiogram of 20% denaturing PAGE, showing the cleavage kinetics of 5'-³²P-labelled target RNA (**14**) by RNase H1 in the native AON (**1**)/RNA (**14**), 15-DPPZ AON (**7**)/RNA (**14**), 15-3T-DPPZ (**8**)/RNA (**14**), 15-3T-Cholest AON (**10**)/RNA (**14**) and 15-2C-Cholest AON (**11**)/RNA (**14**) hybrid duplexes. PDE-Ladder: snake venom PDE ladder. Conditions of cleavage reaction: RNA (0.047 μM) and AONs (5 μM) in buffer, containing 20 mM Tris-HCl (pH 8.0), 20 mM KCl, 10 mM MgCl₂ and 0.1 mM DTT at 21 °C, 0.06 or 0.12 U of RNase H. Total reaction volume is 30 μl. (See Materials and Methods section for full experimental details.)