







In the NMR world however, it makes perfect sense...

The dual ¹3C/¹H CryoProbe™ that fits into your 400 MHz NMR spectrometer

Our new generation of 5mm DCH CryoProbe™ is available now for the first time at 400MHz!

This makes your routine 400MHz spectrometer as sensitive as an 800. Just think what that could mean to you...

- ✓ A ¹³C spectrum on every sample

 ✓ Why not DEPT-135 and ¹³C HSQC spectrum as well?

At last, CryoProbe™ sensitivity for the organic chemist

Only from Bruker



