1104 Notizen

Evidence for Uptake of Plamid DNA into Intact Plants (*Lemna perpusilla*) Proved by an *E. coli* Transformation Assay

M. Frey, S. Rall, A. Roth, and V. Hemleben

Institut für Biologie II der Universität Tübingen, Lehrstuhl für Genetik, Auf der Morgenstelle 28, D-7400 Tübingen

Z. Naturforsch. **35 c**, 1104–1106 (1980); received November 6, 1980

DNA Uptake, Lemna perpusilla, Plasmids

The water plant Lemna perpusilla has been incubated with the E. coli plasmids pMB9 and pBR 325, respectively. Uptake of plasmids has been shown by subsequent transformation of E. coli cells to tetracycline resistance after treatment with Lemna DNA from plasmid-incubated plants. In 7 out of 15 assays we found stable transformants. From the transformation rate an amount of 10^{-4} to 10^{-6} µg plasmid DNA per 10 µg of plant DNA can be calculated.

Introduction

E. coli plasmids can be used as vectors for the transfer of eucaryotic or procaryotic genes into higher plant cells or intact plants (for review see [1]). The aim of this investigation was to show uptake of plasmids into an intact plant. As recipient we used the water plant Lemna perpusilla cultivated under sterile conditions. The plasmids pMB9 or pBR 325 conferring the tetracycline or the tetracycline, ampicillin and chloramphenicol resistance markers were applied to the intact plants. Whereas previous studies on DNA uptake by plant cells often included the analysis of the DNA preparation by caesium chloride gradient centrifugation [2-4] we tried as a more sensitive method for detection of plasmid DNA the bioassay of E. coli transformation. The plasmids of the resulting transformants were analyzed after a miniscreen procedure to see whether a rearrangement of plasmid occurred during the passage through the plant.

Materials and Methods

Recipient plants

Lemna perpusilla, strain 6746, was cultivated under sterile conditions in the medium described by

Reprint requests to Prof. Dr. Vera Hemleben. 0341-0382/80/1100-1104 \$ 01.00/0

Pirson and Seidel [5] at 27 °C and continous light of 3000 lux.

Plasmids and Bacteria

E. coli plasmids pMB9 and pBR 325 were harvested from E. coli K-12, HB101, kindly given by Prof. A. Pühler, University of Bielefeld, and by Dr. R. B. Meagher, University of Athens, Georgia. E. coli K-12, RRI [6] was used for the transformation assays.

Isolation of plant DNA and plasmids

DNA of Lemna perpusilla was isolated after homogenization in a glass homogenizer following a slightly modified procedure described by Hemleben et al. [3]. The purified and reprecipitated DNA was used in the transformation assay.

Plasmid DNA was isolated from *E. coli* by ethidium bromide-caesium chloride gradient centrifugation [7] after amplification with chloramphenicol ($200 \,\mu\text{g/ml}$ for pMB9) or with spectinomycin ($300 \,\mu\text{g/ml}$ for pBR 325).

Incubation of Lemna with plasmid DNA

Ca. 200 exponentially growing *Lemna* plants were transferred to 10–20 ml 0.01 M KCl-0.001 M sodium citrate buffer pH 6.7 containing sterile prepared plasmid DNA (50–70 μg/ml) and streptomycin (10 μg/ml). The incubation was carried out at 27 °C or 6 °C. The lower temperature improves the stability of the plasmids in the medium. The plants were in contact with the medium with the root-like organs and the lower surface of the leaves. At the end of the incubation time (5 min as control or 20–22 h) the plants were washed three times with distilled water, treated with DNase I in 0.002 M MgCl₂ (20 μg/ml, 15 min, 37 °C) and then with proteinase K (100 μg/ml, 30 min, 37 °C) before isolation of the DNA.

Transformation of E. coli

Competent *E. coli* cells were obtained by the $CaCl_2$ method [8]. $5-30\,\mu g$ *Lemna* DNA from plasmid-incubated plants or pure plasmid DNA $(0.0003-3\,\mu g)$ were added to $200\,\mu l$ of competent *E. coli* K-12, RRI cells.



Dieses Werk wurde im Jahr 2013 vom Verlag Zeitschrift für Naturforschung in Zusammenarbeit mit der Max-Planck-Gesellschaft zur Förderung der Wissenschaften e.V. digitalisiert und unter folgender Lizenz veröffentlicht: Creative Commons Namensnennung-Keine Bearbeitung 3.0 Deutschland Lizenz.

This work has been digitalized and published in 2013 by Verlag Zeitschrift für Naturforschung in cooperation with the Max Planck Society for the Advancement of Science under a Creative Commons Attribution-NoDerivs 3.0 Germany License.

Notizen 1105

Table I. Incubation of Lemna perpusilla with the E. coli plasmids pMB9 (exp. 1-3) and pBR 325 (exp. 4-7) and subsequent transformation of competent E. coli cells to tetracycline resistance with the DNA isolated from the plasmid-incubated plants or pure plasmid as control.

Plasmid concen- tration [µg/ml]	Incubation conditions	Lemna DNA in the transforma- tion assay [µg]	Resulting transfor- mants [ml]	Transformation rate/µg Lemna DNA	Transformation rate/µg pure plasmid	Calculated µg of plasmid DNA/10 µg <i>Lemna</i> DNA
1 a) 57	5 min/27 °C	10	0	0	4.5 × 10 ⁻⁴	0
1 b) 57	20 h/27 °C	10	11.4	3.0×10^{-9}	4.5×10^{-4}	8.9×10^{-5}
2) 65	20 h/27 °C	5	5.0	1.6×10^{-8}	1.7×10^{-3}	2.3×10^{-4}
3 a) 75	5 min/27 °C	30	0	0	8.6×10^{-5}	0
3 b) 75	20 h/27 °C	30	8.6	2.7×10^{-10}	8.6×10^{-5}	3.1×10^{-5}
4) 70	22 h/27 °C	18	5.0	6.5×10^{-11}	6.7×10^{-6}	1.7×10^{-4}
5 a) 57	22 h/27 °C	20	6.7	3.8×10^{-9}	1.9×10^{-3}	2.2×10^{-5}
5 b) 57	22 h/6 °C	20	33.3	1.9×10^{-8}	1.9×10^{-3}	9.6×10^{-5}
6 a) 57	22 h/27 °C	18	5.0	6.6×10^{-10}	5.5×10^{-3}	1.6×10^{-5}
6 b) 57	22 h/6 °C	24	13.3	1.3×10^{-9}	5.5×10^{-3}	3.0×10^{-5}
7) 73	22 h/6 °C	20	1.8	2.1×10^{-10}	5.5×10^{-3}	4.5×10^{-6}

Miniscreen

Plasmids present in the transformed bacteria were isolated as described by Birnboim and Doyle [9], treated with the restriction endonuclease *HindIII* and analyzed by gel electrophoresis on 1% agarose gels.

Results and Discussion

Lemna plants were incubated with pMB9 or pBR 325 plasmid DNA under various conditions for 0 to 22 h at 6 or 27 °C, respectively. After the incubation the plants were thoroughly washed, treated with DNase I twice to remove the plasmids from the plant surface and subsequently with proteinase K and then the plant DNA was isolated. These Lemna DNA preparations were used for transformation of E. coli cells to tetracycline resistance. Table I shows the results. Out of 15 incubation assays 7 Lemna DNA preparations from plasmid-incubated plants transformed E. coli with a transformation rate of $10^{-8}-10^{-10}/\mu g$ DNA. These transformants were not due to contamination of the plants or of the competent E. coli cells. Both possibilities were tested by plating incubation medium on nutrient broth agar plates and competent cells on tetracycline-containing plates. DNA from Lemna plants treated only for 5 min with plasmid DNA never resulted in transformed E. coli cells.

The amount of plasmid DNA present in *Lemna* DNA preparations could be determined by compar-

ing the transformation rates with those obtained after treatment of competent E. coli cells with pure plasmid DNA (Table I). These rates are linear to the amount of plasmid between 3 to 0.0003 µg plasmid DNA applied to E. coli cells. Mixing plant DNA and plasmids in vitro (100:1) this only slightly lowered the transformation rate. Therefore it is possible to calculate the portion of plasmid DNA taken up by the plants. Per 10 µg of Lemna DNA 10⁻⁴ to 10⁻⁶ µg plasmid DNA are estimated. Regaining ca. 0.5 µg DNA per Lemna plant this means that about 100-900 intact plasmid molecules are taken up by one Lemna plant. Using intact plants as recipients one could argue that the plasmids do not actually enter the plant cells or even the cell nucleus but stay in the intracellular spaces. One argument against this possibility is that in several experiments we did not get any transformants.

The transformants obtained after treatment of $E.\ coli$ with DNA from plasmid incubated plants were analyzed in more detail after a miniscreen with the restriction endonuclease HindIII which cuts the normal pMB9 and pBR 325 plasmids only once giving rise to one band on agarose gels. Some of these Lemna-plasmid transformants contained in addition a new band with a molecular weight of ca. 3×10^5 , whereas most of these appeared unchanged. This phenomenon is under further investigation and might be an indication that some kind of rearrangement of the plasmids occur in the plant cells.

1106 Notizen

- C. J. Kado and A. Kleinhofs, Int. Rev. Cyt., Suppl. 11 B, pp. 47-80 (1980).
 P. F. Lurquin and C. J. Kado, Mol. Gen. Genet. 154, 113-121 (1977).
 V. Hendelson, M. Ermisch, D. Kimmick, P. Leber, and
- [3] V. Hemleben, N. Ermisch, D. Kimmich, B. Leber, and G. Peter, Eur. J. Biochem. 56, 403-411 (1975).
 [4] W. Gradmann and V. Hemleben, Z. Naturforsch. 31 c,
- 558-564 (1976).
- [5] A. Pirson and F. Seidel, Planta (Berl.) 38, 431-473 (1950).
- [6] H. W. Boyer and D. J. Roulland Dussoix, J. Mol. Biol. **41**, 459 – 472 (1969). [7] M. Bazarel and R. Helinski, J. Mol. Biol. **36**, 185–194
- (1968).
 [8] S. N. Cohen, A. C. Y. Chang, and L. Hsu, Proc. Nat. Acad. Sci. USA 69, 2110-2114 (1972).
 [9] H. C. Birnboim and J. Doyle, Nucl. Acids Res. 7,
- 1513-1523 (1979).