

**Agilent Technologies  
E444xA Option H26**

**User's and Service Guide**



# **Agilent Technologies**

## **E444xA Option H26**

### **User's and Service Guide**

**This guide applies to firmware revision A.06.01 or greater**

**Use this manual with the following documents:**

**Agilent Technologies PSA Series Spectrum Analyzer User's  
and Programmer's Reference Guide**



**Agilent Technologies**

**Manufacturing Part Number: E4440-90560**

**Printed in USA**

**March 2005**

© Copyright 2003-2005 Agilent Technologies

---

## **Warranty Statement**

THE MATERIAL CONTAINED IN THIS DOCUMENT IS PROVIDED "AS IS," AND IS SUBJECT TO BEING CHANGED, WITHOUT NOTICE, IN FUTURE EDITIONS. FURTHER, TO THE MAXIMUM EXTENT PERMITTED BY APPLICABLE LAW, AGILENT DISCLAIMS ALL WARRANTIES, EITHER EXPRESS OR IMPLIED WITH REGARD TO THIS MANUAL AND ANY INFORMATION CONTAINED HEREIN, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. AGILENT SHALL NOT BE LIABLE FOR ERRORS OR FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH THE FURNISHING, USE, OR PERFORMANCE OF THIS DOCUMENT OR ANY INFORMATION CONTAINED HEREIN. SHOULD AGILENT AND THE USER HAVE A SEPARATE WRITTEN AGREEMENT WITH WARRANTY TERMS COVERING THE MATERIAL IN THIS DOCUMENT THAT CONFLICT WITH THESE TERMS, THE WARRANTY TERMS IN THE SEPARATE AGREEMENT WILL CONTROL.

---

## **Safety Notes**

The following safety notes are used throughout this manual. Familiarize yourself with each of the notes and its meaning before operating this instrument. All pertinent safety notes for using this product can be found in the standard manual.

---

<b>WARNING</b>	<b>Warning denotes a hazard. It calls attention to a procedure which, if not correctly performed or adhered to, could result in injury or loss of life. Do not proceed beyond a warning note until the indicated conditions are fully understood and met.</b>
<b>CAUTION</b>	Caution denotes a hazard. It calls attention to a procedure that, if not correctly performed or adhered to, could result in damage to or destruction of the instrument. Do not proceed beyond a caution sign until the indicated conditions are fully understood and met.

---

---

## **Certification**

Agilent Technologies, Inc. certifies that this product met its published specifications at the time of shipment from the factory. Agilent Technologies, Inc. further certifies that its calibration measurements are traceable to the United States National Institute of Standards and Technology, to the extent allowed by the Institute's calibration facility, and to the calibration facilities of other International Standards Organization members.

---

## Definitions

- *Specifications* describe the performance of parameters covered by the product warranty (temperature - 0 to 55 °C, unless otherwise noted.)
- *Typical* describes additional product performance information that is not covered by the product warranty. It is performance beyond specification that 80% of the units exhibit with a 95% confidence level over the temperature range 20 to 30 °C. Typical performance does not include measurement uncertainty.
- *Nominal* values indicate expected performance, or describe product performance that is useful in the application of the product, but is not covered by the product warranty.

---

## Where to Find the Latest Information

Documentation is updated periodically. For the latest information about Agilent Technologies PSA Analyzers, including firmware upgrades and application information, please visit the following Internet URL:

<http://www.agilent.com/find/psa>

---

## DFARS/Restricted Rights Notice

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.



---

# Contents

## 1. Overview

Description .....	1-2
Characteristics .....	1-4
Maximum Input Level .....	1-5
VSWR .....	1-6
Displayed Average Noise Level (Preamp Off) .....	1-7
Noise Figure Measurement Personality and Option H26 .....	1-8
Performance Verification .....	1-12
Block Diagram .....	1-13
Replaceable Parts .....	1-14
Contacting Agilent .....	1-16

## 2. Example Measurements

What is in This Chapter .....	2-2
Equipment for Test Set-ups .....	2-3
Analog Test Set-up .....	2-4
Digital Test Set-up .....	2-6

## 3. User Interface

What is in This Chapter .....	3-2
Front Panel Key Select .....	3-3
AMPLITUDE Y Scale .....	3-3
GPIB Control of Preamp Gain .....	3-4
Preamp Gain Command Interface .....	3-4
Flatness Compensation Table .....	3-5
Accessing the Flatness Compensation Table .....	3-5
GPIB Control of the Flatness Compensation Table .....	3-5
Reloading the Flatness Correction Table .....	3-6
Recreating the Flatness Correction Table .....	3-7

## 4. Appendix

DANL Manual Performance Verification .....	4-2
Procedure .....	4-2

---

## **Contents**

# 1 Overview

## Description

The Agilent Technologies E444xA Option H26 consists of a high band preamp with a nominal gain of 27 dB at the front end of the spectrum analyzer. The performance and ease of use come from the high gain of the preamplifier, the ability to switch the preamp in and out and the use of the instruments flatness correction tables to reduce the amplitude errors.

The preamp has a frequency range from 3 GHz to the upper frequency limit of the analyzer. Refer to [Table 1-1](#). The preamp can be enabled or disabled via remote SCPI commands or from the front panel interface. The lowband preamp (Option 1DS) operates below 3 GHz and the highband preamp (Option H26) operates above 3 GHz.

**Table 1-1**      **Preamp frequency range**

Model	Preamp Frequency range
E4440A	3 GHz to 26.5 GHz
E4443A	3 GHz to 6.7 GHz
E4445A	3 GHz to 13.2 GHz
E4446A	3 GHz to 44 GHz
E4448A	3 GHz to 50 GHz

---

**NOTE** All instrument specifications are considered to be **NOMINAL** values of performance when the highband preamp is engaged. For example, specifications for absolute amplitude accuracy are now **NOMINAL** in performance and may still be useful but are not specified.

---

**CAUTION** **Excessive Relay Switch wear can occur if the instrument is sweeping from below 2.85 GHz and the Stop Frequency is above 3.05 GHz.**

---

**CAUTION** **The Damage Level is reduced to +20 dBm when the Option H26 is enabled.**

Flatness characteristics of the preamp will be stored in the User memory file titled “Other.” This can be used to compensate for the variations in amplitude of the High Band PreAmp. Refer to [“Accessing the Flatness Compensation Table” on page 3-5](#).

When not using the highband preamp, displayed average noise level may still be degraded due to the extra length of cable in the option. Refer to [“Displayed Average Noise Level \(Preamp Off\)” on page 1-7](#).

If the instrument preamp is “ON” and sweeping from below 3 GHz to above 3 GHz, the instrument will automatically switch the highband preamp to the “ON” state above 3 GHz. Refer to [Figure 1-1](#).

**Figure 1-1      Lowband vs. Highband Frequency Diagram for E4440A<sup>1</sup>**



\* The “3.05 GHz” switch point occurs when sweeping from below 3 GHz to above 3 GHz. The highband preamp low frequency point is extended down to 3.0 GHz if the start frequency is set to 3.0 GHz.

---

**NOTE** Please note that in Full Band Sweep, Retrace time has been increased to reduce relay wear at band-cross. To reduce the retrace time, sweep completely in either side of the band-cross frequency.

---

---

**CAUTION** **Excessive Relay Switch wear can occur if the instrument is sweeping from below 2.85 GHz and the Stop Frequency is above 3.05 GHz.**

---

To find descriptions of specific analyzer functions refer to the *Agilent Technologies PSA Spectrum Analyzers User’s and Programming Reference Guide*.

---

1. See Table 1-1, “Preamp frequency range” for other model numbers.

## Characteristics

**Table 1-2 Option H26 Nominal Performance when the Preamp is On**

Description	Nominal Performance		Supplemental Information
	<b>E4440A, E4443A, E4445A</b>	<b>E4446A, E4448A</b>	
<b>Preamp Gain</b>			
3 GHz to 30 GHz	27 dB	27 dB	
Above 30 GHz	n/a	24 dB	
<b>Amplitude</b>			
Measurement Range	DANL to +16 dBm	DANL to +20 dBm	
<b>Max Safe Input Level</b>			<a href="#">page 1-5</a>
Avg Total Power	+20 dBm	+20 dBm	See Caution below
<b>Gain compression</b>			
1 dB Gain Comp Point <sup>a</sup>	-6 dBm	-6 dBm	<a href="#">Table 1-3 on page 1-5</a>
<b>VSWR</b>			<a href="#">page 1-6</a>
<b>DANL degradation<sup>b</sup></b>			<a href="#">page 1-7</a>
3 GHz	0.4 dB	0.4 dB	
12 GHz	0.7 dB	0.8 dB	
26 GHz	1.1 dB	1.4 dB	
50 GHz	n/a	2.2 dB	
<b>Noise Figure</b>			<a href="#">page 1-8</a>

a. at the Preamp RF input

b. with Option H26 preamp off.

---

**CAUTION**

**The Damage Level is reduced to +20 dBm when the Option H26 is enabled.**

---

## Maximum Input Level

The maximum allowable input level without damage is +20 dBm when the preamp is on. The input level calculation varies for different PSA models. See below.

### **E4440A, E4443A and E4445A:**

As mentioned, the maximum allowable input level without damage is +20 dBm when the preamp is on. This is not the power at the input connector minus the on-screen input attenuator setting. The maximum allowable input level is the power at the input connector minus the actual input attenuation. This is because the preamp is located after the 0-4 dB attenuator stages but before the 6, 10, 20, 30 dB attenuator stages. Therefore, the attenuation value entered on the instrument is in most cases, not the input attenuation to the preamp. Refer to [Table 1-3](#) for the actual input attenuation.

### **E4446A and E4448A:**

Here too, the maximum allowable input level without damage is +20 dBm when the preamp is on. For the E4446A and E4448A, the on-screen attenuation is the actual input attenuation to the preamp, therefore the maximum allowable input level without damage is the power at the input connector minus the on-screen attenuation. Refer to [Table 1-3](#) for the actual input attenuation.


**Table 1-3      Input Attenuator vs. Actual Input Attenuation**

<b>Input Attenuator Setting &gt; 3 GHz</b>	<b>Actual Input Attenuation at Preamp Input for E4440A, 43A, and 45A</b>	<b>Actual Input Attenuation at Preamp Input for E4446A, and E4448A</b>
0 dB	0 dB	0 dB
2 dB	2 dB	2 dB
4 dB	4 dB	4 dB
6 dB	0 dB	6 dB
8 dB	2 dB	8 dB
10 dB	0 dB	10 dB
12 dB	2 dB	12 dB


## VSWR

The VSWR graphs show the nominal instrument VSWR from 3 GHz to 26.5 GHz for the E4440A (see [Figure 1-2](#)), and from 3 GHz to 50 GHz for the E4448A (see [Figure 1-3](#)) with the Option H26 Preamp On, and 0 dB attenuation. The graphs show the behavior of two instruments, an E4440A and an E4448A. The E4443A and E4445A will have similar behavior to the E4440A and likewise the E4446A will have similar behavior to the E4448A.

**Figure 1-2**      **Nominal E4440A Instrument VSWR 3 GHz to 26.5 GHz;  
High Band Preamp On, Attenuation = 0 dB**



**Figure 1-3**      **Nominal E4448A Instrument VSWR 3 GHz to 50 GHz;  
High Band Preamp On, Attenuation = 0 dB**



## Displayed Average Noise Level (Preamp Off)

Due to the additional length of RF cable, and insertion loss of the RF switch by the inclusion of option H26, the standard instrument displayed average noise level (DANL) (H26 Off) may be nominally degraded. [Equation 1-1, “DANL degradation for E4440A, E4443A, and E4445A H26”](#) and [Equation 1-2, “DANL degradation for E4446A and E4448A H26”](#) are provided in calculating the degradation for specific frequencies other than those given in [Table 1-2 on page 1-4](#).

[“Appendix” on page 4-1](#) is provided to manually verify DANL performance.

### Equation 1-1 DANL degradation for E4440A, E4443A, and E4445A H26

$$DH26 = D + 0.2 + ((0.45/26.5) \times F) + 7.5 \times \left( \frac{5.758 + 4.89 \times F - 0.146 \times F^2 + 2.433 \times 10^{-3} \times F^3}{1200} \right)$$

### Equation 1-2 DANL degradation for E4446A and E4448A H26

$$DH26 = D + 0.15 + (0.02 \times F) + 7 \times \left( \frac{8.473 + 6.86 \times F - 0.159 \times F^2 + 1.831 \times 10^{-3} \times F^3}{1200} \right)$$

For all equations,  $F$  is the Frequency in GHz.

$D$  is the term used for representing the instrument’s DANL specification found in the PSA Series Specification Guide.

$DH26$  is the term used to represent the resultant instrument DANL from the inclusion of Option H26.

## Noise Figure Measurement Personality and Option H26

This section provides nominal performance information for the PSA series, Option 219, Noise Figure Measurement as it relates to the Option H26, High Band Preamp.

---

<b>CAUTION</b>	<b>Excessive Relay Switch wear can occur if the instrument is sweeping from a start frequency of 3.05 GHz and below to above 3.05 GHz. This is different than when operating in the SA Mode.</b>
----------------	--

---

You need the following equipment to use this personality/hardware combination.

**Table 1-4      Hardware, Firmware and Software Requirements**

Firmware	Software	Hardware	
Revision Number	Noise Figure Measurement Personality	Front End Driver Board	Option 1DS Internal Pre-Amp
≥ A.06.01	Option 219	Rev b or later	Required

---

<b>NOTE</b>	The Noise Figure Measurement personality (Option 219) requires Revision "b" or later of the Front End Driver assembly. This supplies the +28 V output (labelled "NOISE SOURCE DRIVE OUT +28 V (PULSED)" on the rear panel), which is needed to drive the noise source. To see which version is installed on your PSA, press <b>System</b> , <b>Show Hdwr</b> . If you have an earlier revision than Revision "b," contact your Agilent Technologies representative. Refer to <a href="http://www.agilent.com/find/psa">http://www.agilent.com/find/psa</a> for further information.
-------------	--

---


**Table 1-5      Instrument Uncertainty for Noise Figure with Option H26**

Description	Specifications	Supplemental Information
Instrument Uncertainty		Nominally the same as for the 10 MHz to 3 GHz range
3 to 10 GHz	There are no warranted specifications for Option H26. Refer to <a href="#">Table 1-2 on page 1-4</a> for Nominal Performance	H26 Preamp caution <sup>a</sup> Band Crossing caution <sup>b</sup>
10 to 20 GHz		Well-controlled preselector <sup>c</sup>
20 to 26.5 GHz		Good preselector stability <sup>d</sup> Preselector Drift Effects <sup>e</sup>


- a. The Special Option H26 preamp can reduce the total NF measurement uncertainty substantially above 3 GHz because it will reduce the effective noise figure of the measurement system, and thus it will reduce the sensitivity of the total NF uncertainty to the Instrument Gain Uncertainty. But if the signal levels into the preamp are large enough, the preamp may experience some compression. The compression differences between the noise-source-on and noise-source-off states causes an error that must be added to Instrument Noise Figure Uncertainty for use in the Noise Figure Uncertainty Calculator. Such signal levels are quite likely for the case where the DUT has some combination of high gain, high noise figure and wide bandwidth. Here's an example: The measurement will be made at 18 GHz. The typical preamp gain is 25 dB and the noise figure is 7 dB. We will assume the DUT has 20 dB gain, a 10 dB NF, and a passband from 5 to 30 GHz. We will use a noise source with 17 dB ENR. When the noise source is on, the DUT output can be computed by starting with  $kTB$  (-174 dBm/Hz) and adding  $10 \times \log(30 \text{ GHz} - 5 \text{ GHz})$  or 104 dB, giving -70 dBm for the thermal noise. Add to this the ENR of the noise source (17 dB) combined with the NF of the DUT (10 dB) to give an equivalent input ENR of 18 dB, thus -52 dBm input noise power. Add the gain of the DUT (20 dB) to find the DUT output power to be -32 dBm. The noise figure of the H26 preamp may be neglected. The H26 preamplifier gain of 25 dB adds, giving a preamplifier output power of -7 dBm. The typical 1 dB compression point of this amplifier at its output is +19 dBm. Therefore, the output noise is 26 dB below the 1 dB compression point. This amplifier will have negligible compression. As a rule of thumb, the compression of a noise signal is under 0.1 dB if the average noise power is kept 7 dB below the 1 dB CW compression point. The compression in decibels will usually double for every 3 dB increase in noise power. Use cases with higher gain DUTs could be compressed, leading to additional errors.
- b. The band 0 to band 1 crossing should be avoided. In addition to the wear-out mechanisms (see [Caution on page 1-8](#)) involved in measurements that overlap the 3 GHz band crossing, there will also be performance degradations. There will be thermal instabilities in such measurements that will add nominally 0.2 dB Instrument Uncertainty. The uncertainty of some NF or Gain measurements are greatly multiplied from the Instrument Uncertainty. See the Uncertainty Calculator included with the Noise Figure Measurement for details.
- c. In this frequency range, the preselector is well-controlled and there should be no need for special measurement techniques.
- d. In this frequency range, the preselector usually requires no special measurement techniques in a lab environment. But if the temperature changes by a few degrees, or the analyzer frequency is swept or changed across many gigahertz, there is a small risk that the preselector will not be centered well enough for good measurements.
- e. In this frequency range, the preselector behavior is not warranted. There is a modest risk that the preselector will not be centered well enough for good measurements. This risk may be reduced but not eliminated by using the analyzer at room temperature, limiting the span swept to a few gigahertz, and not changing the operating frequency range for many minutes.

Nominal values of Noise Figure are given in [Figure 1-4](#) and [Figure 1-5](#).

**Figure 1-4      Nominal E4440A Instrument Noise Figure;  
3 GHz to 26.5 GHz with H26 Preamp On**



**Figure 1-5      Nominal E4448A Instrument Noise Figure;  
3 GHz to 50 GHz with H26 Preamp On**




**Preamp ON:**

**Computed Measurement NF Uncertainty vs. DUT Gain, > 3 GHz  
(Non-warranted Frequency Range)**

Assumptions: Measurement Frequency 12 GHz, Instrument NF = 6.4 dB<sup>1</sup>, Instrument VSWR = 1.5<sup>2</sup>, Instrument Gain Uncertainty = 2.2 dB, Instrument NF Uncertainty = 0.05 dB, Agilent 346B Noise Source with Uncertainty = 0.2 dB, Source VSWR = 1.25, DUT input/output VSWR = 1.5.

**Figure 1-6 NF Uncertainty vs. DUT Gains > 3 GHz, Preamp On**



1. The preamp/analyizer combination NF is 6.4 dB; the internal preamp alone has a gain of 27 dB and a NF of 5 dB.
2. Instrument VSWR is now that of the internal preamp; VSWR = 1.5

## Performance Verification

To verify the performance of the E444xA Option H26, refer to [Figure 1-7](#). In this example, an 83650B Synthesized Sweeper was used, but a similar RF source could also be used.

---

**CAUTION**      **The Damage Level is reduced to +20 dBm when the Option H26 is enabled.**

---

1. Set the Signal Source to 5 GHz at -30 dBm. Verify the RF Power using a power meter. Connect the source to the PSA RF Input.
2. Measure the amplitude on E444xA H26 with preamp off.
3. Turn on the preamp, insuring that the Amp Corr is also on. Refer “[Front Panel Key Select](#)” on page 3-3 and “[Flatness Compensation Table](#)” on page 3-5 for more information.
4. Measure the amplitude on the PSA screen.

The measured amplitude should be  $\pm 5\text{dB}$  of the measured input signal from [step 2](#).


5. Repeat at 15 GHz, 20 GHz, 35 GHz, and 40 GHz.

---

**NOTE**      Even with option H26 disabled, the instrument may fail the DANL test when performing the *PSA Series Performance Verification Tests*. This is due to the insertion loss in the switches and cables used with H26. If a failure occurs, refer to “[Appendix](#)” on page 4-1 to manually check DANL. This check can be used to determine if the instrument needs to be troubleshooted for a higher than normal H26 DANL level.

---


**Figure 1-7      Performance Verification Set-up**



## Block Diagram


Refer to [Figure 1-8](#). The preamp and high frequency transfer switch are located between the first attenuator and the second attenuator.

**Figure 1-8 E4440A, E4443A, and E4445A H26 Block Diagram**



Refer to [Figure 1-9](#). The preamp and high frequency transfer switch are located between the attenuator and the highband mixer/preselector.

**Figure 1-9 E4446A and E4448A H26 Block Diagram**



## Replaceable Parts

**Table 1-6      E4440A, E4443A, and E4445A Option H26 Replaceable Parts**

Description	Agilent Part Number	Quantity
Microwave amplifier	0955-1663	1
Bracket, Option H26	E4440-00051	1
Semi-rigid cable assembly (Switch In)	E4440-20302	1
Semi-rigid cable assembly (Amp In)	E4440-20303	1
Semi-rigid cable assembly (Amp Out)	E4440-20304	1
Semi-rigid cable assembly (Switch Out)	E4440-20305	1
Cable assembly, switch control	E4440-60427	1
Board Assembly, PSA H26	E4440-60358	1
Cable assembly, amp bias	E4446-60067	1
Transfer switch 26 GHz	N1811TL-CFG002 <sup>a</sup>	1

a. N1811TL-CFG002 has standard options 124-026-201-302-403.

**Table 1-7**

**E4446A, and E4448A Option H26 Replaceable Parts**

Description	Agilent Part Number	Quantity
Microwave amplifier	0955-1617	1
Transfer switch 50 GHz	87222-60015	1
Attenuator bracket	E4440-00020	1
Cable Assembly, ribbon, 10 conductor	E4440-60071	1
Board Assembly, PSA H26	E4440-60358	1
Bracket, Option H26	E4446-00008	1
Semi-rigid cable assembly (Switch In)	E4446-20060	1
Semi-rigid cable assembly (Amp In)	E4446-20061	1
Semi-rigid cable assembly (Amp Out)	E4446-20062	1
Semi-rigid cable assembly (Switch Out)	E4446-20063	1
Wire harness cable assembly (bias)	E4446-60067	1

## Contacting Agilent

By internet, phone, or fax, get assistance with all your test and measurement needs.

This information supersedes all prior HP contact information.			
<b>Online assistance:</b> <a href="http://www.agilent.com/find/assist">www.agilent.com/find/assist</a>			
<b>Americas</b>			
<b>Brazil</b> (tel) (+55) 11 3351 7012 (fax) (+55) 11 3351 7024	<b>Canada</b> (tel) +1 877 894 4414 (fax) +1 303 662 3369	<b>Mexico</b> (tel) 1 800 254 2440 (fax) 1 800 254 4222	<b>United States</b> (tel) 800 829 4444 (alt) (+1) 303 662 3998 (fax) 800 829 4433
<b>Asia Pacific and Japan</b>			
<b>Australia</b> (tel) 1 800 225 574 (fax) 1 800 681 776 (fax) 1 800 225 539	<b>China</b> (tel) 800 810 0508 (alt) 800 810 0510 (fax) 800 810 0507 (fax) 800 810 0362	<b>Hong Kong</b> (tel) 800 933 229 (fax) 800 900 701	<b>India</b> (tel) 1600 112 626 (fax) 1600 112 727 (fax) 1600 113 040
<b>Japan (Bench)</b> (tel) 0120 32 0119 (alt) (+81) 426 56 7799 (fax) 0120 01 2144	<b>Japan (On-Site)</b> (tel) 0120 802 363 (alt) (+81) 426 56 7498 (fax) (+81) 426 60 8953	<b>Singapore</b> (tel) 1 800 275 0880 (fax) (+65) 6755 1235 (fax) (+65) 6755 1214	<b>South Korea</b> (tel) 080 778 0011 (fax) 080 778 0013
<b>Taiwan</b> (tel) 0800 047 669 (fax) 0800 047 667 (fax) 886 3492 0779	<b>Thailand</b> (tel) 1 800 2758 5822 (alt) (+66) 2267 5913 (fax) 1 800 656 336	<b>Malaysia</b> (tel) 1800 880 399 (fax) 1800 801 054	
<b>Europe</b>			
<b>Austria</b> (tel) 0820 87 44 11* (fax) 0820 87 44 22	<b>Belgium</b> (tel) (+32) (0)2 404 9340 (alt) (+32) (0)2 404 9000 (fax) (+32) (0)2 404 9395	<b>Denmark</b> (tel) (+45) 7013 1515 (alt) (+45) 7013 7313 (fax) (+45) 7013 1555	<b>Finland</b> (tel) (+358) 10 855 2100 (fax) (+358) (0) 10 855 2923
<b>France</b> (tel) 0825 010 700* (alt) (+33) (0)1 6453 5623 (fax) 0825 010 701*	<b>Germany</b> (tel) 01805 24 6333* (alt) 01805 24 6330* (fax) 01805 24 6336*	<b>Ireland</b> (tel) (+353) (0)1 890 924 204 (alt) (+353) (0)1 890 924 206 (fax) (+353) (0)1 890 924 024	<b>Israel</b> (tel) (+972) 3 9288 500 (fax) (+972) 3 9288 501
<b>Italy</b> (tel) (+39) (0)2 9260 8484 (fax) (+39) (0)2 9544 1175	<b>Luxemburg</b> (tel) (+32) (0)2 404 9340 (alt) (+32) (0)2 404 9000 (fax) (+32) (0)2 404 9395	<b>Netherlands</b> (tel) (+31) (0)20 547 2111 (alt) (+31) (0)20 547 2000 (fax) (+31) (0)20 547 2190	<b>Russia</b> (tel) (+7) 095 797 3963 (alt) (+7) 095 797 3900 (fax) (+7) 095 797 3901
<b>Spain</b> (tel) (+34) 91 631 3300 (alt) (+34) 91 631 3000 (fax) (+34) 91 631 3301	<b>Sweden</b> (tel) 0200 88 22 55* (alt) (+46) (0)8 5064 8686 (fax) 020 120 2266*	<b>Switzerland (French)</b> (tel) 0800 80 5353 opt. 2* (alt) (+33) (0)1 6453 5623 (fax) (+41) (0)22 567 5313	<b>Switzerland (German)</b> (tel) 0800 80 5353 opt. 1* (alt) (+49) (0)7031 464 6333 (fax) (+41) (0)1 272 7373
<b>Switzerland (Italian)</b> (tel) 0800 80 5353 opt. 3* (alt) (+39) (0)2 9260 8484 (fax) (+41) (0)22 567 5314	<b>United Kingdom</b> (tel) (+44) (0)7004 666666 (alt) (+44) (0)7004 123123 (fax) (+44) (0)7004 444555		
(tel) = primary telephone number; (alt) = alternate telephone number; (fax) = FAX number; * = in country number      11/16/04			

2

# Example Measurements

## What is in This Chapter

This chapter demonstrates analog and digital analyzer measurements with the Agilent E444xA Option H26 and how it relates to a signal with the preamp “ON” or “OFF.”

- “[Equipment for Test Set-ups](#)” on page 2-3.
- “[Analog Test Set-up](#)” on page 2-4.
- “[Digital Test Set-up](#)” on page 2-6.

To find descriptions of specific analyzer functions, refer to the *Agilent Technologies PSA Series Spectrum Analyzers User’s and Programmer’s Reference Guide*.

## Equipment for Test Set-ups

Use the following equipment to set-up the Option H26 for an analog or digital measurement.

**Table 2-1** Required Test Equipment


Test Equipment	Quantity	Recommended Model
Various Adapters	n/a	n/a
Various cables	n/a	n/a
Synthesizer Sweeper 10 MHz to 50 GHz	1	83650B <sup>1</sup>
Attenuator 20 dB or greater	1	8490D

1. An 83630B can be used with the E4440A, E4443A, or E4445A.

## Analog Test Set-up

Using the Agilent Technologies E444xA Option H26 will allow the user to view low amplitude analog signals. Use the following instrument set-ups for viewing an analog signal. Use the “[Front Panel Key Select](#)” on page 3-3 and “[GPIB Control of Preamp Gain](#)” on page 3-4 to control the instrument while making measurements.


**Figure 2-1**      **Analog Test Set-up**



The setup shown in [Figure 2-1, “Analog Test Set-up”](#) can be used to look for low amplitude signals that would normally be within the noise floor of the spectrum analyzer.


In the example, a low amplitude, 26 GHz signal is generated with the Agilent 83650B. The amplitude was further reduced by using an 8490D coaxial attenuator, which produced a signal that was hidden in the noise floor of the spectrum analyzer in normal operation. Refer to [Figure 2-2, “Analog Input Signal with Preamp OFF.”](#)

**Figure 2-2**      **Analog Input Signal with Preamp OFF**



When the internal preamp is turned on, the noise floor of the spectrum analyzer is reduced, revealing the signal of interest. Refer to [Figure 2-3, “Analog Input Signal with Preamp ON.”](#)


**Figure 2-3**      **Analog Input Signal with Preamp ON**



## Digital Test Set-up

Using the Agilent Technologies E444xA Option H26 will allow the user to view low amplitude digital signals. Use the following instrument set-up to view an example of a digital signal. Use the “[Front Panel Key Select](#)” on page 3-3 and “[GPIB Control of Preamp Gain](#)” on page 3-4 to control the instrument while making measurements.

**Figure 2-4**      **Digital Test Set-up**



The test set-up shown in [Figure 2-4, “Digital Test Set-up”](#) can also be used to look for a low amplitude modulated signals. In the example, we try to see the carrier suppression of a multi-tone digital signal. The 20 GHz signal was generated by an Agilent E8267C Vector Signal Generator.


The E8267C was configured using the following keystrokes:

**[AMPLITUDE] > -60 > dBm [MODE] > Multitone > On**

The default settings of the source were used to produce the signal as shown. The amplitude was reduced even further by using a coaxial attenuator, model number 8490D.


This produced the multi-tone signal as shown in [Figure 2-5, “Digital Input Signal with Preamp OFF.”](#) Notice that we cannot detect the carrier signal at 20 GHz.

**Figure 2-5      Digital Input Signal with Preamp OFF**



When the Option H26 highband preamplifier is turned on, the noise floor of the spectrum analyzer is reduced thus revealing the carrier signal at 20 GHz. Refer to [Figure 2-6, “Digital Input Signal with Preamp ON.”](#)

**Figure 2-6      Digital Input Signal with Preamp ON**



**Example Measurements**  
**Equipment for Test Set-ups**

---

---

**3**

## **User Interface**

## What is in This Chapter

This chapter provides instructions on using the preamp via the front panel, through GPIB commands, accessing the flatness compensation table reloading and recreating the flatness correction table for your instrument.

- “Front Panel Key Select” on page 3-3.
- “GPIB Control of Preamp Gain” on page 3-4
- “Flatness Compensation Table” on page 3-5.

To find descriptions of specific analyzer functions, refer to the *Agilent Technologies PSA Series Spectrum Analyzers User’s and Programmer’s Reference Guide*.

## Front Panel Key Select

### AMPLITUDE Y Scale

Activates the reference level function and accesses the amplitude menu keys only while the instrument is in Spectrum Analysis mode. Amplitude menu keys allow you to set functions that affect the way data on the vertical axis is displayed or corrected.

#### Int Preamp

On Off

*Agilent E444xA with Option 1DS and H26* turns the internal preamps on and off. Pressing **Int Preamp (On)** results in a correction being applied to compensate for the gain of the 1DS preamp so that amplitude readings show the value at the input connector. To apply corrections to the H26 Highband Preamp amplitude readings, corrections are manually applied and are explained in the section titled “[Accessing the Flatness Compensation Table](#)” on page 3-5. When the preamps are on, a PA indication appears on the left side of the display.

The lowband preamp is specified to operate over a 1 MHz to 3 GHz frequency range and the highband preamp is specified to operate from 3 GHz to the upper frequency limit of the instrument<sup>1</sup> depending on the model. Refer to [Figure 1-1, “Lowband vs. Highband Frequency Diagram for E4440A”](#) on page 1-3.

Key Access: **AMPLITUDE Y Scale > More 1 of 3**

---

1. Refer to [Table 1-1 on page 1-2](#) for other model numbers and frequency ranges.

## GPIB Control of Preamp Gain

The following SCPI commands apply to the Agilent E444xA Option 1DS and H26.

Refer to *Agilent Technologies PSA Spectrum Analyzer User's and Programmer's Reference Guide*, to locate SCPI command subsystems and subsections that apply to the standard functions of the Agilent Technologies E444xA PSA-Series Spectrum Analyzer.

### Preamp Gain Command Interface

```
[:SENSe]:POWeR[:RF]:GAIN[:STATe] OFF|ON|0|1  
[:SENSe]:POWeR[:RF]:GAIN[:STATe]?
```

The above two commands turns the internal lowband and highband preamp on or off.

Factory Preset  
and \*RST: Off

Remarks: This command is available only with Option 1DS installed.

Front Panel  
Access: **AMPLITUDE/Y Scale > More 1 of 3 > Int Preamp On Off**

Refer to [Figure 2-3, “Analog Input Signal with Preamp ON” on page 2-5](#) for an example of an analog "ON" state of the highband preamp.

Refer to [Figure 2-6, “Digital Input Signal with Preamp ON” on page 2-7](#) for an example of a digital "ON" state of the highband preamp.

## Flatness Compensation Table

The following sections will cover the following:

- “Accessing the Flatness Compensation Table” on page 3-5.
- “GPIB Control of the Flatness Compensation Table” on page 3-5
- “Reloading the Flatness Correction Table” on page 3-6.
- “Recreating the Flatness Correction Table” on page 3-7.

### Accessing the Flatness Compensation Table

During the Swept Flatness test routine at the factory, correction data points are collected with the highband preamp in the OFF state. The test is then run again with the highband preamp ON. The corrections are then applied to get the best flatness response.

The flatness compensation table is available for the user to apply corrections for an improved flatness response while the highband preamp is on. The following keystrokes will aid the user in accessing the flatness compensation table.

**[AMPLITUDE] > More 1 of 3 > Corrections > Other > Correction > On | Off**

---

**NOTE** When sweeping from below 2.85 GHz to above 3.05 GHz, the corrections table will cause amplitude errors in the 2.85 GHz – 3 GHz band. To avoid these errors, set the start frequency above 2.85 GHz.

---

**NOTE** The instrument does not automatically enable the Amplitude Corrections, they must be enabled by the user where appropriate. Unlike the option 1DS lowband preamp where the lowband preamp flatness is corrected automatically, the Amplitude Correction “**OTHER**” must be turned “ON” when the Preamp is turned “ON” and turned “OFF” when the Preamp is turned “OFF.” Measurement errors will result if this is not done.

### GPIB Control of the Flatness Compensation Table

The following SCPI commands related to the Agilent E444xA Option 1DS and H26.

Refer to *Agilent Technologies PSA Spectrum Analyzer User’s and Programmer’s Reference Guide*, to locate SCPI command subsystems and subsections that apply to the standard functions of the Agilent Technologies E444xA PSA-Series Spectrum Analyzer.

## Preamp Compensation Table Command Interface

The following command turns on or off the correction system.

[**:SENSe**] :CORRection:CSET:ALL[:STATe] OFF|ON|0|1

The following command turns on or off the correction file to be used<sup>1</sup>.

[**:SENSe**] :CORRection:CSET3[:STATe] OFF|ON|0|1

## Reloading the Flatness Correction Table

As a cautionary measure, the flatness corrections stored in the “Other” users file have been provided on a matching serialized disk to be reloaded into the specific spectrum analyzer in the event the flatness corrections are deleted from the instrument.

Insert the serialized disk into the instrument. Use the following keystrokes to reload the serialized flatness corrections data from the disk to the specific instrument.

**[FILE] > Load > Type > More 1 of 2 > Corrections > Dir Select > [↑ | ↓] > Drive A > [↑ | ↓] > "H26CORR" > Load Now**

## Verification of the Flatness Correction Data

To verify the download, use the following keystrokes then verify that the data is visible in the table.

**[AMPLITUDE] > More 1 of 3 > Corrections > Other > Edit**

Press **[RETURN]** so as not to disturb the data.

---

1. where CSET3 is the correction set value for "Other".

## Recreating the Flatness Correction Table

This section is provided in the event that both the disk and the internal correction tables become corrupt or lost. This procedure can be used to re-create the flatness correction table.

**Table 3-1** Required Test Equipment (recreating)

Test Equipment	Quantity	Recommended Model
Various Adapters	n/a	n/a
Various cables	n/a	n/a
Synthesizer Sweeper 10 MHz to 50 GHz	1	83650B <sup>1</sup>

1. An 83630B can be used for the E4440A, E4443A or E4445A.

## Connections

1. Connect the 10 MHz Reference Output from the PSA to the Reference In of the source.
2. Turn on the Reference Output on the PSA.

[System] > Reference > 10 MHz Out > On

3. Use [Table 3-2](#) for the PSA settings.

**Table 3-2** PSA Instrument Settings

Setting	Value
Span	2 MHz
Attenuation	10 dB <sup>1</sup>
Ref Level	-30 dBm
Scale/Div	5 dB

1. **MUST be in Manual Mode:**  
 [AMPLITUDE] > Atten > Manual > 10 dB

4. Turn off the corrections on the PSA.

**[AMPLITUDE] > More 1 of 3 > Corrections > Apply Corrections > Off**

5. Use [Table 3-3](#) for the Source settings.

**Table 3-3**      **Source Instrument Settings**

Setting	Value
Frequency	follows PSA frequency
Power Level	-35 dBm

### Procedure

1. Select 199 frequency points for use in the correction table. The first point should be 2.85 GHz and the final point should be the upper frequency limit of the PSA under test.

The PSA and the source will be set to each frequency chosen.

Repeat the following sequence of steps for each frequency on the PSA and Source in generating the replacement correction table. The PSA and the source will be set to each frequency chosen.

- **[AMPLITUDE] > More 1 of 3 > Int Preamp Off**

Adjust the source to display the input signal on the PSA screen. This should be at least 20 dB above the noise floor.

- **[Peak Search]**
- **[Marker] > Delta**
- **[AMPLITUDE] > More 1 of 3 > Int Preamp On**

Record the amplitude delta. Refer to [Table 3-4 on page 3-9](#) for a sample table. The use of a spreadsheet is helpful in recording the points.

---

**NOTE** Any adjustment points  $\leq$  3 GHz will need to be measured with the PSA start frequency of 2.85 GHz and a stop frequency of 3.06 GHz in order to insure that the instrument is in band 1.

---

**Table 3-4**

**A Sample Frequency Amplitude Sheet**

Point	Frequency (GHz)	Delta Amplitude (dB)	Correction (dBm <sup>1</sup> )
1	2.849 <sup>2</sup>	0 dB	0 dB
2	2.85	___ dB	- (___ dB)
...	___ GHz	___ dB	- (___ dB)
200	___ GHz <sup>3</sup>	___ dB	- (___ dB)

1. The correction field will contain the negative, or opposite of the Amplitude field. The correction value is what the instrument uses rather than the Amplitude (or error) value.
2. This point is needed so that the corrections are not applied to Band 0.
3. This should be the upper limit of the PSA.

### Entering the Flatness Correction Data

To manually enter these recorded values into the PSA instrument, perform the following sequence for each point.

- [AMPLITUDE] > More 1 of 3 > Corrections > Other > Edit
- Select [FREQUENCY] 2.849 GHz
- Select [AMPLITUDE] 0 dB

After each entry, the Point counter automatically increments for entry of the next point.

- Select [FREQUENCY] 2.85 GHz
- Select [AMPLITUDE] \_\_\_ dB (Refer to Table 3-4)
- Continue entering data for each point in the table.

Press [RETURN], when finished.

To use the new correction table, insure that the corrections are turned on.

---

**CAUTION**

It is recommended that the correction set entered be saved on the internal drive or a floppy disk for future references. See the *PSA User's and Programming Reference Guide* for information on saving correction values.

---

User Interface  
**Flatness Compensation Table**

4

## Appendix

## **DANL Manual Performance Verification**

The Displayed Average Noise Level manual performance verification measures the preamp “off” DANL of Option H26. The test measures the noise in zero span with a 1 kHz resolution bandwidth, and then normalizes the amplitude to a 1 Hz bandwidth.

DANL is defined as the average of the displayed trace. There is no practical method for manually reading the average of the trace. This procedure averages the trace 100 times and then the operator scrolls the display line to the middle of the trace. The display line reading is considered the trace average. The reading is normalized to a 1 Hz RBW by subtracting 30 dB ( $10 \times \log(1000)$ ) from the display line value.

[Table 4-7 on page 4-7](#) directs the user to a specific Test Record for that instrument. [Table 4-8 on page 4-8](#) is provided for traceability of the instruments used in the testing.

### **Procedure**

1. Connect the appropriate termination to the PSA RF input connector.  
Press **[System] > Alignments > Align All Now**.
2. Press **Preset**. Set the analyzer controls by pressing the following keys:

**SPAN, Span, 0 Hz**  
**Sweep, Sweep Time, 20 ms**  
**Amplitude, Ref Level, -100 dBm**  
**Amplitude, Attenuation, Man, 0 dB**  
**BW/AVG, Res BW, 1 kHz**  
**BW/AVG, Average, On**  
**Single**

3. Press **[AMPLITUDE] > More 1 of 3 > In Preamp > On**. (This is the lowband preamp.)
4. Press **[FREQUENCY]** and enter the first frequency listed in [Table 4-1 on page 4-3](#).

**Table 4-1**

**DANL (PSA with Option 1DS)**

Frequency	Normalized DANL
100 kHz	1)
199 kHz	2)
201 kHz	3)
499 kHz	4)
501 kHz	5)
9.9 MHz	6)
10.1 MHz	7)
1.0 GHz	8)
1.2 GHz	9)
2.4 GHz	10)
2.6 GHz	11)
3.0 GHz	12)

5. Press **Single** and wait for 100 averages.
6. Press **Display, Display Line, On**.
7. Scroll the display line so that it bisects the trace. Read the display line amplitude and subtract 30 from the value. Record the result under “Normalized DANL” in the table.
8. Repeat **step 5** through **step 8** for all frequencies listed in [Table 4-1](#).
9. Record the results of [Table 4-1](#) in the test record table for your model instrument later in this chapter. Refer to [Table 4-7 on page 4-7](#).
10. Press **[AMPLITUDE], More 1 of 3, Int Preamp, Off**.
11. Press **[FREQUENCY]** and enter the first frequency listed in [Table 4-2](#).

**Table 4-2**

**DANL (All PSA Instruments)**

Frequency	Normalized DANL
10 kHz	13)
99 kHz	14)
101 kHz	15)
990 kHz	16)
1.01 MHz	17)
500 MHz	18)

**Table 4-2 DANL (All PSA Instruments)**

Frequency	Normalized DANL
1.1 GHz	19)
1.3 GHz	20)
2.0 GHz	21)
2.4 GHz	22)
2.6 GHz	23)
3.0 GHz	24)
3.1 GHz	25)
4.0 GHz	26)
5.0 GHz	27)
6.5 GHz	28)
6.7 GHz	29)

12. Press **Single** and wait for 100 averages.
13. Press **Display, Display Line, On**.
14. Scroll the display line so that it bisects the trace. Read the display line amplitude and subtract 30 from the value. Record the result under “Normalized DANL” in the table.
15. Repeat step 12 through step 15 for all frequencies listed in [Table 4-2](#).
16. Record the results of [Table 4-2](#) in the test record table for your model instrument later in this chapter. Refer to [Table 4-7 on page 4-7](#). If the analyzer is an E4443A stop here.
17. Repeat step 12 through step 15 for all frequencies listed in [Table 4-3](#).

**Table 4-3 DANL (PSA E4445A, E4440A, E4446A, E4448A)**

Frequency	Normalized DANL
8.0 GHz	30)
9.0 GHz	31)
10.0 GHz	32)
11.0 GHz	33)
12.0 GHz	34)
13.1 GHz	35)

18. Record the results of [Table 4-3](#) in the test record table for your model instrument later in this chapter. Refer to [Table 4-7 on page 4-7](#). If the analyzer is an E4445A stop here.

19. Repeat [step 12](#) through [step 15](#) for all frequencies listed in [Table 4-4](#).

**Table 4-4 DANL (PSA E4440A, E4446A, E4448A)**

Frequency	Normalized DANL
13.3 GHz	<a href="#">36)</a>
14.0 GHz	<a href="#">37)</a>
15.0 GHz	<a href="#">38)</a>
16.0 GHz	<a href="#">39)</a>
17.0 GHz	<a href="#">40)</a>
18.0 GHz	<a href="#">41)</a>
19.9 GHz	<a href="#">42)</a>
20.1 GHz	<a href="#">43)</a>
21.0 GHz	<a href="#">44)</a>
22.4 GHz	<a href="#">45)</a>
22.6 GHz	<a href="#">46)</a>
24.0 GHz	<a href="#">47)</a>
25.0 GHz	<a href="#">48)</a>
26.4 GHz	<a href="#">49)</a>

20. Record the results of [Table 4-4](#) in the test record table for your model instrument later in this chapter. Refer to [Table 4-7 on page 4-7](#). If the analyzer is an E4440A stop here.

21. Repeat [step 12](#) through [step 15](#) for all frequencies listed in [Table 4-5](#).

**Table 4-5 DANL (PSA E4446A, E4448A)**

Frequency	Normalized DANL
27.0 GHz	<a href="#">50)</a>
28.0 GHz	<a href="#">51)</a>
29.0 GHz	<a href="#">52)</a>
30.0 GHz	<a href="#">53)</a>
31.0 GHz	<a href="#">54)</a>
32.0 GHz	<a href="#">55)</a>
33.0 GHz	<a href="#">56)</a>
34.0 GHz	<a href="#">57)</a>
35.5 GHz	<a href="#">58)</a>

**Table 4-5 DANL (PSA E4446A, E4448A)**

Frequency	Normalized DANL
36.5 GHz	<b>59)</b>
37.5 GHz	<b>60)</b>
38.5 GHz	<b>61)</b>
39.0 GHz	<b>62)</b>
40.0 GHz	<b>63)</b>
41.0 GHz	<b>64)</b>
42.0 GHz	<b>65)</b>
43.0 GHz	<b>66)</b>
43.9 GHz	<b>67)</b>

22. Record the results of [Table 4-5](#) in the test record table for your model instrument later in this chapter. Refer to [Table 4-7 on page 4-7](#). If the analyzer is an E4446A stop here.  
23. Repeat step 12 through step 15 for all frequencies listed in [Table 4-6](#).

**Table 4-6 DANL (PSA E4448A)**

Frequency	Normalized DANL
44.1 GHz	<b>68)</b>
45.0 GHz	<b>69)</b>
46.0 GHz	<b>70)</b>
47.0 GHz	<b>71)</b>
48.0 GHz	<b>72)</b>
49.0 GHz	<b>73)</b>
50.0 GHz	<b>74)</b>

24. Record the results of [Table 4-6](#) in the test record table for your model instrument later in this chapter. Refer to [Table 4-7 on page 4-7](#).

**Table 4-7**

**Model Number Verification Test Record**

<b>Model Number</b>	<b>Table Number</b>
E4440A	<a href="#">Table 4-9 on page 4-10</a>
E4443A	<a href="#">Table 4-10 on page 4-13</a>
E4445A	<a href="#">Table 4-11 on page 4-15</a>
E4446A	<a href="#">Table 4-12 on page 4-17</a>
E4448A	<a href="#">Table 4-13 on page 4-20</a>

**Table 4-8 Agilent E444xA Performance Verification Test Record**

Agilent Technologies			
Address: _____	Report No. _____		
	Date _____		
Model E4440A			
Serial No. _____	Ambient temperature _____ °C		
Options _____	Relative humidity _____ %		
Firmware Revision _____	Power mains line frequency _____ Hz (nominal)		
Customer _____	Tested by _____		
<b>Test Equipment Used:</b>			
Description	Model No.	Trace No.	Cal Due Date
Synthesized Signal Generator	_____	_____	_____
Synthesized Signal Source	_____	_____	_____
Synthesized Sweeper #1	_____	_____	_____
Synthesized Sweeper #2	_____	_____	_____
Function Generator	_____	_____	_____
Power Meter, Dual-Channel	_____	_____	_____
RF Power Sensor #1	_____	_____	_____
RF Power Sensor #2	_____	_____	_____
Microwave Power Sensor	_____	_____	_____
Digital Multimeter	_____	_____	_____
Universal Counter	_____	_____	_____
Frequency Standard	_____	_____	_____
RF Power Splitter	_____	_____	_____
Microwave Power Splitter	_____	_____	_____
50 Ω Termination	_____	_____	_____
1 dB Step Attenuator	_____	_____	_____
10 dB Step Attenuator	_____	_____	_____
6 dB Fixed Attenuator	_____	_____	_____
10 dB Fixed Attenuator	_____	_____	_____

**Table 4-8 Agilent E444xA Performance Verification Test Record**

20 dB Fixed Attenuator	_____	_____	_____
Directional Bridge	_____	_____	_____
Directional Coupler	_____	_____	_____
Notes/comments:	_____		
	_____		

**Table 4-9 Agilent E4440A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4440A</b>		Report No. _____		
Serial No. _____		Date _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
<b>Displayed Average Noise Level</b>		Note: Enter results with preamp on in the appropriate section based upon the ambient temperature when the test was performed.		
Preamp On	100 kHz	(1)_____	-159 dBm	±0.07 dB
	199 kHz	(2)_____	-159 dBm	±0.07 dB
	201 kHz	(3)_____	-159 dBm	±0.07 dB
	499 kHz	(4)_____	-159 dBm	±0.07 dB
	501 kHz	(5)_____	-163 dBm	±0.07 dB
	9.9 MHz	(6)_____	-166 dBm	±0.07 dB
	10.1 MHz	(7)_____	-169 dBm	±0.07 dB
	1.0 GHz	(8)_____	-168 dBm	±0.07 dB
	1.2 GHz	(9)_____	-167 dBm	±0.07 dB
	2.4 GHz	(10)_____	-165 dBm	±0.07 dB
	2.6 GHz	(11)_____	-165 dBm	±0.07 dB
	3.0 GHz	(12)_____	-165 dBm	±0.07 dB
Preamp Off	10 kHz	(13)_____	-134 dBm	±0.07 dB
	99 kHz	(14)_____	-134 dBm	±0.07 dB
	101 kHz	(15)_____	-144 dBm	±0.07 dB
	990 kHz	(16)_____	-144 dBm	±0.07 dB
	1.01 MHz	(17)_____	-149 dBm	±0.07 dB
	500 MHz	(18)_____	-154 dBm	±0.07 dB
	1.1 GHz	(19)_____	-154 dBm	±0.07 dB
	1.3 GHz	(20)_____	-153 dBm	±0.07 dB
	2.0 GHz	(21)_____	-153 dBm	±0.07 dB
	2.4 GHz	(22)_____	-153 dBm	±0.07 dB

**Table 4-9 Agilent E4440A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4440A</b>		Report No. _____		
Serial No. _____		Date _____		
Test Description	Minimum	Results Measured	Maximum	Measurement Uncertainty
2.6 GHz		(23)_____	-152 dBm	±0.07 dB
3.0 GHz		(24)_____	-152 dBm	±0.07 dB
3.1 GHz		(25)_____	-151 dBm	±0.07 dB
4.0 GHz		(26)_____	-151 dBm	±0.07 dB
5.0 GHz		(27)_____	-151 dBm	±0.07 dB
6.5 GHz		(28)_____	-151 dBm	±0.07 dB
6.7 GHz		(29)_____	-149 dBm	±0.07 dB
8.0 GHz		(30)_____	-149 dBm	±0.07 dB
9.0 GHz		(31)_____	-149 dBm	±0.07 dB
10.0 GHz		(32)_____	-149 dBm	±0.07 dB
11.0 GHz		(33)_____	-149 dBm	±0.07 dB
12.0 GHz		(34)_____	-149 dBm	±0.07 dB
13.1 GHz		(35)_____	-149 dBm	±0.07 dB
13.3 GHz		(36)_____	-146 dBm	±0.07 dB
14.0 GHz		(37)_____	-146 dBm	±0.07 dB
15.0 GHz		(38)_____	-146 dBm	±0.07 dB
16.0 GHz		(39)_____	-146 dBm	±0.07 dB
17.0 GHz		(40)_____	-146 dBm	±0.07 dB
18.0 GHz		(41)_____	-146 dBm	±0.07 dB
19.9 GHz		(42)_____	-146 dBm	±0.07 dB
20.1 GHz		(43)_____	-146 dBm	±0.07 dB
21.0 GHz		(44)_____	-142 dBm	±0.07 dB
22.4 GHz		(45)_____	-142 dBm	±0.07 dB
22.6 GHz		(46)_____	-141 dBm	±0.07 dB
24.0 GHz		(47)_____	-141 dBm	±0.07 dB

Appendix  
DANL Manual Performance Verification

**Table 4-9 Agilent E4440A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4440A</b>		Report No. _____		
Serial No. _____		Date _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
	25.0 GHz	(48)_____	-141 dBm	±0.07 dB
	26.4 GHz	(49)_____	-141 dBm	±0.07 dB

**Table 4-10 Agilent E4443A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4443A</b>		Report No. _____		
<b>Serial No.</b> _____		Date _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
<b>Displayed Average Noise Level</b>		Note: Enter results with preamp on in the appropriate section based upon the ambient temperature when the test was performed.		
Preamp On				
	100 kHz	(1)_____	-159 dBm	±0.07 dB
	199 kHz	(2)_____	-159 dBm	±0.07 dB
	201 kHz	(3)_____	-159 dBm	±0.07 dB
	499 kHz	(4)_____	-159 dBm	±0.07 dB
	501 kHz	(5)_____	-163 dBm	±0.07 dB
	9.9 MHz	(6)_____	-166 dBm	±0.07 dB
	10.1 MHz	(7)_____	-169 dBm	±0.07 dB
	1.0 GHz	(8)_____	-168 dBm	±0.07 dB
	1.2 GHz	(9)_____	-167 dBm	±0.07 dB
	2.4 GHz	(10)_____	-165 dBm	±0.07 dB
	2.6 GHz	(11)_____	-165 dBm	±0.07 dB
	3.0 GHz	(12)_____	-165 dBm	±0.07 dB
Preamp Off				
	10 kHz	(13)_____	-134 dBm	±0.07 dB
	99 kHz	(14)_____	-134 dBm	±0.07 dB
	101 kHz	(15)_____	-144 dBm	±0.07 dB
	990 kHz	(16)_____	-144 dBm	±0.07 dB
	1.01 MHz	(17)_____	-149 dBm	±0.07 dB
	500 MHz	(18)_____	-154 dBm	±0.07 dB
	1.1 GHz	(19)_____	-154 dBm	±0.07 dB
	1.3 GHz	(20)_____	-153 dBm	±0.07 dB
	2.0 GHz	(21)_____	-153 dBm	±0.07 dB
	2.4 GHz	(22)_____	-153 dBm	±0.07 dB

**Appendix**  
**DANL Manual Performance Verification**

**Table 4-10 Agilent E4443A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4443A</b>		<b>Report No.</b> _____		
<b>Serial No.</b> _____		<b>Date</b> _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
2.6 GHz		(23)_____	-152 dBm	±0.07 dB
3.0 GHz		(24)_____	-152 dBm	±0.07 dB
3.1 GHz		(25)_____	-151 dBm	±0.07 dB
4.0 GHz		(26)_____	-151 dBm	±0.07 dB
5.0 GHz		(27)_____	-151 dBm	±0.07 dB
6.5 GHz		(28)_____	-151 dBm	±0.07 dB
6.7 GHz		(29)_____	-149 dBm	±0.07 dB

**Table 4-11 Agilent E4445A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4445A</b>		Report No. _____		
Serial No. _____		Date _____		
Test Description	Minimum	Results Measured	Maximum	Measurement Uncertainty
<b>Displayed Average Noise Level</b>		Note: Enter results with preamp on in the appropriate section based upon the ambient temperature when the test was performed.		
Preamp On				
	100 kHz	(1)_____	-159 dBm	±0.07 dB
	199 kHz	(2)_____	-159 dBm	±0.07 dB
	201 kHz	(3)_____	-159 dBm	±0.07 dB
	499 kHz	(4)_____	-159 dBm	±0.07 dB
	501 kHz	(5)_____	-163 dBm	±0.07 dB
	9.9 MHz	(6)_____	-166 dBm	±0.07 dB
	10.1 MHz	(7)_____	-169 dBm	±0.07 dB
	1.0 GHz	(8)_____	-168 dBm	±0.07 dB
	1.2 GHz	(9)_____	-167 dBm	±0.07 dB
	2.4 GHz	(10)_____	-165 dBm	±0.07 dB
	2.6 GHz	(11)_____	-165 dBm	±0.07 dB
	3.0 GHz	(12)_____	-165 dBm	±0.07 dB
Preamp Off				
	10 kHz	(13)_____	-134 dBm	±0.07 dB
	99 kHz	(14)_____	-134 dBm	±0.07 dB
	101 kHz	(15)_____	-144 dBm	±0.07 dB
	990 kHz	(16)_____	-144 dBm	±0.07 dB
	1.01 MHz	(17)_____	-149 dBm	±0.07 dB
	500 MHz	(18)_____	-154 dBm	±0.07 dB
	1.1 GHz	(19)_____	-154 dBm	±0.07 dB
	1.3 GHz	(20)_____	-153 dBm	±0.07 dB
	2.0 GHz	(21)_____	-153 dBm	±0.07 dB
	2.4 GHz	(22)_____	-153 dBm	±0.07 dB

**Appendix**  
**DANL Manual Performance Verification**

**Table 4-11 Agilent E4445A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4445A</b>		<b>Report No.</b> _____		
<b>Serial No.</b> _____		<b>Date</b> _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
2.6 GHz		(23)_____	-152 dBm	±0.07 dB
3.0 GHz		(24)_____	-152 dBm	±0.07 dB
3.1 GHz		(25)_____	-151 dBm	±0.07 dB
4.0 GHz		(26)_____	-151 dBm	±0.07 dB
5.0 GHz		(27)_____	-151 dBm	±0.07 dB
6.5 GHz		(28)_____	-151 dBm	±0.07 dB
6.7 GHz		(29)_____	-149 dBm	±0.07 dB
8.0 GHz		(30)_____	-149 dBm	±0.07 dB
9.0 GHz		(31)_____	-149 dBm	±0.07 dB
10.0 GHz		(32)_____	-149 dBm	±0.07 dB
11.0 GHz		(33)_____	-149 dBm	±0.07 dB
12.0 GHz		(34)_____	-149 dBm	±0.07 dB
13.1 GHz		(35)_____	-149 dBm	±0.07 dB

**Table 4-12 Agilent E4446A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4446A</b>		Report No. _____		
Serial No. _____		Date _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
<b>Displayed Average Noise Level</b>		Note: Enter results with preamp on in the appropriate section based upon the ambient temperature when the test was performed.		
Preamp On	100 kHz	(1)_____	-158 dBm	±0.07 dB
	199 kHz	(2)_____	-158 dBm	±0.07 dB
	201 kHz	(3)_____	-158 dBm	±0.07 dB
	499 kHz	(4)_____	-159 dBm	±0.07 dB
	501 kHz	(5)_____	-161 dBm	±0.07 dB
	9.9 MHz	(6)_____	-167 dBm	±0.07 dB
	10.1 MHz	(7)_____	-167 dBm	±0.07 dB
	1.0 GHz	(8)_____	-166 dBm	±0.07 dB
	1.2 GHz	(9)_____	-166 dBm	±0.07 dB
	2.4 GHz	(10)_____	-163 dBm	±0.07 dB
	2.6 GHz	(11)_____	-163 dBm	±0.07 dB
	3.0 GHz	(12)_____	-163 dBm	±0.07 dB
	Preamp Off			
Preamp Off	10 kHz	(13)_____	-136 dBm	±0.07 dB
	99 kHz	(14)_____	-136 dBm	±0.07 dB
	101 kHz	(15)_____	-144 dBm	±0.07 dB
	990 kHz	(16)_____	-144 dBm	±0.07 dB
	1.01 MHz	(17)_____	-149 dBm	±0.07 dB
	500 MHz	(18)_____	-153 dBm	±0.07 dB
	1.1 GHz	(19)_____	-153 dBm	±0.07 dB
	1.3 GHz	(20)_____	-152 dBm	±0.07 dB
	2.0 GHz	(21)_____	-152 dBm	±0.07 dB
	2.4 GHz	(22)_____	-151 dBm	±0.07 dB
	2.6 GHz	(23)_____	-151 dBm	±0.07 dB

**Appendix**  
**DANL Manual Performance Verification**

**Table 4-12 Agilent E4446A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4446A</b>		<b>Report No.</b> _____		
<b>Serial No.</b> _____		<b>Date</b> _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
3.0 GHz		(24)_____	-151 dBm	±0.07 dB
3.1 GHz		(25)_____	-150 dBm	±0.07 dB
4.0 GHz		(26)_____	-150 dBm	±0.07 dB
5.0 GHz		(27)_____	-150 dBm	±0.07 dB
6.5 GHz		(28)_____	-150 dBm	±0.07 dB
6.7 GHz		(29)_____	-145 dBm	±0.07 dB
8.0 GHz		(30)_____	-145 dBm	±0.07 dB
9.0 GHz		(31)_____	-145 dBm	±0.07 dB
10.0 GHz		(32)_____	-145 dBm	±0.07 dB
11.0 GHz		(33)_____	-145 dBm	±0.07 dB
12.0 GHz		(34)_____	-145 dBm	±0.07 dB
13.1 GHz		(35)_____	-145 dBm	±0.07 dB
13.3 GHz		(36)_____	-144 dBm	±0.07 dB
14.0 GHz		(37)_____	-144 dBm	±0.07 dB
15.0 GHz		(38)_____	-144 dBm	±0.07 dB
16.0 GHz		(39)_____	-144 dBm	±0.07 dB
17.0 GHz		(40)_____	-143 dBm	±0.07 dB
18.0 GHz		(41)_____	-143 dBm	±0.07 dB
19.9 GHz		(42)_____	-143 dBm	±0.07 dB
20.1 GHz		(43)_____	-141 dBm	±0.07 dB
21.0 GHz		(44)_____	-141 dBm	±0.07 dB
22.4 GHz		(45)_____	-141 dBm	±0.07 dB
22.6 GHz		(46)_____	-138 dBm	±0.07 dB
24.0 GHz		(47)_____	-138 dBm	±0.07 dB
25.0 GHz		(48)_____	-138 dBm	±0.07 dB
26.4 GHz		(49)_____	-138 dBm	±0.07 dB
27.0 GHz		(50)_____	-140 dBm	±0.07 dB

**Table 4-12 Agilent E4446A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4446A</b>		Report No. _____		
<b>Serial No.</b> _____		Date _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
	28.0 GHz	(51)_____	-140 dBm	±0.07 dB
	29.0 GHz	(52)_____	-140 dBm	±0.07 dB
	30.0 GHz	(53)_____	-140 dBm	±0.07 dB
	31.0 GHz	(54)_____	-140 dBm	±0.07 dB
	32.0 GHz	(55)_____	-132 dBm	±0.07 dB
	33.0 GHz	(56)_____	-132 dBm	±0.07 dB
	34.0 GHz	(57)_____	-132 dBm	±0.07 dB
	35.5 GHz	(58)_____	-132 dBm	±0.07 dB
	36.5 GHz	(59)_____	-127 dBm	±0.07 dB
	37.5 GHz	(60)_____	-127 dBm	±0.07 dB
	38.5 GHz	(61)_____	-129 dBm	±0.07 dB
	39.0 GHz	(62)_____	-129 dBm	±0.07 dB
	40.0 GHz	(63)_____	-129 dBm	±0.07 dB
	41.0 GHz	(64)_____	-129 dBm	±0.07 dB
	42.0 GHz	(65)_____	-129 dBm	±0.07 dB
	43.0 GHz	(66)_____	-129 dBm	±0.07 dB
	43.9 GHz	(67)_____	-129 dBm	±0.07 dB

**Table 4-13 Agilent E4448A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4448A</b>		Report No. _____		
Serial No. _____		Date _____		
Test Description	Minimum	Results Measured	Maximum	Measurement Uncertainty
<b>Displayed Average Noise Level</b>		Note: Enter results with preamp on in the appropriate section based upon the ambient temperature when the test was performed.		
Preamp On	100 kHz	(1)_____	-158 dBm	±0.07 dB
	199 kHz	(2)_____	-158 dBm	±0.07 dB
	201 kHz	(3)_____	-158 dBm	±0.07 dB
	499 kHz	(4)_____	-159 dBm	±0.07 dB
	501 kHz	(5)_____	-161 dBm	±0.07 dB
	9.9 MHz	(6)_____	-167 dBm	±0.07 dB
	10.1 MHz	(7)_____	-167 dBm	±0.07 dB
	1.0 GHz	(8)_____	-166 dBm	±0.07 dB
	1.2 GHz	(9)_____	-166 dBm	±0.07 dB
	2.4 GHz	(10)_____	-163 dBm	±0.07 dB
	2.6 GHz	(11)_____	-163 dBm	±0.07 dB
	3.0 GHz	(12)_____	-163 dBm	±0.07 dB
	Preamp Off			
Preamp Off	10 kHz	(13)_____	-136 dBm	±0.07 dB
	99 kHz	(14)_____	-136 dBm	±0.07 dB
	101 kHz	(15)_____	-144 dBm	±0.07 dB
	990 kHz	(16)_____	-144 dBm	±0.07 dB
	1.01 MHz	(17)_____	-149 dBm	±0.07 dB
	500 MHz	(18)_____	-153 dBm	±0.07 dB
	1.1 GHz	(19)_____	-153 dBm	±0.07 dB
	1.3 GHz	(20)_____	-152 dBm	±0.07 dB
	2.0 GHz	(21)_____	-152 dBm	±0.07 dB
	2.4 GHz	(22)_____	-151 dBm	±0.07 dB
	2.6 GHz	(23)_____	-151 dBm	±0.07 dB

**Table 4-13 Agilent E4448A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4448A</b>		Report No. _____		
<b>Serial No.</b> _____		Date _____		
<b>Test Description</b>	<b>Minimum</b>	<b>Results Measured</b>	<b>Maximum</b>	<b>Measurement Uncertainty</b>
	3.0 GHz	(24)_____	-151 dBm	±0.07 dB
	3.1 GHz	(25)_____	-150 dBm	±0.07 dB
	4.0 GHz	(26)_____	-150 dBm	±0.07 dB
	5.0 GHz	(27)_____	-150 dBm	±0.07 dB
	6.5 GHz	(28)_____	-150 dBm	±0.07 dB
	6.7 GHz	(29)_____	-145 dBm	±0.07 dB
	8.0 GHz	(30)_____	-145 dBm	±0.07 dB
	9.0 GHz	(31)_____	-145 dBm	±0.07 dB
	10.0 GHz	(32)_____	-145 dBm	±0.07 dB
	11.0 GHz	(33)_____	-145 dBm	±0.07 dB
	12.0 GHz	(34)_____	-145 dBm	±0.07 dB
	13.1 GHz	(35)_____	-145 dBm	±0.07 dB
	13.3 GHz	(36)_____	-144 dBm	±0.07 dB
	14.0 GHz	(37)_____	-144 dBm	±0.07 dB
	15.0 GHz	(38)_____	-144 dBm	±0.07 dB
	16.0 GHz	(39)_____	-144 dBm	±0.07 dB
	17.0 GHz	(40)_____	-143 dBm	±0.07 dB
	18.0 GHz	(41)_____	-143 dBm	±0.07 dB
	19.9 GHz	(42)_____	-143 dBm	±0.07 dB
	20.1 GHz	(43)_____	-141 dBm	±0.07 dB
	21.0 GHz	(44)_____	-141 dBm	±0.07 dB
	22.4 GHz	(45)_____	-141 dBm	±0.07 dB
	22.6 GHz	(46)_____	-138 dBm	±0.07 dB
	24.0 GHz	(47)_____	-138 dBm	±0.07 dB
	25.0 GHz	(48)_____	-138 dBm	±0.07 dB
	26.4 GHz	(49)_____	-138 dBm	±0.07 dB
	27.0 GHz	(50)_____	-140 dBm	±0.07 dB

**Table 4-13 Agilent E4448A Performance Verification Test Record**

<b>Agilent Technologies</b>				
<b>Model E4448A</b>		Report No. _____		
Serial No. _____		Date _____		
Test Description	Minimum	Results Measured	Maximum	Measurement Uncertainty
28.0 GHz		(51)_____	-140 dBm	±0.07 dB
29.0 GHz		(52)_____	-140 dBm	±0.07 dB
30.0 GHz		(53)_____	-140 dBm	±0.07 dB
31.0 GHz		(54)_____	-140 dBm	±0.07 dB
32.0 GHz		(55)_____	-132 dBm	±0.07 dB
33.0 GHz		(56)_____	-132 dBm	±0.07 dB
34.0 GHz		(57)_____	-132 dBm	±0.07 dB
35.5 GHz		(58)_____	-132 dBm	±0.07 dB
36.5 GHz		(59)_____	-127 dBm	±0.07 dB
37.5 GHz		(60)_____	-127 dBm	±0.07 dB
38.5 GHz		(61)_____	-129 dBm	±0.07 dB
39.0 GHz		(62)_____	-129 dBm	±0.07 dB
40.0 GHz		(63)_____	-129 dBm	±0.07 dB
41.0 GHz		(64)_____	-129 dBm	±0.07 dB
42.0 GHz		(65)_____	-129 dBm	±0.07 dB
43.0 GHz		(66)_____	-129 dBm	±0.07 dB
43.9 GHz		(67)_____	-129 dBm	±0.07 dB
44.1 GHz		(68)_____	-126 dBm	±0.07 dB
45.0 GHz		(69)_____	-125 dBm	±0.07 dB
46.0 GHz		(70)_____	-125 dBm	±0.07 dB
47.0 GHz		(71)_____	-125 dBm	±0.07 dB
48.0 GHz		(72)_____	-125 dBm	±0.07 dB
49.0 GHz		(73)_____	-125 dBm	±0.07 dB
50.0 GHz		(74)_____	-124 dBm	±0.07 dB

**A**  
amplitude y scale 3-3  
analog example 2-5  
analog test setup 2-4

**B**  
block diagram 1-13

**C**  
characteristics 1-4

**D**  
DANL 1-4, 1-7  
    manual verification 4-2  
description 1-2  
digital example 2-7  
digital test setup 2-6  
displayed average noise level 1-7

**E**  
equipment for test setup 2-3  
example  
    analog measurement 2-5  
    digital measurement 2-7

**F**  
flatness compensation table  
    front panel access 3-5  
    GPIB control 3-5  
    recreating 3-7  
    reloading via front panel 3-6  
    verifying 3-6  
frequency range, preamp 1-2  
front panel access  
    preamp 3-3  
functional testing  
    test descriptions 2-3  
    test list 2-3

**G**  
gain, preamp 1-4  
GPIB  
    flatness compensation table on/off 3-5  
    preamp on/off 3-4

**I**  
Int Preamp On Off key 3-3  
int preamp on off key 3-3

**M**  
max safe input level 1-4

**N**  
noise figure personality 1-8

**O**  
option 219 1-8

**P**  
parts, replaceable 1-14  
performance verification 1-12  
preamp  
    flatness compensation table  
        front panel 3-5  
        gpib control 3-6  
        recreating 3-7  
        reloading 3-6  
        verifying 3-6  
        front panel access 3-3  
        gpib control on/off 3-4  
    preamp frequency range 1-2  
    preamp gain 1-4  
    preamplifier on/off via gpib 3-4

**R**  
recreating flatness compensation table 3-7  
reloading flatness compensation table 3-6  
replaceable parts 1-14

**S**  
SCPI  
    command subsystems, locating 3-4

**T**  
test equipment list  
    functional testing 2-3

**V**  
verifying flatness compensation table 3-6  
verifying performance 1-12